

INSTRUCTIONS

JNM-PMX60si

HIGH RESOLUTION NMR INSTRUMENT

No. I-E061-4
(NM105179)





**JNM-PMX60SI HIGH RESOLUTION
NMR INSTRUMENT**

TABLE OF CONTENTS

CHAPTER 1 GENERAL DESCRIPTION

1.1	Outline	1- 1
1.2	Specifications	1- 2
1.2.1	Instrument Specifications	1- 2
1.2.2	Installation Requirements	1- 3
1.3	Composition	1- 5
1.3.1	Standard Composition	1- 5
1.3.2	Options	1- 6
1.4	Standard Samples and Consumables	1- 7
1.4.1	Standard PMX60SI Samples	1- 7
1.4.2	Special Order Consumables and Standard Samples	1- 8
1.5	Connecting up Oscilloscope (Optional)	1- 9

CHAPTER 2 DESCRIPTION OF PANELS

2.1	Spectrometer Front Panel	2- 1
2.2	Spectrometer Table Panel	2- 6
2.3	Spectrometer Rear Panel	2- 9
2.4	Magnet Console: Top	2-10
2.5	Magnet Console: Rear	2-11
2.6	Oscilloscope (Optional)	2-12

CHAPTER 3 OPERATION

3.1	Start-up	3- 1
3.2	Sample Tube Insertion	3- 2
3.3	Initial Settings and NMR Signal Observation	3- 5
3.4	Resolution Adjust	3- 8
3.4.1	Recorder Method	3- 8
3.4.2	Oscilloscope Method	3-11
3.4.3	SIGNAL INTENSITY Meter Method	3-13
3.5	Spinning Adjustment	3-15
3.6	Spinning Side Bands (SSB)	3-16
3.7	Waveform Adjustment	3-17
3.8	Signal Expansion and Offset	3-18
3.8.1	Signal Expansion	3-18
3.8.2	Signal Offset	3-21
3.9	Operation Modes	3-23
3.9.1	TMS SET	3-23
3.9.2	LOCK	3-24
3.9.3	RESO	3-26
3.9.4	INTEG	3-28
3.9.5	AMPL SET	3-32
3.9.6	AUTO	3-34
3.9.7	MANUAL	3-35
3.9.8	CRT	3-35

3.10 NMR Lock (Internal Lock)	3-36
3.11 Spin Decoupling	3-40
3.12 Automatic TMS Settings	3-44

CHAPTER 4 MAINTENANCE

4.1 Table of Maintenance Checks	4- 1
4.2 Parts to be Cleaned	4- 1
4.2.1 Console Cover	4- 1
4.2.2 Panel	4- 2
4.2.3 Top of Magnet Console	4- 2
4.2.4 Sample Tube and Rotor	4- 2
4.2.5 Probe	4- 2
4.2.6 Air Filters	4- 2
4.2.7 Recorder Cartridge Pen	4- 3
4.3 Recorder Maintenance	4- 4
4.4 Installation Room Cleanliness	4- 4
4.5 Sweep Width Correction	4- 5

CHAPTER 5 TROUBLESHOOTING AND COUNTERMEASURES

5.1 Absence of Signals	5- 1
5.2 Deterioration of Resolution	5- 1
5.3 Large Spinning Sidebands	5- 1
5.4 Poor Signal Reproducibility (Horizontal Drift)	5- 2
5.4.1 Single Direction Drift	5- 2
5.4.2 Random Drift	5- 2

CHAPTER 6 PERFORMANCE TESTS

6.1 Resolution Test	6- 1
6.2 Stability Test (Internal Lock)	6- 4
6.3 Sensitivity Test	6- 7
6.4 Integration Test	6-10
6.5 Spin Decoupling Test	6-13

SERVICE NOTES

CHAPTER 1 GENERAL DESCRIPTION

1.1 OUTLINE

The JNM-PMX60SI is a compact high resolution ^1H NMR instrument that even a beginner will find easy to operate.

- o Especially designed for measurements of hydrogen nuclei, both the field and the resolution of this instrument are very stable. This is achieved through the adoption of a magnetically shielded permanent magnet whose temperature is precision-controlled.
- o Due to the adoption of a large-sized X-Y recorder and chart paper that is calibrated into Hz and ppm, spectrum analysis can be easily carried out.
- o Almost all operations are performed by feather-touch key or push button. Since an automatic mode (interlocking mechanism) is provided, even a beginner can easily carry out measurements.
- o With the incorporation of all solid-state circuitry, this instrument provides both high reliability and high stability. Moreover, maintenance is simple.
- o The instrument has low power consumption. Only 150 W is required when no attachments are used. 400 W is required with all the attachments in use.
- o Spin decoupling and NMR internal lock are built in.

1.2 SPECIFICATIONS

1.2.1 INSTRUMENT SPECIFICATIONS

1. Applicable nuclei	^1H
2. Standard frequency	60 MHz
3. Standard field	1.4092 T (14,092 G)
4. Resolution	0.4 Hz
5. Sensitivity	$S/\text{Nav} \geq 40$ (Ratio between maximum peak (S) of 1% ethylbenzene quartet and average noise (Nav).)
6. Integration accuracy	2% or less (Average deviation when integration waveform of 5% ethylbenzene is successively repeated five times and recorded.)
7. Sweep	
Method	Magnetic field sweep
Width	60, 120, 180, 240, 300, 360, 480, 600, 1200 Hz
Time	25, 50, 100, 250, 500, 1000 sec
8. Magnetic field shift	
Offset	$\pm 60, \pm 600$ Hz
Rough shift	4 kHz (1 kHz steps)
Continuous shift	1.5 kHz (Continuously variable)
9. Sample tube dia. (Outer)	5 mm ϕ
10. Sample volume	0.3 ml (Approx.)
11. Internal lock	
Reproducibility	0.2 Hz (Average deviation when successively repeated 5 times and recorded at sweep time of 250 sec.)
Sweep method	Magnetic field/frequency sweep
12. Spin decoupling	
Frequency range	2.2 kHz
Offset	Rough adjustment: -600, 0, 600 Hz Fine adjustment: 1,000 Hz, continuously variable

Frequency stability	0.2 Hz
Irradiation method	Magnetic field/frequency synchro-track
13. Recorder	Plane-surface type X-Y recorder
Chart size	280 mm (vert.) × 420 mm (horiz.)
14. Power	AC 100/115 V, single phase 50/60 Hz, 400 VA
15. Dimensions/weight	

Console	Dimensions (cm)			Weight (kg)
	Width	Depth	Height	
Spectrometer inc. stand	80	73	97	64
Magnet console	84	69	87	350

1.2.2 INSTALLATION REQUIREMENTS

The following requirements must be met before the PMX60SI can be installed.

1. Room area

Floor space: 2 m × 2.5 m

Room entrance: 80 cm wide or over




Fig. 1-1 Installation layout

- | | |
|----------------------------|---|
| 2. Floor strength | 350 kg/m ² or over, no vibration |
| 3. Temperature | 18°C ~ 28°C
Permissible fluctuation: ±1°C/hr,
±3°C/24hrs. |
| 4. External magnetic field | Install the PMX60SI in an area where magnetic field fluctuation is minimum. Do not install it near a large transformer, an elevator, a magnetic field generator, or a magnetic substance. If ferromagnetic objects are in motion nearby, poor reproducibility will result. |
| 5. Power | AC 100 V ±10%, 50/60 Hz, 500 VA
<ul style="list-style-type: none">• Use an automatic voltage regulator in areas where line voltage fluctuation is large.• Be sure to ground the green wire of the power cable.• The power outlet should be within 3 m of the instrument.• Do not use the same outlet as an air conditioner. |

Note: These specifications may be changed without notice so as to reflect modifications.

1.3 COMPOSITION

1.3.1 STANDARD COMPOSITION




Fig. 1-2 Standard layout

Item	Pcs.
1. Spectrometer (Inc. stand)	1
2. Magnet console	1
3. Air compressor with 3 filters	1
4. Cables/hoses	1 set
5. Accessories (May be changed without notice.)	
① Chart paper (Standard)	100 sheets
(White paper)	100 sheets
② Sample tubes	10
③ Sample tube plugs	10
④ Standard sample (5% mixed sample)	1

- | | |
|---|--------|
| ⑤ Probe cleaning rods (For bobbin and for stator) | 1 each |
| ⑥ Tools | 1 set |
| ⑦ Recorder accessories | 1 set |
| ⑧ Air filters (5 μ , 10 μ , 60 μ) | 1 each |
| ⑨ Spare parts (fuses, lamps) | 1 set |
| ⑩ Vinyl cover | 1 |
| ⑪ Instruction manual | 1 vol. |




Fig. 1-3 Accessories

1.3.2 OPTIONS


1. Oscilloscope
2. Cartridge pen set

1.4 STANDARD SAMPLES AND CONSUMABLES

1.4.1 STANDARD PMX60SI SAMPLES

Sample	Color	Application	Escaping gas	Standard signal	Solvent
Mixed sample *	Sky blue	- o General adjustment - o Resolution - o Sweep width correction	×	TMS	CDCl_3

* Spectrum of mixed sample



Corresponding signal No. (See above)	Sample	Position where signal appears (Hz) ppm
①	TMS	0
②	Cyclohexane	85.6 1,43
③	Acetone	129 2,15
④	Trichloroethane	164 2,73
⑤	P-Dioxane	221 3,68
⑥	Dichloromethane	318 5,3
⑦	Chloroform	439 7,32

1.4.2 SPECIAL ORDER CONSUMABLES AND STANDARD SAMPLES

Method of ordering: Always indicate part No. and nomenclature.

• Consumable Items

Item	Standard designation	Part No.
Chart paper	PMX60SI	421000368
Recorder pen	HP disposable type Red: No. 5081-1190 Blue: No. 5081-1191	
Sample tube	5 mm ϕ	422000400
Spin rotor	PMX60SI	
Recorder pen	K & E cartridge pen	
Recorder ink for cartridge pen	Black	

• Standard Samples

Sample	Color	Application	Escaping gas	Standard signal	Solvent	Part No.
90% ethylbenzene		Ordinary adjustment	○	TMS	CDCl_3	427000467
5% ethylbenzene		Ordinary adjustment, resolution, integration	○	TMS	CDCl_3	427000521
1% ethylbenzene		Sensitivity		TMS	CDCl_3	427000491
90% crotonaldehyde		Spin decoupling	○	TMS	CDCl_3	427000548
20% acetoaldehyde		Resolution	○	TMS	CDCl_3	427000572
95% chloroform		Spinning sidebands	○	TMS	CDCl_3	427000653

1.5 CONNECTING UP OSCILLOSCOPE (OPTIONAL)




Fig. 1-4

Pr. No. 15084

JEOL

CERTIFICATION

Instrument: HIGH RESOLUTION
NMR SPECTROMETER

Model: JNM - PMX 60 Si

Serial No.: NM 105206-619

This is to certify that the subject instrument has passed the final performance inspection at our factory.

Date of Inspection: 20, JUL. 1987

Inspected by:

M. IWASA

Approved by:

K. OHTAKA

Manager K. OHTAKA



Inspection Department Final Test Division

 JEOL LTD. / 日本電子

Tokyo, Japan

CERTIFICAT D'INSTALLATION

DE **JEOL** LTD
New Tokyo Bldg.
3-2, Marunouchi, Chiyoda-ku, Tokyo, JAPAN

Installation et réglage de l'appareil : NMR SPECTROMETER

Modèle : JNM-PMX 60Si N° d'Identification : NM105206-619

A compléter et signer par l'utilisateur.

Date : 29 septembre 1987

Signature : [Signature]

Nom : LE GUILLANTON

Qualité : Directeur de Recherche au CNRS

Organisation : CNRS / UIC

Adresse : BP 808

(N°

Rue)

49005 ANGERS Cedex

(Département

Ville)

Téléphone : 41 88 33 12 Poste : 450

Telex : _____

Signature de l'Ingénieur JEOL :

[Signature]

La signature de ce certificat signifie le début de la période
de garantie de l'appareil de base et des accessoires y figurant

Accessoires	N°	Accessoires	N°