
INSTRUCTION MANUAL

MANUEL D'INSTRUCTIONS

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CONTENTS

- I - GENERAL
- II - USE OF MTDC WITH SPOOLS FOR GLASS TEST TUBES
 - 1/ Standard spool ref :A21.400 for 200 glass test tubes Ø 13 x 100 mm (75 or 65 mm)
 - 2/ Initial set-up
 - 3/ Operation
- III - USE OF MTDC WITH SPOOLS FOR GILSON SNAP TUBES.
 - 1/ Spools
 - 2/ Initial set-up
 - 3/ Operation
- IV - USE OF MTDC WITH SPOOL FOR PREPARATIVE CHROMATOGRAPHY
 - 1/ Spool
 - 2/ Initial Set-up
 - 3/ Operation
- V - USE OF MTDC WITH SPOOL ref.A21.460 for 120 scintillation mini vials Ø 16 x 55 mm.
- VI - USE OF MTDC WITH SPOOL ref.A21.480 for 27 scintillation vials Ø 28x65 mm.
- VII - USE OF THE MTDC/V WITH A GILSON MINIPULS VOLUMETRIC PUMP.
- VIII - ELECTRONIC CIRCUIT DIAGRAMS

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I - GENERAL

This fraction collector enables the collecting of column effluents during analytical or preparative chromatography by either drop counting or timing. It is also possible to collect precise volume by using the collector in conjunction with a Gilson volumetric peristaltic pump.

Characteristics :

- Capacity : 200 glass test tubes Ø 13 x 100 mm (75 or 65mm)
 - or : 200 polypropylene snap tubes of 15 or 6 ml.
 - or : 49 collecting funnels for preparative chromatography
 - or : 120 scintillation mini-vials Ø 16 x 55 mm
 - or : 27 scintillation vials Ø 28 x 65 mm.
- Collection :
 - Timing from 6 sec to 99 min. in increments of 6 seconds.
 - Drop counting from 1 to 999 drops.
- Integrated circuit electronic drop counter with high counting rate (max. 10 per second.)
- Marker plug delivering a contact closure at each tube change.
- A plug for remote control by contact closure
- A plug for electro-valve
- Auxiliary plug (supplying mains power) which is shut off at the end of collection.
- Preselection of the number of fractions to be collected (0 to 200) by means of a movable stopper.
- Automatic shut-off at the end of the collection with diversion of the effluent to a drain tube.
- Size : 35 X 33 X 33 cm
- Operating temperature : 0 to 50°C.

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Accessories supplied with the instrument

- 1 spool.ref.:A21.400 for 200 glass test tubes Ø 13x100 mm.
- 2x100 test tubes ref.: 450.051 Ø 13 x 100 mm.
- 1 mast ref.:B20.259 (1m X Ø 16 mm)
- 1 female 4 pin socket ref.: 506.064
- 1 bag ref.:A21.060 containing :
 - 2 column support rods
 - 1 movable stopper
 - 1 catheter holder
 - 1 teflon tube Ø 1.5x500 mm
 - 1 key for securing the mast

II - OPERATION OF MTDC WITH SPOOL FOR GLASS TEST TUBES

Before mounting the spool, the instrument must be switched off. To do this, the switch (Fig.I.1) must be flicked downwards. Then check that the mains voltage corresponds with that needed for the instrument (220 V - 50 Hz, 115 V - 60 Hz, 100 V - 50/60 Hz) and switch on the power. The circuit breaker (Fig. I,3) must be on.

1/ Standard spool ref.:A21.400 for 200 glass test tubes Ø 13 x 100 mm.

This standard spool can hold 200 glass test tubes Ø 13 x 100 mm but can also be supplied to hold tubes of Ø 13 x 75 mm and Ø 13 x 65 mm (to be specified at the time of the order).

Whatever the type of spool used, the operation procedure is the same.

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By means of the small bar (Fig. II.1) move the test tube supporting base into its raised position (Fig. II). Load the spool with tubes and then move the base back into its lower position (Fig. III).

b) Mounting the spool

Lift up the small lock (Fig.I,10) with the index finger while at the same time pulling out the spool holder handle (Fig.I.11). Mount the loaded spool, and then push in the handle, (Fig.I.11) as far as it will go.

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2/ Initial Set-up

a) Connection between MTDC and recorder

The fraction collector delivers a contact closure at each tube change which enables one to locate the fractions collected on the recorded curve of the absorption spectra during continuous chromatography column effluent analysis. Connect the marker plug on the fraction collector (Fig. IV.1) to the one on the recorder.

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b) Connection between chromatography column and fraction collector.

The outlet tube (catheter) of the column must be connected to the detection block of the head of the fraction collector.

Slightly loosen the screw shown in Fig. V.1). The clamp (Fig.V.2) which fixes the catheter in place can now be moved aside allowing the passage of the catheter. Lower the catheter into the top third of the glass protector tube using the small notch as a guide, then tighten up the screw so as to keep the catheter in place.

Warning ! the catheter must never cut across the detection axis (Fig.V,3) situated half-way up the glass protector tube.

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b) Method of operation.

Position the head over the first tube taking care that the pin (Fig.VI.1) which moves the spool around is exactly above the second hole of the top disc of the spool.

Select the required collection mode (Fig.I.12)

-switch UP : timing

-switch DOWN : drop counting.

Set the number of drops to be collected or the number of time units (1 unit = 6 secs.) by means of the thumb wheel switch (Fig.I.14).

NOTE : with the adjuster set at 0 the collector moves continuously from tube to tube.

Flick the ON and OFF switch (Fig.I.1) up, and check that the counter (Fig.I.15) is at zero. If not, return to zero by pressing the small switch (Fig.I.13). The collector is now ready for operation.

c) Automatic stop at the end of collection.

When collection is finished, the contactor (Fig.VII.1) touches the fixed stopper (Fig.VII.2) and this stops the collector. The luminous counter (Fig. I.15) will go out, but the power is effectively still on. On the other hand, the auxiliary plug will be shut off.

Attention ! At the end, the head stops moving before it has filled the last 2 tubes. Therefore only 198 tubes are filled.

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When the collector has stopped, the detection head opens (Fig.VIII) and the column effluent flows away down the drain channel (Fig. VIII.1). Consequently always make sure the drain pipe (Fig.VIII.2) is appropriately placed to dispose of the liquid.

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d) Selection of number of fractions to be collected.

If one does not want to use the full capacity of the collector, it is possible to pre-select the number of tubes to be filled. To do this, simply replace one of the test tubes with the movable stopper (Fig.I.6).

Attention ! the collector will still miss the last 2 tubes (example: to fill 50 tubes, place the stopper in the place of the 53rd tube).

As before, when the collector stops, the detection head will open and the column effluent will flow away down the drain pipe. The luminous counter will go out, and the auxiliary plug will be shut off.

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e) Removing spool at end of collection

- . flick the ON and OFF switch (Fig.I.1) down.
- . CLOSE the detection head
- . Remove the spool.

III USE OF MTDC WITH SPOOLS FOR GILSON SNAP TUBES.

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c) Mounting the snap tubes on spools ref.:A21.160 andA21.161

. 'Assembling the snap tubes

The boxes of snap tubes contain bands of 50 pre-assembled snap tubes. Every 10th snap tube is pink ; this enables the tubes to be counted easily.

To assemble the snap tubes, see (Fig. XI). To disassemble, follow the reverse procedure.

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Mounting snap tubes on spools ref.:A21.160 andA21.161 see (Fig. XII).

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The snap tubes must be correctly wound around the spool, and pulled as tightly as possible by using the fixing spring (Fig.XII).

d) Mounting the spool

- Spool ref.:A21.160 for 200 15 ml snap tubes. The procedure is the same as in paragraph II,1,b.

- Spool ref.:A21.161 for 200 6 ml snap tubes. In this case an adapter is fitted under the spool see (Fig. XIII).

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- 2/ Initial Set up
(see paragraph II.2).

- 3/ Operation
(see paragraph II.3)

IV USE OF MTDC WITH SPOOL FOR PREPARATIVE CHROMATOGRAPHY.

- 1/ Spool ref.:A21.500 for preparative chromatography.

This spool enables collection of 49 fractions of any volume. But, in this case the collector can only be used in timing mode.

This spool comprises :

- One fixed lower disc which holds the funnels (Fig.XIV.1) . The liquid flows from these funnels into the receiving bottles through tubing.

- One movable upper disc (Fig. XIV.2) which enables distribution of the column effluent into each funnel successively.

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Check that the ON/OFF switch (Fig.I.1) is in the "stop" position. Lift up the small lock (Fig.I.10) with the index finger while at the same time pulling out the spool holder handle (Fig.I.11). Position the lower disc which is fixed to the spool holder handle by means of a small metal foot (Fig. XIV.3). The movable upper disc can now be placed inside the lower disc.

Turn the upper disc so that the catheter holder (Fig. XIV.4) is on your right, and then push in the spool handle (Fig. I.11) as far as it will go.

2/ Initial Set-up

- a) Connection of MTDC to recorder (see paragraph II.2.a)
- b) Connection of chromatography column to fraction collector.

Join the chromatography column outlet tube onto the catheter holder (Fig. XIV.4).

3/ Operation of MTDC

Rotate the upper disc so that the catheter holder comes to rest exactly above the first funnel which is behind the detection head (Fig. XIV). Position

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a) Mise en place des fioles

voir Fig. XV.

the detection head in such a way that the pin (Fig.XIV.5) which moves the spool around is above the first hole of upper disc.

Flick the collection mode switch (Fig. I.12) up(for timing) and set the required number of time units (1 unit = 6 secs) on the thumb wheel switch (Fig.I.14)

Note : With the adjuster set at zero, the collector moves continuously from tube to tube.

Flick the ON and OFF switch (Fig.I.1) up, and check that the counter display (Fig.I.15) is at zero. If not, return it to zero by pressing the small switch (Fig. I.13). The collector is now ready for operation and will stop automatically when collection has finished, (i.e. when the last funnel has been filled). When it has stopped the detection head will open (Fig.VIII) and the luminous counter (Fig.I,15) will go out, but the power is effectively still on.On the other hand,the auxiliary plug will be shut off.

To begin a new collection, flick the ON and OFF switch to STOP, close the detection head and proceed as before.

V USE OF MTDC WITH SPOOL REF.:A21.460
FOR 120 SCINTILLATION MINI VIALS
Ø 16x55 mm.

See paragraph II. Simply replace the glass test tubes with scintillation mini vials.

VI USE OF MTDC WITH SPOOL REF.A21.480
FOR 27 SCINTILLATION VIALS Ø 28x65mm

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b) Mounting the spool.

This spool requires the use of an adapter (Fig. XVI) to position the detection head directly above the scintillation vials. Make sure the ON and OFF switch (Fig.I.1) is in position STOP, and then lift up the small lock (Fig.I.10) while pulling out the spool holder handle (Fig.I.11). Mount the spool in such a way that the fixed stop (Fig. XV.1) is on your right, then push in the spool handle as far as it will go.

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2) Initial Set-up

- a) Connection of MTDC to recorder

see paragraph II.2.a.
- b) Connection from chromatography column to fraction collector.

see paragraph II.2.b.

3/ Operation.

Rotate the spool so that the fixed stop (Fig. XV.1) is behind the detection head and that the pin which moves the spool around is directly above the first hole in the upper disc.

Select the required collection mode (Fig.I.12) switch UP - timing switch DOWN - drop counting.

Set either the number of drops to be collected or the number of time units (1 unit = 6 secs) on the thumb wheel switch (Fig.I.14)

Note : with the switch set at zero, the collector moves continuously from tube to tube.

Flick the ON and OFF switch (Fig.I.1) up and check that the counter display (Fig. I.15) is at zero. If not return it to zero by pressing the small switch (Fig. I.13).

The collector is now ready for operation and will stop automatically when collection has finished. When it has stopped the luminous counter (Fig.I.15) will go out, but the power is effectively still on. On the other hand, the auxiliary plug will be shut off.

The detection head will open (Fig. VIII) and the column effluent flows away down the drain channel (Fig. VIII.1). Consequently, always make sure the drain pipe (Fig. VIII.2) is appropriately placed to dispose of the liquid.

4/ Removal of spool

- . Flick the ON and OFF switch (Fig.I.1) to "stop".
- . Close the detection head.
- . Remove the spool.

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
VIII SCHEMAS ELECTRIQUES.

VII USE OF THE MTDC/V WITH A GILSON
MINIPULS VOLUMETRIC PUMP.

Unscrew the drop counting head plug of the MTDC/V (Fig. IV.3).
Screw in (Fig.IV.3) the plug of the pump-collector connecting cable (supplied with the pump) and screw again in (Fig.IV.2) the drop counting head plug of the MTDC/V.

The switch (Fig.I.12) of the MTDC/V must be switched down : drop counting.

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oxygraphes, respiromètres.

" GILSON " products

Instruments for liquid sample handling :

adjustable volume pipettes, fixed volume pipettes, disposable polypropylene tips, manual and automatic diluters, manual and automatic dispensers, peristaltic pumps, high pressure metering pump.

Systems for automatic handling of liquid samples :

sample changers, system for automatic preparation of samples.

Instruments for liquid chromatography :

fraction collectors, mini cold room volumetric peristaltic pumps, high pressure metering pump, column effluent analysers, recorders, peak detector, automatic gradient former, peak collector, chromatography columns, chemically inert fittings and valves.

Electrophoresis :

gel fractionators, fraction collector bases.

Instruments for respirometry :

oxygraphs, respirometers.