UCO - 0662 - CH  
NOT 022  
Notice de Spectrophotomé-  
tre UV Jobin-Yvon.  
PATSTEC 2395  
Uco. Angers. IBE 17.01.71

I N S T R U M E N T S   S . A .  
=+=+=+=+=+=+=+=+=+=+=+=+=+=+=

DIVISION JOBIN YVON  
-----

DEPARTEMENT CHROMATOGRAPHIE  
-----

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SPECTROPHOTOMETRE UV HOLOCHROME HMD - P  
=====

## I - GENERALITES

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## I - GENERAL

The Holochrome HMD is a dual beam spectrophotometer which has been designed in particular for the qualitative and quantitative monitoring of column effluents. The Holochrome HMD uses a stabilized deuterium light source, and a holographic grating monochromator to provide continuous calibrated wave length adjustment between 190 and 600 nanometers.

### Accessories supplied with the instrument :

- 1 set of 2 analytical blown quartz cells (ref.: 809.511) for differential measurements, 10 mm light path, 70 µl volume, max pressure : 1 Bar (supplied with holder and 0.5 mm I.D. connection tubes).
- 1 cable ref.: 9.533 for 1 channel Gilson recorder or 2 cables for 2 channel Gilson recorders.
- 1 accessory package ref.: 9.532 which contains :
  - . 1 standard recorder cable
  - . 2 preassembled teflon connecting tubes (30 cm long)
  - . 2 polypropylene tube - end fittings (natural colour) for Teflon tubes, O.D. 1.5 mm.
  - . 2 black polypropylene tube end fittings for Teflon tubes O.D. 1.5 mm.
  - . 1 female luer adaptor
  - . 1 male luer adaptor
  - . 2 standard air reference aperture blocks (diameter 2 mm)
  - . 1 event connector
  - . 1 set of 5 couplings.

### SPECIFICATIONS :

Optical unit : Holographic monochromator.

Wave-length : 190 to 600 nm ( $\pm$  2.5nm)

Half amplitude bandwidth : 10 nm.

Base line drift : less than 0.0005 A.U per hour at constant temperature, after warm-up.

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Recorder outputs : 10 mV - 100 mV. A third, fixed output of 100 mV for recording at 2 A.U. sensitivity is also available.

Event marker : allows for event marking on the recorder either by manual command (front panel of instrument) or by remote control (back panel).

#### Standard Voltages :

240 V - 50 Hz ref.: 9.504  
220 V - 50 Hz ref.: 9.500  
115 V - 60 Hz ref.: 9.505  
100 V - 50/60 Hz ref 9.506

Operating temperature : 0 to 50°C  
(see Gilson cold room, Mini RE/2)

## II INITIAL SET-UP

Place the instrument as far as possible from any light source (eg. lamps, windows etc..), and turn the sensitivity control to zero (Fig. I.1). Open the door giving access to the cells (Fig. I.2) and check that these have been correctly mounted.

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1) For working on differential mode with respect to the buffer, (as in the case of gradients or buffer sensitive to U.V. radiations) two cells are used : reference and measurement.

To make sure they are well in place, slide each cell block into its support, and move it fully to the end of the slot (Fig.II,1).

For measuring on differential mode see Fig.III.

2) For working on semi-differential mode with respect to the light source (as in the case of usual separations in which the buffer does not absorb U.V. light), insert a reference aperture (air) block (Fig. IV) instead of the reference cell block. This is not a true differential mode recording but it compensates for any variations of the light source.


FIG. II



FIG. III

BLOC A ORIFICE CALIBRE  
APERTURE (AIR) BLOCK

FIG. IV

Select the required wavelength by means of the direct reading potentiometer (Fig.I.3). This must then be locked using the small locking lever which is underneath the potentiometer, (Fig.I.4).

Do not forget to push the filter (Fig.II.2) in, if you are working in the visible (360 to 600 nm), or to pull it out if you are working in U.V. (190 to 360 nm).

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Close the door of the cell compartment taking care that the tubes from the cells are well positioned in their corresponding slots at the top and bottom of the door (Fig. II.3).

Be sure that the line voltage corresponds with that needed for the instrument 240 V - 50 Hz or 220 V - 50 Hz or 115 V - 60 Hz or 100 V - 50/60 Hz. Switch the instrument on. The circuit breaker which is at the back must be in the "ON" position (Fig. V.1). Press the "Power" switch (Fig. I.5) on and allow the instrument to warm up to ensure stable operation.

### III CONNECTION BETWEEN HOLOCHROME AND GILSON ONE CHANNEL AND TWO CHANNEL RECORDERS.

The Holochrome allows simultaneous recording at two different sensitivity levels.

- at 2 A.U. full scale, in all cases fixed output 100 mV (Fig. V.2).
- between 0.01 and 2 A.U. full scale, according to the sensitivity chosen. Output 100 mV (Fig. VI.3) or 10 mV (Fig. IV.4).

Connect the recorder cables onto the appropriate output sockets of the Holochrome.

### IV OPERATION

When one is working at very high sensitivity certain precautions are necessary to ensure a drift of no more than 0.0005 A.U/h

- constant ambient temperature
- avoid areas where there is smoke or other U.V. absorbing vapours.
- avoid the presence of any intense light source and direct sunlight.
- avoid air draught.

To balance the Holochrome the cells must be filled with the buffer. To

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
FIG. VII

avoid having trapped air bubbles in the cells, the flow should be from the bottom (lower slot tubing) to the top (upper slot tubing).

#### 1) Setting the baseline

With the sensitivity control at zero, set the pen on your recorder for the desired baseline using the pen position control on the recorder.

#### 2) Setting the Holochrome

Turn the sensitivity control of the Holochrome to 2 A.U. and adjust the coarse and fine balance knob of the Holochrome (Fig.I.6) to return pen to baseline. Repeat this operation for the whole range of sensitivities up to 0.01 A.U. and then choose the desired sensitivity.

### V CALIBRATION

All calibration is carried out at the factory. One can check this calibration using a solution with known absorbance, and correct it if necessary. To do this use the small calibration trim pot (Fig. VI.1) which is accessible through a small hole on the left side of the Holochrome, to adjust the recorded value to the real value of the calibration solution.

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### 1) Assembly

Slide the scanner assembly front plate under the Holochrome, and screw the two knurled thumbscrews into the two holes to join the Holochrome and gel scanner (Fig. VIII.1).

. remove the two black plastic plugs from the cell compartment (Fig.I.7 and Fig.X.1).

. slide the cell-blocks out of position and remove from Holochrome. Insert the reference aperture (air) block in place of the reference cell.

. insert the 245 mm black anodized guide tube into the gel tube block and slide this assembly into the place where the measuring cell block was previously mounted (Fig.IX). The gel

tube guide therefore protrudes from the Holochrome to fit into the gel scanning attachment. Close the door of the cell compartment, and insert the 220 mm black anodized gel tube guide through the opening in the door and into the block.

. Unscrew the plastic knurled knob 1/4 of a turn (Fig.VII.1) and slide it to the position shown in (Fig.VII.2). Connected to a pushing rod, this knob enables the gel tube to be moved into the detector cell.

### 2) Setting-up of the Holochrome

Balancing the signal of the Holochrome should be done on a "pure" part of the gel tube which does not contain any sample material, or this can be accomplished by having a plug of pure gel at one end of the gel tube (see paragraph IV.1.2). When the signal is balanced, choose the desired sensitivity range. Insert the gel tube into the guide as far as it will go. Then, switch on the gel scanning attachment (Fig.VII.3). The gel is pushed through the tube at a speed of 20 mm/min which should correspond to the chart speed of the recorder. At the end the gel stops automatically.

Note: It is necessary to use quartz gel tubes for scanning in the U.V. range.

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With the instrument switched off, open the door to the lamp compartment (Fig.X). Unplug the white plastic connector (Fig. X.2). Remove the two screws (Fig.X.3), and take out the lamp block. DO NOT HANDLE THE NEW DEUTERIUM LAMP DIRECTLY WITH THE FINGERS AS FINGERPRINTS CAN DAMAGE IT. If necessary, clean it with alcohol. Position the lamp in place and plug in the white connector.

Note: The deuterium lamp has a counter which allows one to check the length of its life. The lamp is guaranteed for 500 hours, normal lifetime is about 1000 hours. However, lamps which have lasted less than 500 hours can only be replaced free of charge if they are returned with their counter.

### VIII RECOMMENDATIONS

- Do not turn power off and then on again immediately after. Wait until the lamp has cooled down before switching on again so as to prolong the life of the lamp.
- Use the tubing with the screw-on connectors at the inlet and outlet of the cells. A whole range of Gilson tubing, connectors and valves can be supplied on request (see catalogue Chromatography columns)

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. If the lamp does not light up when switched on, check that the circuit breaker on the back of the instrument is on.

. Always make sure that the cells are clean. (Use methanol, and if possible dry them with a draught of dry air or with dry nitrogen).

. If the recorder baseline drifts or is unstable, check the following points:

- block cells should be pushed into position as far as they will go.
- sample and reference cells should be free from any solid particles
- no bubbles are blocking the optical paths of the cells
- the instrument should not be exposed to excessive temperature changes, direct sunlight or draughts.
- the U.V. lamp should work correctly.

---

## HOLOCHROME SPECIFICATIONS

---


|                           |                                                                                                                                                                                                                                                                                                                                                                                       |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Opto-electronics:         | Holographic Monochromator and Deuterium light source (lamp life expectancy is approx. 1000 hours).                                                                                                                                                                                                                                                                                    |
| Wavelength:               | 190-600 nm, ( $\pm 2.5$ nm). (For recording at 360 nm and above, the broad band cut-on filter must be used).                                                                                                                                                                                                                                                                          |
| Half-amplitude bandwidth: | 10 nm                                                                                                                                                                                                                                                                                                                                                                                 |
| Base Line Drift:          | Less than .0005 OD/hour at constant temperature, after warm-up.                                                                                                                                                                                                                                                                                                                       |
| Sensitivity Range:        | .01, .02, .05, 0.1, 0.2, 0.5, 1.0, and 2.0 O. D. full scale (with 10 mm lightpath flow cells).<br>*1.0 O. D. full scale in visible range. (360 to 600 nm.)                                                                                                                                                                                                                            |
| Balance Range:            | With Reference Aperture (air), background density up to 0.5 O. D. can be balanced out.<br><br>When using buffer solution in the Reference Cell, the buffer O. D. should not exceed 1 O. D. (with buffer in Reference Cell, noise and drift may be greater than when using Reference Aperture (air).)                                                                                  |
| Cuvettes                  | Blown Quartz - 10 mm lightpath 70 $\mu$ l volume, STANDARD (Pressure rated to 35 Kg/cm <sup>2</sup> -500 p. s. i.)<br><br>3 mm lightpath, 32 $\mu$ l <u>volume-OPTIONAL</u> (blown quartz)<br><br>KEL-F, with Quartz Windows, for use with degassed liquids- <u>OPTIONAL</u> 10 mm lightpath 32 $\mu$ l volume or 3 mm lightpath, 10 $\mu$ l volume. (Pressure rated to 500 p. s. i.) |
| Recorder Outputs:         | 10 mV - 100 mV; A third, fixed output is included for recording at 2 OD sensitivity (with 10 mm lightpath flow-cell), independent of range setting. (See Step 5 under Initial Set-up, for recorder hookup information.) (Fixed 2 O. D. sensitivity is for use with 10 mV Recorder only, unless HM was ordered with a special output).                                                 |

---


## MODEL HM - REPLACEMENT PARTS

| <u>PART<br/>NUMBER:</u> | <u>DESCRIPTION:</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|-------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 0585                    | Replacement Deuterium Lamp Assembly, mounted in a block and pre-aligned for optimum efficiency.                                                                                                                                                                                                                                                                                                                                                                                                      |
| 0586                    | HM Kel-F cuvette assembly, 3 mm light path, 10 $\mu$ l volume (for use with degassed liquids). *                                                                                                                                                                                                                                                                                                                                                                                                     |
| 0587                    | HM Kel-F cuvette assembly, 10 mm light path, 32 $\mu$ l volume (for use with degassed liquids). *                                                                                                                                                                                                                                                                                                                                                                                                    |
| 0588                    | Replacement Blown Quartz Cuvette Assembly, 10 mm light path, 70 $\mu$ l volume (standard flow cell). *                                                                                                                                                                                                                                                                                                                                                                                               |
| 05471                   | Event Connector, mating, Switchcraft 09GM4M. (Switch Closure).                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| 05546                   | Reference Aperture (air). (Two required per unit; to allow checking deuterium lamp and HM electronics independent of flow cells).                                                                                                                                                                                                                                                                                                                                                                    |
| 05559                   | Recorder Cable, 1.52 meters long (5-feet).                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| 13287                   | HM Amplifier printed circuit board, complete.                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| 13288                   | HM Lamp/Power Supply printed circuit board, complete.                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| 0528                    | GEL SCANNING ATTACHMENT FOR HM. Intended for scanning 8 mm O.D. quartz (UV range: Glass gel tubes may be used for scanning in the visible range) gel tubes such as those which can be used with the Gilson Gel Fractionators (B-100 or B-200 base unit with GMA/GCB Gel Mincing Assembly). Contains motor driven advancing mechanism which moves the gel through the light path at a rate of 2 cm/min. Includes gel tube block, mounting hardware for Model HM, and 5 quartz gel tubes, 265 mm long. |
| 27267                   | Quartz gel tubes for gel scanning attachment, 165 mm long, 8 mm O.D. X 6 mm I.D. (Pkg. of 10).                                                                                                                                                                                                                                                                                                                                                                                                       |
| 27268                   | Quartz gel tubes for gel scanning attachment, 265 mm long, 8 mm O.D. X 6 mm I.D. (Pkg. of 10).                                                                                                                                                                                                                                                                                                                                                                                                       |

\* If liquid reference is required for sample monitoring (Dual Beam), two cuvette assemblies are needed (one Sample, one Reference). If using solvent that does not absorb UV light, only a Sample flow cell is required, along with one 05546 Reference Aperture (air) positioned in the reference side of the HM.

PC Board #13288

REGULATOR-DEUTERIUM LAMP P.S.

PART No

DATE 8-18-68 OLSON APP

REVISED

MECHANICAL ELECTRONICS