

Spectra-Physics



SP4200 Operators Manual

Part Number A0099-091 (insert only)
Part Number A0099-092 (W/BASIC and Binder)

SP4200
TABLE OF CONTENTS

SECTION 1: INSTALLATION

1.1	UNPACKING.....	1-1
1.2	CHECKOUT PROCEDURE.....	1-2
1.2.1	Electrical Connections.....	1-2
1.2.2	Printer/Plotter Paper Installation.....	1-3
1.2.3	Battery Backup of the SP4200.....	1-4
1.2.4	Internal Adjustments.....	1-4
1.2.5	Test Signal.....	1-4
1.2.6	Connection to a Detector.....	1-5
1.2.7	Input Selector SHUNT Position.....	1-6
1.2.8	External Connections.....	1-7
1.2.9	LABNET Unit Number Selector Switches.....	1-7
1.2.10	Calculating Memory Requirements.....	1-7

SECTION 2: BRIEF OPERATION

2.1	INTRODUCTION.....	2-1
2.2	CHECKLIST FOR OPERATION AND DATA ACQUISITION.....	2-1
2.3	USING OTHER FUNCTION KEYS.....	2-4

SECTION 3: OPERATION

3.0	INTRODUCTION.....	3-1
3.1	KEYBOARD.....	3-2
3.2	FUNCTION KEYS (Figure 3.3).....	3-4
3.2.1	INJECT/END (A or B) Keys.....	3-4
3.2.2	Paper Advance.....	3-4
3.2.3	A <-> B Keys.....	3-5
3.2.4	Soft Reset.....	3-6
3.2.5	STATUS Key.....	3-6
3.2.6	PT EVAL Key.....	3-7
3.2.7	LEVEL Key.....	3-7
3.2.8	RUN TIME Key.....	3-7
3.2.9	TIME Key and Time and Date Settings.....	3-8
3.2.10	Keyboard and File Lock/Unlock.....	3-8

3.3	DISPLAY AND INDICATOR LED'S.....	3-8
3.4	RUNNING A CHROMATOGRAM.....	3-9
3.5	USING DIALOG TO SET UP A FILE OF OPERATING PARAMETERS.....	3-10
3.5.1	File Identification.....	3-16
3.5.2	Replay Option.....	3-16
3.5.3	File Name.....	3-19
3.5.4	Time Functions.....	3-20
3.5.5	Calculation Method Entry.....	3-21
3.5.6	Number of Levels.....	3-21
3.5.7	Injections/Level: INJ/LEV.....	3-22
3.5.8	Component Table.....	3-22
3.5.9	RRT REF PEAK.....	3-23
3.5.10	RF REF PEAK: IP=.....	3-23
3.5.11	Internal Standard.....	3-24
3.5.12	SAMPLE TABLE.....	3-24
3.6	FILE LISTING.....	3-26
3.6.1	File Identification.....	3-26
3.6.2	Report Parameters.....	3-30
3.6.3	Time Functions.....	3-30
3.6.4	Component Table.....	3-31
3.6.5	Sample Table.....	3-32
3.7	EDITING A FILE.....	3-33
3.7.1	Deleting A File.....	3-33
3.7.2	Changing One or More Parameters.....	3-33
3.7.3	Editing Through Dialog.....	3-33
3.8	COPYING FILES.....	3-34
3.9	MONITORING VALUES.....	3-35
3.9.1	To Print a Parameter.....	3-35
3.9.2	To Display a Parameter Value.....	3-35
3.9.3	Printing or Displaying Names.....	3-35
3.11	SELECTING PEAK HEIGHT INSTEAD OF AREA.....	3-36
3.12	PEAK OVERRANGE MESSAGE.....	3-37

SECTION 4: DATA REDUCTION

4.1	GENERAL.....	4-1
4.2	CALCULATION OF COMPONENT CONCENTRATIONS.....	4-2
4.2.1	Response Factors (RF).....	4-3
4.2.2	Response Factor Variances (CV).....	4-3
4.2.3	Method 0 - Area % (or Height %).....	4-4
4.2.4	Method 1 - Normalization Using Response Factors..	4-5

4.2.5	Method 2 - Internal Standard.....	4-7
4.2.6	Method 5 - External Standard.....	4-10
4.3	COMPONENT IDENTIFICATION AND REFERENCE PEAKS.....	4-13
4.3.1	Absolute Retention Time (RT).....	4-13
4.3.2	Reference Peaks and Relative Retention Times (RRT).....	4-15
4.3.3	Multiple Reference Peaks.....	4-16
4.3.4	Dead Time and Relative Retention Times (DT).....	4-16
4.4	MULTILEVEL CALIBRATIONS.....	4-17
4.4.1	Calculations Used for Linear and Non-Linear Fits..	4-19
4.4.2	Method Selection with Multilevel Calibration/Analysis.....	4-21
4.4.3	Dialog for Multilevel Calibration.....	4-24
4.4.4	Parameters Used for Multilevel Calibration.....	4-27
4.5	REPORT OF CONCENTRATION OF UNKNOWNNS.....	4-27
4.6	IDENTIFICATION OF VALUES BEYOND EXPECTED CONCENTRATION RANGES.....	4-29
4.7	STATISTICS REPORTS.....	4-30
4.7.1	Statistics Calculations.....	4-31
4.7.2	Importance of the IX (INDEX) Parameter.....	4-32
4.8	ELIMINATING PEAKS FROM REPORTS (TI, MA and MC).....	4-32
4.8.1	Solvent and Other Early Peaks (TI).....	4-32
4.8.2	Minimum Concentration (MC).....	4-33
4.8.3	Minimum Area (MA).....	4-35
4.8.4	Weighted Update (WU).....	4-35
4.9	NORMALIZATION OF CONCENTRATIONS IF MN = 2 or 5.....	4-36
4.10	REPORT OF COMPONENT TABLE VALUES (SUMMARY REPORT).....	4-36
4.11	PEAK GROUPING.....	4-38
4.12	STORAGE OF COMPONENT DATA.....	4-39
4.12.1	Peaks Identified in Component Table.....	4-39
4.12.2	Peaks in the Chromatogram.....	4-40
4.13	OFFLINE DATA ENTRY.....	4-41
4.14	REPLOTTING COMPONENT FROM PREVIOUS CALIBRATION RUN.....	4-42

SECTION 5: TIME FUNCTIONS

5.1	GENERAL DESCRIPTION.....	5-1
5.2	TIME FUNCTIONS THROUGH DIALOG.....	5-1
5.3	IMMEDIATE EXECUTION OF TIME FUNCTIONS.....	5-2
5.4	EDITING OF FILE TIME FUNCTIONS.....	5-3
5.4.1	Deleting a Time Function.....	5-3
5.5	PLOTTING AND REPORTING TIME FUNCTIONS.....	5-4

5.5.1	Peak Markers (PM).....	5-4
5.5.2	Printing of Retention Times (NT).....	5-4
5.5.3	Printing of Time Function Event Times (NE).....	5-5
5.5.4	Auto Zero (AZ).....	5-6
5.5.5	End Run (ER).....	5-7
5.5.6	Start Run (SR).....	5-7
5.5.7	BCD/Autosampler Vial Decoding (XD).....	5-7
5.6	CONTROL OF EXTERNAL EVENTS THROUGH TTL CONTACT CLOSURES.....	5-11
5.7	CONTROL OF BASELINE CORRECTION.....	5-11
5.7.1	Integration Inhibit (II).....	5-11
5.7.2	Integration Inhibit Mode (IM).....	5-11
5.7.3	Forced Tangent Skimming (FT).....	5-13
5.7.4	Disabled Tailing Peak Logic (TP).....	5-13
5.7.5	Forced Forward Horizontal Baselines (FH).....	5-16
5.7.6	Forced Backward Horizontal Baselines (BH).....	5-17
5.7.7	Negative Peaks (NP).....	5-18
5.7.8	Valley-to-Valley Baselines (BL).....	5-20
5.7.9	Prevention of Baseline (NB).....	5-21
5.7.10	Maximum Baseline Level (MB).....	5-22
5.8	CONTROL OF THE PRINTER/PLOTTER.....	5-22
5.8.1	Attenuation (AT).....	5-22
5.8.2	Control of the Chart Speed (CS).....	5-23
5.9	CONTROL OF GROUPED PEAKS.....	5-23
5.9.1	Maximum Peak Width Multiple (MP).....	5-24
5.9.2	Peak Onset (PO).....	5-24
5.10	CONTROL OF THE "GO" COMMAND (GO).....	5-25

SECTION 6: PRINCIPLES OF INTEGRATION

6.1	GENERAL.....	6-1
6.2	INTEGRATION: DATA SAMPLES AND DATA BUNCHES.....	6-2
6.2.1	Slope Calculation.....	6-3
6.3	PEAK DETECTION.....	6-5
6.3.1	Peak Thresholds.....	6-5
6.3.2	Effective Slope Sensitivity.....	6-7
6.4	AUTOMATIC PARAMETERS: PARALLEL PEAK DETECTORS.....	6-7
6.4.1	Maximum Peak Width Multiple (MP).....	6-8
6.4.2	Changing Number of Peak Detectors.....	6-9
6.5	PEAK CREST IDENTIFICATION.....	6-9
6.6	PEAK ENDINGS.....	6-10
6.6.1	Baseline Detection.....	6-10

6.6.2	Tangent Fit.....	6-11
6.6.3	Fused Peak Splitting.....	6-12
6.6.4	Shoulder Peak Splitting.....	6-12
6.7	TAILING (SKEWED) PEAKS.....	6-14
6.8	BASILINE CORRECTION AND AREA ALLOCATION.....	6-15
6.9	PLATEAU DETECTION.....	6-17
6.10	PEAK HEIGHTS.....	6-17
6.11	PEAK DETECTION AND BASELINE CORRECTION TROUBLESHOOTING.....	6-19
6.11.1	Review of Role of PW and PT in Peak Detection.....	6-19
6.11.2	Selecting a PW Value.....	6-20
6.11.3	How to Determine Your Peak Detection and Baseline Correction.....	6-22
6.11.4	Small Peaks or Noise Being Integrated.....	6-22
6.11.5	Noise Spikes.....	6-23
6.11.6	Optimizing Detector Input.....	6-23
6.11.7	Peaks Not Detected.....	6-24
6.11.8	Detection of Fused and Shoulder Peaks.....	6-25
6.12	WHEN AUTOMATIC PARAMETERS ARE NOT ENOUGH.....	6-26
6.12.1	Automatic Parameters- Changing PW Before Run Beginning.....	6-26
6.12.2	Resetting PW and PT During Run (Special Dialog)...	6-26
6.12.3	Automatically Doubling the PW and PT Values (T2)...	6-27
6.12.4	Changing PW and PT via Time Functions.....	6-28
6.13	DATA SLICE MODE.....	6-29
6.13.1	Enabling and Disabling Data Slice Mode.....	6-29

SECTION 7: PERIPHERAL DEVICES

7.1	DESCRIPTION.....	7-1
7.2	EXTERNAL STORAGE.....	7-2
7.2.1	Current Loop/Cassette Player Cable.....	7-2
7.2.2	Recording and Playing Back Raw Data.....	7-3
7.2.3	Saving and Loading Files and BASIC Programs.....	7-3
7.2.4	Storing Data on the Kerr Minifile.....	7-4
7-3	EXTERNAL CONTROLS CABLE.....	7-4
7-4	LABNET CONNECTOR.....	7-5
7-5	AUTOSAMPLERS.....	7-6
7.5.1	Hewlett-Packard 7670-7671.....	7-6
7.5.2	Varian 8000.....	7-7
7.5.3	Micromeritics Model 725.....	7-8
7.5.4	WISP Autosampler.....	7-8

TABLE OF ILLUSTRATIONS

SECTION 1 INSTALLATION

1.1	SP4200 Showing Power Switch and Channel Selector Switches.....	1-1
1.2	Location of the SP4200 Serial Tag.....	1-2
1.3	Example of Chromatogram and Report Using the Test Signal.....	1-5
1.4	Detector Input Connections to the SP4200.....	1-6

SECTION 2 BRIEF OPERATION

2.1	The SP4200 Computing Integrator.....	2-1
2.2	Test Chromatogram and Method 0 Report.....	2-3
2.3	Example of a STATUS Printout.....	2-5

SECTION 3 OPERATION

3.1	SP4200 Controls and Indicators.....	3-1
3.2	SP4200 Input Switches, Channel A and B.....	3-2
3.3	SP4200 Alphanumeric Keyboard and Function Keys.....	3-3
3.4	STATUS Printout.....	3-6
3.5	Display and Indicator LEDs.....	3-9
3.6	Example of a Chromatogram and Report Using Default Parameters.....	3-11
3.7	Examples of Baseline Correction (BC) Codes.....	3-11
3.8	Example of a Replay Dialog Printout.....	3-14
3.9	Example of Replay Option Responses.....	3-16
3.10	Example of Response to File Identification Query.....	3-20
3.11	Example of a Time Function Table.....	3-20
3.12	Example of Response to the Method Number: MN= Query.....	3-21
3.13	Example of Component Table Dialog.....	3-22
3.14	Example of Sample Table Dialog.....	3-24
3.15	Example of File Listing	3-27
3.16	Example of the Report Parameter Listing.....	3-30
3.17	Example of a Time Function Table Listing.....	3-31
3.18	Example of Component Table Listing.....	3-31
3.19	Example of Edited Component Table Listing.....	3-32
3.20	Example of a Sample Table Listing.....	3-33
3.21	Editing the Time Function File Through Dialog.....	3-34
3.22	Example of a Report with both Peak Areas and Heights.....	3-36
3.23	Report Showing Peak Overrange Condition.....	3-37

SECTION 4 DATA REDUCTION

4.1	Example of a Chromatogram and Report Using Method 0.....	4-1
4.2	Example of a Chromatogram and Report Using Method 1.....	4-6
4.3	Example of a Chromatogram and Report Using Method 2.....	4-9
4.4	Example of a Chromatogram and Report Using Method 5.....	4-12
4.5	The "Window" for Component Peak Identification.....	4-13
4.6	Criteria for Component Peak Identification Windows.....	4-14
4.7	Example of "NEW FILE" Report.....	4-15
4.8	The Largest Peak in a Reference Peak Window.....	4-16
4.9	Example of a Multilevel Calibration Requiring a Linear Fit....	4-17
4.10	Example of a Multilevel Calibration Requiring a Nonlinear Fit..	4-17
4.11	Graph of a Multilevel Calibration Using a Linear Function....	4-20
4.12	Example of a Report of the Coefficients from a Nonlinear Calibration.....	4-22
4.13	Example of a Multilevel Calibration Using a Nonlinear Function.....	4-22
4.14	Example of Dialog for Multilevel Calibration of One Component and an Internal Standard.....	4-23
4.15	Example Graph for nitrogen multilevel calibration.....	4-25
4.16	Graph showing the relationship between VA,VB, VC, etc., and Concentration for a single component.....	4-26
4.17	Response Factor of Default Peak Applied to all Unidentified Peaks.....	4-28
4.18	Asterisk Identifies Components Whose Concentrations Exceed Range of Expected Concentrations.....	4-29
4.19	Example of a Statistics Report.....	4-30
4.20	Example of a BAD DATA Printout.....	4-32
4.21	Examples of Reports with and without T1.....	4-33
4.22	Examples of Reports with and without MC.....	4-34
4.23	Example of a Report with an MA being Used.....	4-35
4.24	Example of a Report with Unnormalized Total Reported.....	4-36
4.25	Example of Report with LC=3.....	4-37
4.26	An Example of Peak Grouping.....	4-38
4.27	Example Offline Data Entry and Report.....	4-41

SECTION 5 USE OF TIME FUNCTIONS IN REPORTS AND BASELINE CORRECTION

5.1	Example of a Time Function Dialog.....	5-2
5.2	Example of a Time Function Listing.....	5-2
5.3	How to Interpret the Location of a Peak Marker.....	5-4
5.4	Examples of Chromatograms with and without Time Function Times Printed.....	5-5
5.5	Example of Use of AZ to Automatically Rezero Baseline.....	5-6

5.6	Example of the End Run Time Function.....	5-7
5.7	Location of Connector J110.....	5-9
5.8	Effects of II and IM on Baseline Correction.....	5-12
5.9	Forcing Tangential Baseline Correction Using FT.....	5-14
5.10	Disabling Tailing Peak Identification with TP.....	5-15
5.11	Forward Horizontal Baseline Correction with FH.....	5-16
5.12	Backward Horizontal Baseline Correction with BH.....	5-17
5.13	Examples of Ways to Deal with Negative Peaks.....	5-19
5.14	Force Baseline at Valley Point.....	5-20
5.15	NB Time Function used to Prevent Finding Baseline Between Peaks.....	5-21

SECTION 6 PRINCIPLES OF INTEGRATION

6.1	Graphic Example of Segmented Data Sample.....	6-2
6.2	Equation for Calculating Slope Using Data Bunch Integrals....	6-3
6.3	First Derivative of Slope Calculation Plotted.....	6-4
6.4	Effect of Data Bunching on Peak Detection Sensitivity.....	6-4
6.5	Peak is Detected when Slope Exceeds Peak Threshold Value.....	6-5
6.6	Step in Peak Detection and Confirmation.....	6-6
6.7	Automatic Parameters Using Multiple Peak Detectors.....	6-8
6.8	Example of Peak Crest Identification.....	6-9
6.9	Examples of Using Retention Time Location to Determine Appropriateness of PW Parameter Value.....	6-10
6.10	Baseline Test.....	6-11
6.11	Criteria for Peaks Ending with Tangent Fit.....	6-11
6.12	Fused Peaks Splitting.....	6-12
6.13	Splitting of Front Shoulders.....	6-13
6.14	Splitting of Rear Shoulders.....	6-13
6.15	Tailing Peak Test.....	6-14
6.16	Nonhorizontal Baseline Correction.....	6-15
6.17	Forward Horizontal and Backward Horizontal Baseline Correction.....	6-15
6.18	Tangential Baseline Correction.....	6-16
6.19	Riders on Tailing Peaks are Baseline Corrected Using Tangential Baselines.....	6-16
6.20	The Plateau Test.....	6-17
6.21	Peak Height is Taken as Value of Largest Bunch.....	6-19
6.22	Best PW Value for a Peak.....	6-20
6.23	Using Retention Time Placement to Determine Appropriateness of Integration Bunch Size.....	6-21
6.24	Using Peak Markers to Identify Peak Beginnings and Endings...	6-21
6.25	Examples of Baseline Corrections and Codes.....	6-22
6.26	Example of Noise Spike Rejection using PW Value.....	6-23

6.27	Example of Small Peak Detection using PW and PT Values.....	6-24
6.28	Placement of Peak Markers at Valley Points.....	6-25
6.29	Changing PW Value to Change Peak Ending.....	6-25
6.30	PW and PT Parameter Doubling using T2.....	6-27
6.31	Examples of Dialog and Listing of Time Functions for PW and PT Parameters.....	6-28

SECTION 7

7.3	LABNET Cable Connection.....	7-5
7.1	Current Loop/Cassette Cable Connection.....	7-2
7.2	External Controls Cable Connector.....	7-5
7.4	Cable 2135-010 Diagram.....	7-7
7.5	Modified Cable 2135-010 Diagram.....	7-7

APPENDIX A

A.1	Location of the Line Power Connector.....	A-1
A.2	Position of the Line Voltage Selector for 120 Volts.....	A-1
A.3	Location of Power Switch and Input Selector Switch.....	A-2
A.4	Location of SP4200 Cover Mounting Screws.....	A-3
A.5	Analog Board Voltage Range Selector Switch.....	A-3
A.6	LABNET Interface Board and Switch Sl.....	A-4

TABLE OF TABLES

SECTION 1 INSTALLATION AND CHECKOUT

1.1	Accessory Parts List.....	1-2
-----	---------------------------	-----

SECTION 3 OPERATION

3.1	Description of Function Keys.....	3-5
3.2	Description of SHIFT and Function Keys.....	3-6
3.3	Time Function Codes and Their Meanings.....	3-14
3.4	File Variables Used in Reports.....	3-24

SECTION 4 DATA REDUCTION

4.1	Guidelines for Selecting a Calculation Method.....	4-2
4.2	Parameters Used in Multilevel Calibrations.....	4-27
4.3	Parameters Associated with Components in the Component Table.....	4-39
4.4	Parameters Used to Access Storage.....	4-40

SECTION 5 USE OF TIME FUNCTIONS

5.1	Time Function Codes and their Meanings.....	5-3.1
5.2	J110 T-Function Pin-Out.....	5-10
5.3	Non-Time Function Pin-Outs of J110.....	5-10

APPENDIX A

Line Voltage Fuse Ratings.....	A-2
Memory Availability vs System Configuration.....	A-5

APPENDIX B

Table of ASCII Characters.....	B-1
SP4200 Keyboard Function Keys.....	B-5
SP4200/SP7100/LABNET Fault Codes.....	B-6

APPENDIX C

Nuances of SP4200 LABNET.....	C-1
-------------------------------	-----

FEDERAL COMMUNICATIONS COMMISSION RADIO FREQUENCY INTERFERENCE STATEMENT

Warning: This equipment generates, uses and can radiate radio frequency energy and if not installed and used in accordance with the instruction manual, may cause interference to radio communications. As temporarily permitted by regulation it has not been tested for compliance with the limits for Class A computing devices pursuant to Subpart J of Part 15 of FCC Rules, which are designed to provide reasonable protection against such interference. Operation of this equipment in a residential area is likely to cause interference in which case the user at his own expense will be required to take whatever measures may be required to correct the interference.

A<->B Key3-5
 Absolute Retention Time4-13
 AC (Statistics Enable)4-30
 Active File3-11, 3-13
 Additional Reports4-36
 AN query3-24
 Area
 Allocation.....6-15
 Baseline Corrected5-11
 Percent Method4-4
 Percent Report4-4
 Relative Peak %4-9
 Sum of
 Attenuation (AT)4-22
 Auto ZERO (AZ).....5-6
 Automatic Parameters6-7, 6-26
 Autosampler
 use with Vial Decoding5-7
 Connections.....7-6 to 7-9

 Baseline
 Area Allocation6-15
 Correction Codes5-11
 Correction Control.....5-11
 Correction Troubleshooting .6-19
 Detection6-10
 Forced Backward Horizontal .5-17
 Forced Forward Horizontal ..5-16
 Maximum level (MB).....5-22
 Noise6-22
 Prevention of (NB).....5-21
 Subtraction, Replay Option..3-18
 Setting level3-7
 Valley-to-Valley (BL).....5-20
 Battery Backup Test.....A-2
 BASIC Programs, Saving.....7-3
 BCD/Autosampler Vial Decoding.....5-7
 BYTES IN BASIC message3-6
 BYTES IN STORAGE message3-14

 Calculating Memory Requirements....A-5
 Calculation Method Entry3-18
 Calculations
 Linear/Non-Linear Fits4-19
 Method 0, Area % or Height...4-4
 Method 1, Normalization4-5
 Method 2, Internal Standard .4-7
 Method 5, External Standard 4-10
 Slope6-3
 Statistics4-31
 CALIB Key3-5
 Calibration
 Flag (CA)3-22
 Index (CI)3-22
 Multilevel ...3-19, 4-17 to 4-27

Calibration, Con't
 Normal3-19
 Samples Between3-24
 Channel
 Assignment1-1
 Chart (See Printer/Plotter)
 Speed5-23
 CI (Samples Between Calib.)3-24
 Codes
 Component Table4-39
 FaultB-6
 Report3-28, 3-29
 Sample Table3-29
 Time Function3-15, 5-3
 Component Concentration (CC)3-23
 Component Identification4-13
 Component Table3-22, 3-31, 3-32
 Codes4-39
 Component Window (CW)3-30
 Concentration
 For Calibrations (CC)3-23
 Units (CU)3-28
 CONTROL key3-3
 CONTROL key Functions.....3-3 to 3-5, 3-8
 Copy Files3-30
 Current Loop7-2

 Data
 Bunches6-2
 Integration (see Integration)
 Samples6-2
 Slice Method6-29
 Date/Time Setting3-8
 Dead Time (DT)4-16
 DELETE Key3-2
 Deleting
 Characters3-2
 Files3-33
 Parameters3-35
 Time Functions5-3
 Detector
 Changing number of Peak6-9
 Connection to SP42001-5
 Input1-5
 Level AdjustmentA-4
 Noise1-5
 Optimizing Input6-23
 Output1-5
 DIALOG Key3-11
 Dialog
 Editing3-33
 Exiting3-25
 Step-by-Step3-16 to 3-32
 Disable Tailing Peak Logic (TP)5-13

Edit File (FE)	3-33, 5-3	Index Parameter (IX)	4-32
Editing		INJ Key (A or B)	3-5
Component Table	3-31	Inject/End Keys	3-4
Dialog	3-33	Injection/Level query (RC).....	3-22
File	3-33	Integration	
Time Functions	5-3	Data Bunches	6-2
Electrical Connections	1-2, A-1	Data Samples	6-2
End Run Time Function (ER)	5-7	Inhibit (II)	5-11
ENTER Key	3-3	Inhibit Mode (IM)	5-11
ERASE Key	3-2	Time (see Run Time)	
ESCAPE Key	3-4	Internal Standard	
Exiting Dialog	3-4	Amount (IS)	3-24
Expected Concentration Ranges	3-23, 4-29	Method (MN=2)	3-21
External Controls Option	5-9	Peak (IP)	3-23, 3-28
External Standard	4-10	Internal Switches.....	A-3
External Storage.....	7-2 to 7-4	INT INH Key	3-5, 5-11
Fault Codes	Appendix 1-6	KA - KC	4-27
File		Keyboard	
Active File (FI)	3-12, 3-16	Function Key	3-4
Copying	3-34	Operation	3-2
Deleting	3-33	Overview	3-1
Edit File (FE)	3-33		
Editing	3-33, 5-3	Last Concentration Switch (LC)	4-36
Identification	3-16, 3-26	LABNET	
Listing	3-26	Connector	7-5
Name	3-19	Interface Board Switch	A-4
Fits		Unit Number Setup	1-7
Least Squares	4-19	Least Squares Fits	4-19
Linear/Non-Linear	4-19	LEDs (lights)	
Tangent	6-11	Display	3-8, 3-9
Forced Tangent Skimming (FT)	5-13	Indicator	3-9
Form Feed Sensor	B-2	RUN	3-9
Formatting Reports	4-36	LEVEL Key	3-7
Free Memory Message	3-6	Linear/Non-Linear	
Function Keys	3-4	Fits	4-19
Fused		Query	4-18
Peak Splitting	6-12	Listing (see Printing)	
Peaks	6-25	Lock/Unlock	
 		File (LK)	3-8
GO command	5-25	Keyboard	3-8
Graphics Control	4-24, 4-42	LS (Multilevel Calibration Switch) ..	4-17
Grouped Peaks		Lower Case	3-3
Assigning Names to	4-38		
Baseline Correction Code	5-11	Maximum	
Using Time Function File	4-38	Baseline Level	5-22
Control of	5-23	Peak Width Multiple (MP).....	5-24, 6-8
 		Messages	
Horizontal Baselines		Fault	Appendix 1-6
Forced Forward (FH).....	5-16	Peak Overrange	3-37
Forced Backward (BH).....	5-17	Method	
		Area % or Height	4-4

Method con't		Eliminating From Reports	4-32
Data Slice	6-29	Endings	6-10
External Standard	4-10	Fused	6-12
Internal Standard	4-7	Grouping	4-38, 5-23
Normalization	4-5	Identified in Component Table	4-39
Number Query (MN)	3-21	Indicator LED	3-9
Selection W/Multilevel Calib..	4-21	In Chromatogram	4-40
Minimum		Internal Standard (IP)	3-24
Area (MA)	4-35	Markers (PM).....	5-4
Concentration (MC)	4-33	Maximum Peak Width (MP)	5-24
Multilevel Calibration.....	4-17	Not Detected	6-24
Method Selection	4-21	Number	
Parameters used	4-27	Onset	5-24
Multiple Reference Peaks	4-16	Overrange Message	3-28
		Parallel	6-7
Name of File Query (NM)	3-26	Reference (IP)	3-21, 4-15
Negative Peaks (NP).....	5-18	Relative Ret. Time Reference ..	3-23
Noise		Retention Time, Absolute (RT) ..	4-13
Integration of	6-22	Shoulder	6-12
Spikes	6-23	Small	6-22
Normal Calibration	3-21	Solvent (TI).....	4-32
Normalization		Storage File (PS)	4-39, 4-40
Method	4-5	Tailing (skewed)	6-14
of Concentrations	4-36	Threshold (PT)	3-5, 6-5, 6-19
Number of Levels			6-26
NV=0	3-22, 4-18	Width (PW)	3-30, 6-7, 6-8, 6-10,
NV=1	3-22, 4-18		6-19, 6-20, 6-26 to 6-28
NV 1	3-22, 4-18	Peak Height	
Query	3-22, 4-18	Reporting	6-17
		Selecting (PH)	3-36, 6-17
Offline Data Entry	4-41	Peak Splitting	
Offset - chart (OF)	5-6	Fused	6-12, 6-24
Optimizing Detector Input	6-23	Shoulder	6-12, 6-24
Options		Plateau Detection	6-17
External Control	5-9	Prevention of Baseline (NB)	5-21
Replay RX+	3-16	Printer/Plotter	
Parameters (also see Codes)		Attenuation	5-22
Changing	3-18, 3-33	Offset	5-6
Printing	3-35	Paper Installation	1-3
for Multilevel Calibration ...	4-27	Control	5-22
Report	3-30	Speed	5-23
Display Value	3-35	Printing	
Index (IX)	4-32	Names	3-35
Peak		Parameters	3-35
Area	3-36	Retention Times	5-4
Crest Identification	6-9	System Status	3-6
Detection	6-5, 6-19	Time Function Times (NE)	5-5
Detection Indicator (see RUN light)		PT EVAL Key	3-7
Detection Threshold (PT)	6-19, 6-26		
		Recording Raw Data	7-3
		Reference Peak Window (RW)	3-29
		Reference Peaks (RF)	3-23

Relative Retention Time	3-23	SHIFT key	3-2
Reference Peak Index (RP)	3-23	Shoulder Peak Splitting	6-12
Reference Peak Window (RW)	3-29	Slope	
Replicate		Calculation	6-3
Analysis (RA)	3-24	Negative	6-4
Calibrations (RC)	3-29	Positive	6-4
Replotting	4-42	Sensitivity	6-7
REPORT Key		Small Peaks	6-22
Report		Soft Reset	3-6, 2-14
Codes	3-28, 3-29	Solvent Peaks (TI)	4-32
Component Value Tables	4-36	Standard Deviation (Statistics).....	6-14
Concentrations of Unknowns	4-27	Start Run (SR)	5-7
Formatting	4-36	Statistics	
Parameters in File Listing	3-30	Calculations	4-31
Using Area Slice Method	6-29	Report	4-30
Using Area% Method	4-4	Selecting of	4-30
Using External Standard Method.....	4-10	STATUS Key	3-6
Using Internal Standard Method.....	4-7	Status listing	3-6
Using Normalization Method	4-5	Storage	
Response Factor (RF)	4-3	of Component Data	4-39
Response Factor Variances (CV)	4-3	of Peak Data	4-39, 4-40
Retention Time Absolute (RT)	4-13	Summary Reports	4-36
Run			
Begin/End	5-7	Tailing Peak Logic (disable) (TP) ..	5-13
Control Keys	3-7	Tangent	
LED (light)	3-9	Fit	6-11
Number (RN)	3-29	Skimming	5-13
Time	3-7	Time	
RUN TIME Key	3-7	Absolute Retention	4-13
RX+ Replay Option	3-16	Base (TB)	3-7
Graphics Control	4-24, 4-42	Relative Retention	3-21, 4-13
Menu	3-16	Run	3-7
Play Chromatogram	3-17	TIME key	3-8
Record Chromatogram	3-17	Time Related Keys	3-8
Single Chan Baseline Subtract.....	3-18	Time Functions (TF)	
Storage Bytes	3-17	Adding	3-20
Subtract Stored Chromatogram.....	3-17	Codes	3-15
		Deleting	5-3
Sa - Sg	5-7	Dialog Mode	5-1
Sample		Direct Mode	3-15
Amount (SA)	3-25	Editing	3-20
Between Calibration (CI)	3-24	Plotting And Reporting	5-4
Index (SI)	3-25	Times (TT)	3-20, 3-30, 5-1
Injection/	3-24	Values (TV)	3-20, 3-35, 5-1
Name (SN)	3-25		
Report (see Report)		Unknown Concentration (DP)	4-27
Table	3-24, 3-32		
Table Codes	3-29	Values, Monitoring	3-35
SAMPLE Key	3-24		
Saving Files	7-3	Weighted Update (WU)	4-35
Scale Factor	3-25		

INDEX OF FILE DIALOG AND TIME FUNCTION PARAMETERS

AC 3-28, 4-30, 31, 32, 39
AN 3-24, 28
AT 1-2, 3-9, 15, 5-22
AZ 3-15, 5-6
BC 1-3, 3-11, 6-22
BH 3-15, 5-17, 6-15
BL 3-15, 5-20
CA 3-27, 28
CC 3-23, 27, 4-39
CH 1-3, 3-5, 27
CI 3-14, 24, 27, 28, 4-3
CM 3-27
CN 3-14, 22, 23, 27
CS 1-2, 3-9, 15, 5-22, 23
CU 3-24, 25
CV 3-28, 4-3, 4
CW 3-28, 30, 4-13, 16
CZ 3-27, 28
DT 3-28, 4-16
DP 3-28
EC 3-14, 23, 28, 4-29, 39
EH 4-29
EL 4-29
ER 3-15, 5-7
ET 3-15
FE 3-14, 16, 20
FH 3-15, 5-16, 6-15
FI 3-14, 16, 26
FILE NAME 3-19, 3-20
FT 3-15, 5-13, 14
GO 3-15, 5-25
GP 3-15, 4-38, 5-23
II 3-15, 5-11, 12
INJ/LEV 3-22
IM 3-15, 5-11, 12, 13
IP 3-14, 23, 24, 28, 30
IS 3-14, 3-25
IX 3-28, 4-3, 32
KA 4-19, 21, 27, 39
KB 4-19, 21, 27, 39
KC 4-19, 21, 22, 27, 39
LC 3-28, 4-36, 37, 39
LK 3-28
LS 3-28 4-17, 18
MA 4-32, 35, 40, 6-22
MB 3-15, 5-22
MC 3-29, 4-32, 33, 34, 6-22
MP 3-15, 5-24, 6-8, 9
MN 3-16, 21, 26, 27, 20, 29, 4-4, 36
NB 3-15, 5-21
NE 3-15, 5-5
NM 3-26

NO 3-29, 4-36
 NP 3-15, 5-18, 19
 NT 3-15, 5-4
 NV 3-14, 21, 22, 29, 4-3, 18, 24, 27
 OD 3-29
 PH 3-27, 31, 4-4, 6-17
 PM 3-15, 5-4, 6-21
 PO 3-15, 5-24, 24
 PS 3-27, 4-40
 PSA 4-40
 PSF 4-40
 PSH 4-40
 PSR 4-40
 PST 4-40
 PT 3-15, 27, 30, 6-5, 6, 7, 8, 17, 19, 20, 22, 24, 25, 26, 27, 28
 PW 1-2, 3-15, 27, 30, 6-7, 8, 10, 18, 19, 20, 21, 23, 24, 25, 26, 27, 28
 RA 3-14, 24, 29, 4-32
 RC 3-14, 22, 29, 4-3, 24, 27, 32
 RF 3-22, 23, 4-3, 4, 39
 RN 3-29
 RP 3-14, 23, 27, 29,
 RRT 3-23, 4-3, 15
 RT 3-14, 21, 22, 23, 27, 4-3, 13, 15, 24, 29, 39
 RW 3-27, 29, 32, 4-15, 16
 SA 3-24, 25
 Sa 4-27
 Sb 4-27
 Sc 4-27
 Sd 4-27
 Sd 4-27
 Se 4-27
 Sf 4-27
 Sg 4-27
 SC 4-31
 SI 3-14, 24, 25, 27, 29
 SM 3-27, 33
 SN 3-24, 25, 27, 33
 SP 3-27, 29
 SR 3-15, 5-7
 SZ 3-27, 29
 TB 3-27, 5-1
 TF 3-14, 20, 4-38, 5-1, 3
 TFN 1-2, 5-2, 3, 3.1
 TP 3-15, 5-13, 15
 TT 3-14, 20, 30, 4-38, 5-1, 2, 3
 TV 3-14, 20, 30, 4-38, 5-1, 3
 T1-T8 3-15, 29, 4-32, 5-9, 10, 18, 6-26, 27,
 UC 4-39
 VA 4-24, 25, 27
 VB 4-24, 25 27,
 VC 4-25
 WU 3-29, 4-35
 XD 3-15, 5-7, 8
 XF 3-14, 24, 25, 27

SECTION 1 INSTALLATION

1.1	UNPACKING.....	1-1
1.2	CHECKOUT PROCEDURE.....	1-2
1.2.1	Electrical Connections.....	1-2
1.2.2	Printer/Plotter Paper Installation.....	1-3
1.2.3	Battery Backup of the SP4200.....	1-4
1.2.4	Internal Adjustments.....	1-4
1.2.5	Test Signal.....	1-4
1.2.6	Connection to a Detector.....	1-5
1.2.7	Input Selector SHUNT Position.....	1-6
1.2.8	External Connections.....	1-7
1.2.9	LABNET Unit Number Selector Switches.....	1-7
1.2.10	Calculating Memory Requirements.....	1-7

SECTION 1 INSTALLATION AND CHECKOUT



Figure 1.1 SP4200 Showing Power Switch and
Channel Selector Switches

1.1 UNPACKING

1. Upon receipt of the instrument inspect the packaging for evidence of damage. If the packaging appears to be damaged, have the damage noted by the carrier on both the delivery receipt and the freight bill. The carrier is responsible for damage incurred during shipping.
2. Unpack the instrument and inspect it for damage. If damage is found notify the carrier and your Spectra-Physics representative.

A list of accessories and their part numbers is included on the next page.

NOTE

It is the buyer's responsibility to note any damage at the time of delivery, and report any damage to the carrier and your Spectra-Physics representative immediately. Do not return goods to Spectra-Physics without authorization from a Spectra-Physics representative.

Table 1.1 Accessory Parts List.

P/N	QUANTITIES		DESCRIPTION
	Ver. 030/070 (120V)	Ver. 040/080 (240V)	
A0099-092	1	1	Instruction Manual SP4200
A0022-010	1	1	Z-Fold Paper
A0099-093	1	1	Instruction Card
A1366-010	1	0	Power Cord, UL tested
6003-0330	0	1	Power Cord, CEE Coded
5401-0620	1	1	Alignment tool
5101-0060	2	0	Fuse, SloBlo, 2.0 amp.
5101-0550	2	4	Fuse, SloBlo, 1.0 amp.
5101-0660	6	6	Fuse, SloBlo, 5.0 amp.
5101-0740	1	1	Fuse, SloBlo, 6.25 amp.
A1249-010	1	1	Hood, Connector, Mod.
2114-0380	1	1	Hood, Connector
2105-2110	1	1	Connector, 18/36 pin
2105-0910	1	1	Connector, 6/12 pin
A2116-010	1	1	Cover, Connector, External
A2156-010	1	1	Installation Option Cables Connection

1.2 CHECKOUT PROCEDURE

Locate the SP4200 near the detector. Make sure the placement provides access to the rear of the instrument.

1.2.1 Electrical Connections

The SP4200 operates on AC line frequencies of either 50 or 60 Hz, and is preset at the factory for either 120 VAC or 240 VAC as described on the Serial Tag (Figure 1.2). If the AC line voltages in your area do not come within the range of +10 to -20% of these ratings, change the line voltage input selector to accommodate your line voltage (see Appendix 1 for instructions on changing the line voltage selection). Line voltage inputs of 100, 120, 220, or 240 can be selected.





Figure 1.2 Location of the SP4200 Serial Tag

1.2.2 Printer/Plotter Paper Installation


A pack of thermographic paper is included in the accessory kit. If this pack is the only one you have, you should order more immediately (P/N A0022-010).

To Install Paper

1. If there is a sheet still in the instrument, hold the ENTER key down to let it run out.
2. Remove the paper you are installing from its plastic jacket.
3. Lift the top page; make sure the number "001" is on the top right corner. Fold page 200 back over page 001 so page 199 is on the bottom.
4. Push back the clear plastic writing tray to expose the paper bay entrance.




5. Place the paper pack in the paper bay, and slide the paper down and back. A slight up and down movement may be necessary to make the paper slide.




6. Place the top of the first page in the paper receiver slot. Press the ENTER key to start the paper roller.

7. Hold the ENTER key down and guide the paper into the slot until the sprocket keyholes in the paper are engaged by the platen sprockets and the top of the page is driven up to about 7 cm above the top of the print head.



> d



8. Close the plastic writing tray and feed the page into the slot in the tray. If the top of the page is not parallel with the slot, pull the upper edge of the lower corner while feeding paper with the ENTER key until it is parallel.

1.2.3 Battery Backup of the SP4200

The SP4200 is equipped with a battery backup feature to prevent loss of data stored in read/write memory (RAM) in the event of a power failure. Appendix A describes a test of the battery backup system.

1.2.4 Internal Adjustments

It is assumed that internal controls and switches are properly set at the factory and no modification of these settings is necessary for initial operation of the SP4200. If changes become necessary, refer to Appendix A for directions on adjusting the analog board voltage range selector switch and the data system analog board adjust control.

NOTE: The Analog board voltage range selector switch is factory preset at 1.0 V, which is suitable for most applications.

1.2.5 Test Signal

Turn on the SP4200 by pressing the power switch shown in Figure 1.1. Set both channel A and B (if present) input selector switches to the TEST position. Press INJECT/END A. A "peak" should be plotted at approximately 0.33 minute intervals at 60 Hz and at 0.4 minute intervals at 50 Hz line power. After five or six peaks, press INJECT/END again. A report of retention times and areas for the peaks should be given, with the time between peaks and their areas consistent. Figure 1.3 shows an example of the test signal and its report. Be sure to try both channels A and B (if present). Enter PLOT"B" AUTO to plot the channel B signal.




Figure 1.3 Example of Chromatogram and Report Using the Test Signal

2.2.6 Connection to a Detector

It is usually advisable to use the "recorder" outputs from the detector rather than the "computer" or "integrator" outputs. The latter provide an unattenuated signal to the SP4200 which may include unnecessary electronic noise or detection of small peaks that are not of interest.

NOTE

Do not change the detector attenuator or range settings during integration or between integrations that are to be compared. Doing this changes the signal input to the SP4200. The SP4200 plot can be attenuated if it is necessary to keep peaks on scale. (Section 3.3 discusses printer/plotter attenuation.)

The banana jacks are color coded for the interconnect cable P/N A1463-010. Red is the positive lead, black is the signal return, and the clear lead is the shield. Try connecting the detector both with the shield grounded through the A2156-010 assembly, and then NOT grounded at all. Use the connection that is the least noisy.

Figure 1.4 below shows the detector input banana jacks. There is an "Installation Option Cables Connection" assembly (P/N A2156-010) in the Accessory Kit. The assembly is for use in a high electronic noise environment. See **NOTES**, below. Push the assembly onto the grounding bracket that is just to the right of the input banana jacks. Connect the grounding shield of the cable to the ground of the Installation Option. This ground is a Philips head screw near the top outer edge.



Figure 1.4 Detector Input Connections (banana jacks) to the SP4200

To test the connections:

1. Set the input selector switch to the NORM position and press the INJECT/END key that corresponds to the channel of interest.
2. Rotate the detector "balance" or "zero" control to provide what would normally be a positive pen deflection. If the printer/plotter print head is deflected in a negative direction, reverse the input lines or change the detector polarity switch, if one is available.
3. Press the LEVEL key. When the signal input to the SP4200 is zeroed, the display of LEVEL reads 1000 ± 2 . One unit in the display is equivalent to 10 microvolts, thus the properly adjusted SP4200 has an offset of 10 millivolts. See Section 6 for a complete description of the necessity for this baseline offset.

NOTES

Use of P/N A2156-010 may be required in some installations to avoid EMI interference in other communications equipment.

A signal of less than -10 mV causes a condition where no integration can take place. When this occurs, there is no display of level and plotting stops. The display is blank.

1.2.7 Input Selector SHUNT Position

Place the input selector switches in the SHUNT position. Press the LEVEL key. The display should show 1000 ± 2 . If the display does not show 1000 ± 2 , refer to Appendix A.3.2 for adjustment instructions. Press the A $\leftrightarrow$ B key and repeat the above action for the other channel, if present.

1.2.8 External Connections

External connections and cables, including LABNET connections, current loop/cassette player cable, and the external controls cable, are covered in Section 7 of this manual, "Peripheral Devices and Storage."

1.2.9 LABNET Unit Number Selector Switches (Range 0 to 15)

The LABNET unit number for either channel is selected on the LABNET thumbwheel selector switches located below the right side of the keyboard. The left hand switch is channel A. The right hand switch is channel B. If LABNET is NOT going to be used, set both switches to zero; or if only one channel is going to be used with LABNET, set the other channel to zero.

NOTE

Setting only one channel to zero disables that channel and the unit becomes, in effect, a single channel integrator. This results in about 800 bytes of additional memory for use by the active channel.

The actual unit numbers selected can be shown in the LED display. First, set the desired unit number or numbers on the thumbwheels. They are readout in the LED display BUT have not been executed. Now, execute a "soft reset" by pressing the CTRL and SHIFT keys simultaneously, holding them down and pressing RPRT. You also may execute this by typing RESET.

NOTE

A soft reset stops integration on BOTH channels, stops the execution of BASIC Programs on both channels, sets Time Base (TB) to zero, and much more. Appendix E-20 further describes the uses of soft reset.

Switch S1 (see Appendix A-4) on the LABNET interface board selects the range of the LABNET unit number selectors. This is preset at the factory for unit numbers less than 16. To set LABNET unit numbers in the range of 16 to 31, refer to Appendix C, "Nuances of LABNET" in this manual.

1.2.10 Calculating Memory Requirements

The SP4200 provides 16K bytes of RAM memory. Of this, several thousand bytes are required for overhead and reserve, depending on system configuration. If you are writing large or complex programs, you may wish to maximize the amount of memory available to you. This can be accomplished by configuring a single channel, not connected to LABNET. The procedure is covered in Appendix A.