
POLAROGRAPHIE

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type

PRG 3

Three-electrode
DC
and
superimposed AC
POLAROGRAPHY
SYSTEM
for research
and analysis
applications

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This instrument is specially intended for three-electrode polarography, with or without superimposed AC signal ; its particularly broad capabilities allow it to be used in research as well as in analytical fields. In many cases, the results obtained are equivalent to those of pulse polarography ; in addition, new possibilities are opened such as studies of adsorption and electrode kinetics.

The equipment is composed of a rack and a shielded measuring cell. In the rack are mounted :

- A PRT 30-01 potentiostat ;
- An UAP 3 adapting unit ;
- An EPL 2 potentiometric recorder.

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1/ Potentiostat

The PRT 30-01 potentiostat has been specially designed for polarography with superimposed AC, tensammetry and cyclic voltammetry. New circuits, entirely fitted with silicon semiconductors (transistors, diodes, linear integrated circuits) produce a unique combination of high gain (10 000), wide frequency response (5 MHz unity gain cross-over point), negligible drift ($\pm 200 \mu\text{V}$ long term) and particularly low background noise (less than $100 \mu\text{V p.p.}$) ; the potentiostat phase characteristic is adjustable with precision by means of a resistance/capacitance network. The instrument can be controlled by an externally programmed DC signal, with or without superimposition of an AC sinewave. The output voltage ($\pm 30 \text{ V}$) allows measurements to be made in relatively high resistance media (non aqueous solvents) ; the available output current ($\pm 100 \text{ mA}$) is suitable both for potentiostatic and potentiokinetic techniques, and for electrolysis processes with the use of a diaphragm (supporting electrolyte purification, coulometry, etc.).

It is possible to use other potentiostats of the PRT series, but their use will cause a slight reduction in the overall performance of the instrument.

2/ Recorder

The EPL 2 potentiometric recorder is a high precision instrument, which can be used either as an integral part of the PRG 3 polarography system, or independently. The preamplifier circuits are housed in a plug-in module, which can be selected among nearly twenty types, according to the electrical parameter to be recorded ; for polarographic purposes, the plug-in to be used may be a "V" type with a 12,5 mV range (TV 11 GD, TV 1 GR, TV 1 G) or a "I" type such as TI 20 G or TI 1 G. The chart is driven by a pulse operated stepping motor : an electronic gear box allows twelve speeds to be selected, from 0.15 to 600 mm/minute, in forward or in backward direction, with the possibility for the motor to be controlled by an external signal, thus allowing X-Y operation.

Detailed characteristics of the plug-in modules can be found in leaflet 6-TEPL, those of the recorder in leaflet 8-EPL 2.

3/ Adapting unit

The UAP 3 unit includes the following elements :

● DC voltage programming circuit

Two ten-turn potentiometers allow the presetting of the voltage excursion limits, which can be located, at will, in the anodic or in the cathodic zone ; the scanning is effected by a precision potentiometer, mechanically linked to the recorder chart drum shaft by means of an electromechanical clutch and a gear train. The voltage span, on the chart, can be of 20, 50 or 100 mV/cm.

The following excursion modes can be selected by a switch :

- Single excursion, with automatic stop ;
- Double excursion, with automatic stop ;
- Multiple-cycle excursion.

In the two latter cases, recorder acts as an X-Y tracer.

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
● AC voltage source

The superimposed AC signal is generated by a high stability phase-shift oscillator, which also produces the adjustable phase reference signal applied to the demodulating network (patent pending). The frequency is continuously

adjustable from 3 Hz to 3 kHz in three ranges, covered by a ganged potentiometer with precision dial ; the amplitude is adjustable between 0 and 500 mV in nine ranges, covered by coarse and fine potentiometers ; the exact value is monitored by a built-in meter.

● Phase demodulation

The demodulation phase angle is adjustable from - 45 to + 135 ° (- 50 to + 150 grades) in three ranges, covered by a ten-turn potentiometer with digital display (1000 reading points). This allows both in-phase and out-of-phase detection of the AC response, thereby eliminating capacitive current and enhancing the actual faradaic response of the system. Ordinary (non phase-sensitive) demodulation is also possible.

● Drop-time control

The drop time of the mercury electrode is controllable between 0.01 and 100 seconds by means of an electromagnetic hammer, laterally striking the capillary with an adjustable strength. The hammer is driven by an electronic amplifier associated with a frequency divider ; the imposed drop time frequency is always an integral submultiple of the signal frequency, in order to avoid beat phenomena.

● Current amplifiers and sampling circuit

The AC amplifier of the current measuring circuit has high gain and wide frequency response ; it has been specially designed in order to reduce to a minimum undesirable phase shifts. The most sensitive current ranges are obtained by increasing the DC amplifier gain, after demodulation ; twenty ranges are thus obtained, from 5 nA to 5 mA (AC or DC) for the recorder f. s. d. (250 mm = 10 ") ; currents of higher values can be measured by fitting the recorder with a "multirange" plug-in.

● Damping circuit

Damping is obtained by an operational amplifier connected in a low-pass filter configuration ; the time constant is adjustable between approximately 0.1 and 10 seconds, in two ranges adjustable by a potentiometer.

● Complementary circuits

The DC polarography section includes a five range offset current source, each range being covered by a ten-turn potentiometer with digital display (offset adjustable from 0 to ± 10 mA), and a residual current compensation circuit (slope correction) adjustable between 0 and 5 μ A/V, in two ranges covered by a ten-turn calibrated potentiometer, with digital display.

Several pilot lights and a meter with function switch allow the correct operation of the instrument to be checked. Control knobs and pilot lights corresponding to a specific function are identified by the same color.

4/ Measuring stand type CPR 3 B

This stand, specially designed for high sensitivity AC measurements, is fitted with a removable shield in the form of two half-cylinders of perforated steel sheet. An electronic temperature controller (reference RTP-POL) and a pump, which feeds the cell jacket circuit with thermostating fluid, are located in the stand base. The latter also includes a pressure reducing valve, connected to the distributing circuit of the inert gas (N_2) used for removing oxygen from the solution.

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
Thanks to the elimination of capacitive current, the phase-sensitive demodulation allows quantitative measurements to be made at concentrations as low as 10^{-7} M and semi-quantitative ones at 10^{-8} M, in reversible systems ; the sensitivity is of the order of 10^{-6} M for non-reversible systems (organic species).

Due to the fact that the level of superimposed AC signal can be set to a very low value, it is possible to achieve separation of substances whose peak potentials are only 30 to 50 mV apart, depending upon the system reversibility and concentration ratio of the studied species. Proper adjustment of the demodulation phase angle and of the signal frequency makes it possible to enhance or minimize the peaks corresponding to electrochemically active species, according to their kinetic behaviour at the electrode ; the same is true for adsorbable species, giving tensammetric peaks.



Cyclic voltammetry techniques, for which the PRG 3 systems is particularly suitable, allows corrosion and microcorrosion studies to be made (either in DC conditions or with superimposed AC signal) as well as fundamental studies of reversible and non-reversible systems. DC and even AC amperometric titrations are to be particularly noted, among the numerous potential capabilities of this instrument. Moreover, the elimination of capacitive current allows superimposed AC signal measurements in the kilohertz range, which looks promising for non-reversible reactions and non-aqueous solvent studies.

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
Left and below : Curves obtained at various concentrations of electroactive species. Concentrations of Pb^{++} and Cd^{++} are indicated for each curve ; supporting electrolyte : HCl , 10^{-1} M.

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The PRG 3 polarograph is so designed that it can be associated with a pulse polarograph UAP 4 and a switching unit type CM 34. The apparatus thus assembled is suitable for every polarographic technique (pulse ; superimposed AC ; DC) as well as numerous other electrochemical methods, such as programmed or cyclic voltammetry.

Detailed information and leaflet upon request.

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