
Enregistreur XY BD 90/91

KIPP & ZONEN

Enregistreur XY - BD 90/91

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DELFT-HOLLAND

BD 90/91 (8103)

DIRECTIONS FOR USE
(preliminary issue)
BD 90/ BD 91 RECORDER

KIPP & ZONEN

DELFT-HOLLAND

1. General Information.

1.1. Introduction.

The Kipp & Zonen BD 90/91 recorder is meant to visualise graphically the relation of supplied voltages.

This recorder could be given an outstanding physical presentation by combining the latest mechanical and electrical techniques.

A special feature of this design is the revolutionary carriage construction driven by a minimum of moving parts. Plug-in units facilitate exchange without recalibration. Electronic pen travel end stops and automatic power cut-off in case of mechanical overload are standard.

When driven off scale, the pen is lifted, and a flashing led gives an indication.

A carrying handle in front of the plug-in units also allows a horizontal or slant chart table position.

2. Inspection and Installation.


2.1. Unpacking.

The recorder is packed in a specially developed shipping carton and if damage to this carton is evident, ask the carrier's agent to be present when the instrument is unpacked.

Inspect the instrument for mechanical damage and if there is damage, contact your local supplier.

Check before throwing away the carton whether the standard accessories are delivered.

2.2.1.

KIPP & ZONEN
DELFT-HOLLAND

Standard accessories:

- 6 fibre pens short
- 6 fibre pens long (BD 91 only)
- 100 sheets of graph paper DIN A3
- 1 power cord
- 10 fuses
- 1 instruction manual
- 1 dust cover

KIPP & ZONEN

DELFT-HOLLAND

2.2.2.

1. Power switch.
2. Mains connector with fuses.
3. Input terminals X-amplifier.
4. Remote control connector.
5. Input terminals Y-amplifier.
6. Reposition lever carrying handle.
7. Input modules.
- 7.1. Overrange indicator.
- 7.2. Zero suppression switch (option).
- 7.3. Polarity switch.
- 7.4. Variable / fixed range mode switch.
- 7.5. Input / zero volt selector switch.
- 7.6. Zero setting knob.
- 7.7. Range selector.
- 7.8. Variable range knob.
8. Time base module.
- 8.1. Reset / start time sweep.
- 8.2. Units selector.
- 8.3. Sweep factor selector.
- 8.4. X / time base mode selector.
9. Control unit.
- 9.1. Pen up / down switch.
- 9.2. Record / standby switch.
- 9.3. Paper hold / free switch.
- 9.4. Power on led.
10. Chart positioning led's.

KIPP & ZONEN

DELFT-HOLLAND

2.3. Mains voltage.

For protection of operating personnel and equipment, the recorder has factory installed wiring and fuses for different mains voltages.

If the main supply differs from the indication on the recorder, please notify your local distributor.

2.4. Input signals.

The terminal marked "L" (low) must be connected to that side of the source having the lowest impedance to ground, the other side of the source to "H" (high).

For shielding purposes the "Measure Ground" terminal can be used.

However, this terminal is not meant for mains safety grounding. When the remote control unit (option) is fitted see also 3.3.8

KIPP & ZONEN

DELFT-HOLLAND

3. Operating procedure.

3.1. Writing system.

Reliable fibre pens are used as a standard. After having removed the cap, the pen must be gently inserted into the penholder avoiding sideways pressure.

When the recorder is not in use, the tip of the fibre pens should be covered by the cap, avoiding drying of the pens.

3.2. Paper loading.

On the control unit set "standby/record" switch to standby.
Set power switch to "on".

The pen moves to the far right side of the recorder.

Set paper free/hold switch to free.

Just under the writing table three red leds are located which are positioned parallel to the X-writing direction.

Adjust a sheet of graph paper so the leds illuminate a graticule line, and set the paper free/hold switch to hold.

Return switch standby to record.

It also is possible to change paper in record position, however electrostatically the paper is only given free when both pens are in "up" position.

3.3. Measuring.

- 3.3.1. Set on control unit pen up/down switch to up.
Set on plug-in amplifiers zero/range switches to zero.
Set on time base unit (option) mode switch X/T to X.

3.3.2. Polarity.

Set polarity switch to disired direction.

For a positive voltage on terminal red and the polarity switch in position|positive, the direction of the pen will be from zero to the right (X) or far side(Y).

KIPP & ZONEN

DELFT-HOLLAND

3.3.3. Range.

Set variable/fixed switch to fixed.

Set range switch to desired span.

3.3.4. Zero.

Move with zero knob pen to desired origin coordinates.

Set zero/range switch to range.

First make a dry run. When this is correct, recommence and set on control unit pen to down position.

If necessary, readjust zero knob with pen in down position.

3.3.5. Variable mode.

In this mode, ranges between the fixed spans may be set, down to 40% of the fixed spans.

Set variable/fixed switch to variable.

Set zero/range switch to range.

With maximum input amplitude adjust the variable knob to the desired deflection.

3.3.6. Zero suppression (option).

If, with large input voltage swings, the pen runs off scale on the right side (X) or far side (Y), the zero position can be shifted to the left (X) or front (Y) by turning the zero suppression knob counter clockwise

3.3.7. Time base unit (option).

On control unit, set pen to up position.

On the time base unit set X/T switch to T.

Set reset/start switch to reset.

Set sweep factor and unity switch.

Set with zero knob on X unit pen course to starting position.

KIPP & ZONEN

DELFT-HOLLAND

With pen in down position, readjust fine to starting line.

Set reset/start switch to start.

When a second sweep is to be made, press reset switch;
the pen is lifted during reset cycle.

After pen down, press switch to start position.

KIPP & ZONEN

DELFT-HOLLAND

3.3.8. Remote control (option).

For remote control a connector on the rear is added to which several input/output and control lines are connected.

All line voltage conditions with respect to pin 2, unless otherwise stated.

Loading on outputs; 100 kOhm/V unless otherwise stated.

<u>PIN.</u>	<u>Mnem.</u>	<u>Function.</u>	<u>Condition.</u>															
1	Vs	Voltage supply for external use.	+ 24 V unstabilised to pin 2, 50 mA max.															
2	GC	General common.	Common of X amplifier, Control unit and time base unit.															
3	OPY1	Overrange Y1 servo system.	Output: + 15 V to pin 2; normal operation 0 V; Overrange on either side of scale, pen is lifted															
4	Phd	Paper hold down control	On control unit paper hold/free switch in free position. External: switch open or TTL high: paper free. switch closed or TTL low: paper hold.															
5	U/D2	Pen up/down control Y2 BD 91 only.	<table><tr><th><u>Frontswitch</u></th><th><u>U/D</u></th><th><u>Pen</u></th></tr><tr><td>up</td><td>high</td><td>up</td></tr><tr><td>up</td><td>low</td><td>down</td></tr><tr><td>down</td><td>high</td><td>down</td></tr><tr><td>down</td><td>low</td><td>up</td></tr></table>	<u>Frontswitch</u>	<u>U/D</u>	<u>Pen</u>	up	high	up	up	low	down	down	high	down	down	low	up
<u>Frontswitch</u>	<u>U/D</u>	<u>Pen</u>																
up	high	up																
up	low	down																
down	high	down																
down	low	up																
6	U/D1	Pen up/down control Y1.	TTL or external switch to pin 2.															

TTL or external switch to pin 2.

KIPP & ZONEN
DELFT-HOLLAND

<u>in.</u>	<u>Mnem.</u>	<u>Function.</u>	<u>Condition.</u>
	MU	Muting of servosystems X and Y1.	External switch between pin 39 and pin switch open - servomotors on switch closed - servomotors off
	NC	Not connected.	
	NC	Not connected.	
0	NC	Not connected.	
1	RC/SB	Record/standby.	On control unit rec/standby switch in record position. External: switch open or TTL high; record. switch closed or TTL low; standby.
12	NC	Not connected.	
13	NC	Not connected.	
14	NC	Not connected.	
15	MX1	Monitor output preamplifier X.	50 mV/cm, with respect to pin 43 Load 20 kOhm min.
16	NC	Not connected.	
17	NC	Not connected.	
18	MY1	Monitor output preamplifier Y1.	50 mV/cm, with respect to pin 21 Load 20 kOhm min.
19	NC	Not connected.	
20	NC	Not connected.	
21	Y1C	Common of Y1 main amplifier.	
22	NC	Not connected.	
23	NC	Not connected.	
24	NC	Not connected.	

KIPP & ZONEN

DELFT-HOLLAND

<u>Pin.</u>	<u>Mnem.</u>	<u>Function.</u>	<u>Condition.</u>
25	Y2C	Common of Y2 main amplifier only for BD 91.	
26	HX	High input X amplifier.	} BD 91 only. }
27	LX	Low input X amplifier.	
28	HY1	High input Y1 amplifier.	
29	LY1	Low input Y1 amplifier.	
30	HY2	High input Y2 amplifier.	
31	LY2	Low input Y2 amplifier.	
32	OPY2	Overrange Y2 servosystem.	BD 91 only. Output: + 15V to pin 2: normal operation

Wired to preamplifier or
main amplifier (50 mV/cm).

OV: overrange on either side
of scale and pen is
lifted.

33	TCL	Timebase unit common logic.	Is connected to common X amplifier pin 2.
34	A4	} Control lines time base- generator prom.	} Front switch set fully counter clockwise to 20. Input TTL or switch to pin 33.
35	A3		
36	A2		
37	A1		
38	A0		

A4	A3	A2	A1	A0	Sweepfactor
1	1	1	1	1	20
1	1	1	1	0	10
1	1	1	0	1	5
1	1	1	0	0	2
1	1	0	1	1	1
1	1	0	1	0	0.5
1	1	0	0	1	0.2
1	1	0	0	0	0.1
1	0	1	1	1	0.05
1	0	1	1	0	EXT - see pin 40
0	x	x	x	x	Second part of prom programmed upon request.

**KIPP & ZONEN**
DELFT-HOLLAND

<u>Pin.</u>	<u>Mnem.</u>	<u>Function.</u>	<u>Condition.</u>
39	R/S	Reset/start of time base generator.	On time base unit reset/start switch in reset position. External: switch open or TTL high; reset time base. switch closed or TTL low going; start time base.
40	EXT	External generator input.	Frontswitch time base factor on position EXT. or Setting controllines A0-A4. See pins 34 - 38. T.T.L. level 1 pulse ≈ 0.1 mm
41	S/M	Unit setting: Seconds/cm - Minutes/cm.	Front switch time base in position min./cm. T.T.L. low or switch closed - s/cm. T.T.L. high or switch open - min/cm.
42	D/A	Digital/analog converter output of time base unit.	0.5 mV/pulse 12 bit resolution, 10 bit linearity.
43	TCA	Time base unit common analog.	Is connected to pin 33.
44	OXR	Overrange X servosystem righthand side.	- 12 V to pin 2 : no overrange + 12 V to pin 2 : overrange.
45	OXL	Overrange X servosystem lefthand side.	- 12 V to pin 2 : no overrange + 12 V to pin 2 : overrange.
46	OY1F	Overrange Y1 servosystem farside.	- 12 V to pin 21 : no overrange + 12 V to pin 21 : overrange.
47	OY1C	Overrange Y1 servosystem close side.	- 12 V to pin 21 : no overrange + 12 V to pin 21 : overrange.
48	OY2F	Overrange Y2 servosystem farside.	- 12 V to pin 25 : no overrange + 12 V to pin 25 : overrange.
49	OY2C	Overrange Y2 servosystem close side.	- 12 V to pin 25 : no overrange + 12 V to pin 25 : overrange.
50	GND	Measure ground	not for safety grounding.

1.2.

technical specifications

1.2.1.

Main frame

Paper size	: 420 x 297 mm (DIN A3)
Writing area	: 380 x 250 mm
Writing format	: A3 or A4 Internal switch setting
Paper adjustment	: Indication by three LED's
Paper hold down	: Electrostatic
Pen	: Fibre tip disposable
Pen lift	: Electrical: internal or external (option)
Slewing speed	: X axis: 140 cm/s (BD 90) 140 cm/s (BD 91) Y axis: 140 cm/s
Acceleration	: X axis: 5500 cm/s ² (BD 90) 4000 cm/s ² (BD 91) Y axis: 5000 cm/s ²
Pen response	: X axis: 0,30 s (BD 90) 0,35 s (BD 91) 0-95% fsd Y axis: 0,20 s 0-95% fsd
Input	: Floating, asymmetrical
Input impedance	: 1 MOhm
Range	: 50 mV/cm (X and Y)
Inaccuracy (including non-linearity, deadzone and range)	: 0.2% max.
Overrange	: Electronic limiting of pen travel one mm before mechanical end stops. In overrange, pen is lifted.
Overload	: Automatic power cut-off in case of mechanical overload
Normal mode rejection	: 50 Hz: 50 dB min.

1.2.2.

Input module (identical for X, Y and Y' axis)

Input	: Floating, asymmetrical
Input impedance	: 1 MOhm
Fixed ranges	: 50 μ V/cm - 5 V/cm 16 ranges in 1-2-5 sequence
Inaccuracy	: 0.2% max. of range
Variable mode	: 40 - 100% of fixed ranges
Zero setting	: 2 x 380 (X) or 250 (Y) mm.
Zero drift	: 2 μ V/°C
Polarity	: Direction of pen travel can be reversed; additional range error 0.1% max. in negative direction
Common mode rejection	: 140 dB ac. and dc.
Zero suppression (option)	: 1 to 5 times 200 mm.
Inaccuracy	: 0.4 mm max. for each 200 mm zero suppression

1.2.3.

Rear input/output connector (option)

Remote control of	: All time base functions (option) T.T.L. or switch Pen up-down; T.T.L. or switch Paper hold-free; T.T.L. or switch Record-standby; T.T.L. or switch Muting of X and Y servo systems; switch
Inputs	: X, Y, Y' amplifiers Pulse input time base; 1 pulse $\approx$ 0.1 mm; T.T.L.
Outputs	: 24 V dc - 100 mA max. Overrange indication Time base ramp; 50 mV/cm Monitor X, Y, Y' amplifiers

1.2.4.

Time base module (option)

Panel controls	: XY - tY mode, reset - start of sweep
Sweep factors	: 20-10-5-2-1-0.5-0.2-0.1-0.05-EXT, $\frac{s}{cm} - \frac{min}{cm}$
Sweep nonlinearity	: 0.1% f.s.d. max.