

POLAROGRAPHE modulaire type PRGS

PRGS modular
POLAROGRAPHER

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Three-electrode polarography has many advantages due to the fact that the imposed potential control circuits, and the current measuring circuits are independent. It is possible, in particular, to carry out measurements in highly resistive media (e.g. non aqueous solutions), thus obtaining results which are not affected by the ohmic drop between electrodes.

The availability of recently designed electronic components (field effect transistors, silicon integrated circuits) has helped to create compact, inexpensive and easy to use measuring equipment, nevertheless capable of operating with outstanding performance so far as sensitivity and accuracy are concerned. In addition, the availability of a general purpose plug-in type recorder (EPL series) allows the measurement of other quantities such as potential, pH, conductivities, temperatures, etc. by a simple choice of the appropriate plug-in preamplifier (eighteen units available).

The PRGS modular polarograph is composed of the following :

- A general purpose plug-in recorder of the EPL type ;
- A "TiPol" plug-in, housing all the main and accessory polarographic circuits ;
- A polarographic measuring cell, type CPRS or CPRA, this latter able to accommodate a MPO 3 hammer.
- Furthermore, it is possible to add to the above equipment the following modules :
- GCMS module, for the regulation of the hammer rate controlling mercury drop time, and rapid polarography measurements ;
- ADAPAL adapting unit, for polarography with superimposed AC signal with phase sensitive detection.

Finally, it is possible, by the addition of a syringe-burette or of a coulometric generator, to carry out amperometric titrations with a three-electrode circuit, with or without superimposed AC signal.

The following table gives the various combinations which are possible, depending upon the measurement technique selected.

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
TECHNICAL SPECIFICATIONS

① EPL recorder

The detailed characteristics of the two models of potentiometric recorders will be found in leaflets 12-EPL 1 and 10-EPL 2.

② TiPol plug-in

In this plug-in are grouped : the potentiostat ; the circuits and controls for potential sweeps, current measurement, damping ; various complementary controls : slope adjustment control (capacitive current compensation), recorder control, etc.

Potentiostat :

- Drift : less than 1 mV per day (non cumulative) ;
- Maximum output voltage : ± 10 V ;
- Maximum output current : ± 5 mA ;
- Maximum power : 50 mW ;
- Risettime : less than 50 μ s ;
- Warning light signalling overload, open circuit or erroneous connection.

Potential excursion control :

- Initial potential : adjustable between + 5 and - 5 V (resolution : 5 mV) ;
- Forward and reverse sweep ;
- Sweep speed : 15mV.min⁻¹ to 6000 mV.min⁻¹ in nine ranges (progression : 15 - 30 - 60 - 150, etc.).

Current measuring circuits :

- Provision for setting the zero at any part of the recorder scale, or to offset it from 0 to 50 % beyond the left or right scale limits ;
- Twelve measuring ranges, from 0.25 μ A to 1250 μ A for the full inscription width (250 mm = 10 ") of the recorder ; progression : 0.25 - 0.5 - 1.25 - 2.5, etc. Measurements up to 5 mA by addition of an external, calibrated resistor.

Damping :

- Time constant adjustable from 0 to approximately 5 seconds in six steps.

Slope compensation :

- Progressively adjustable between 0 and 2 μ A.V⁻¹.

Complementary circuits :

- Control for the automatic synchro start between chart drive mechanism and potential sweeps ;
- Output delivering a low-impedance signal of the same value as the potential imposed on the dropping electrode, thus allow-

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
CPRA

③ CPRS measuring cell

This cell, of very simple design, includes :

- A measuring vessel (allowing operation on a typical volume of solution of 10 ml) ;
- A rubber stopper drilled for insertion of the electrodes and the deoxygenation tubing ;
- A mercury reservoir, fitted with a PTFE plunger stopcock ;
- The supporting hardware.

④ CPRA measuring cell

- "All-glass" tight construction, allowing measurements in nonaqueous, corrosive or toxic media ;
- Accurate indication of the mercury height with respect to the capillary tip ;
- Mercury supply stopcock with platinum contact ;
- Gas input tube ;
- Rapid fitting, thermostatically controlled measuring vessel ; permits measurements on volumes of 5 to 50 ml ;
- Possibility to adapt the MPO 3 hammer (see below : "GCMS pulse generator").

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COMPLEMENTARY MODULES

① GCMS* pulse generator for hammer control

- Mercury drop time adjustable between 0.1 and 9.9 sec., by 0.1 sec. steps ;
- Electronically adjustable stroke intensity ;
- Output for the possible synchronization of an oscilloscope to the hammer strokes ;
- Dimensions : 120 x 90 x 250 mm (4 3/4 x 3 3/4 x 10") ;
- Weight : approximately 3 kg (6 3/4 lbs).
- This instrument is intended to be used in conjunction with hammers of type MPO 1 or MPO 3.


② ADAPAL* unit for polarography with superimposed AC signal

The ADAPAL unit is to be used with the PRGS polarograph, the GCMS hammer pulse generator and the CPRA measuring cell fitted with the MPO 3 hammer ; adjustments and controls relative to the potential excursion, the damping, the current origin and the starting of the measurements are effected from the Tipol plug-in ; the ADAPAL unit includes the superimposed voltage source, and the current amplification and phase sensitive detection circuits.

- Frequency of the superimposed signal : 20 Hz ;
- Amplitude of the superimposed signal : 10 mV ;
- Measurement of the AC current through the cell : twelve ranges, between 50 nA and 250 μ A (for the 250 mm or 10"


* This module is designed for 50 Hz operation ; upon special order, it can be fitted with special circuits for 60 Hz mains supply.

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

- Synchronous detection, Demodulation phase angles : 0 and $90^\circ \pm 3^\circ$ (continuous adjustment) ;
- Various outputs for oscilloscopic display ;
- Dimensions : 120 x 130 x 260 mm (4 3/4 x 5 1/4 x 10 ") ;
- Weight : approximately 3.5 kg (7 3/4 lbs.) .

Polarography with superimposed sinewave and synchronous detection marks a noticeable progress, by the increase of both resolution and measuring sensitivity. Considering the last criterion, one can carry out analyses of trace amounts of material at concentrations more than ten times lower than those detectable by conventional polarography, or by the technique of superimposed AC signal without synchronous detection.

$$\begin{aligned} \text{Pb}^{++} & 10^{-4} M ; \\ \text{Cd}^{++} & 2 \cdot 10^{-4} M ; \\ \text{Zn}^{++} & 10^{-4} M. \end{aligned}$$

Cu^{++} $5 \cdot 10^{-5} M$;
 Pb^{++} $10^{-4} M$;
 Tl^{+} $2 \cdot 10^{-4} M$;
 In^{+++} $5 \cdot 10^{-5} M$;
 Cd^{++} $10^{-4} M$;
 Zn^{++} $2 \cdot 10^{-4} M$.

Electrolyte support : $\text{HCl } 10^{-1} N$.

Pb⁺⁺ 10⁻⁵ M ;
Cd⁺⁺ 10⁻⁵ M (en présence de Zn⁺⁺ 10⁻³ M).

Electrolyte support : $\text{HCl } 10^{-1} \text{ N.}$