
POTENTIOSTATS

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Herein is described the complete line of TACUSSEL electronic potentiostats, including improved models and new additions. There is a wide choice of models, allowing applications ranging from the most modern electrochemical relaxation techniques to more classical techniques, such as corrosion study and control, electrochemical compound preparation, controlled potential separation, polarography. The power capabilities offered should meet with most standard requirements ; facilities are available to construct special types of potentiostats for specific purposes.

Electronic potentiostats are versatile instruments, which can be voltage programmed and adapted to various closed loop configurations, thus differing from ordinary power supplies. The speed with which a potentiostat may follow a program voltage source depends upon the model type. TACUSSEL potentiostats are noted for their fast (99 %) rise time and provision for phase shift control, thus allowing a considerable range in electrochemical cell conditions to be adequately handled for a particular type of study or control.

Due to its electronic configuration, a potentiostat may be employed in the potential control, current control or voltage control modes : TACUSSEL potentiostats are particularly adapted to all of these operating conditions. The circuitry is such that control accuracy is maintained, whether the unit is operating at extremely low or near saturation output power.

CONTENTS


Introduction	2
Contents	2
Operating principles :	
1 - Use as a potentiostat	3
2 - Other uses :	
- use as a current source (galvanostatic operation)	3
- use as a voltage source	3
PRT and PIT series	
General features and recapitulative table	4
Common characteristics	5
Particular characteristics :	
- Types with interchangeable plug-ins	
(PRT...X and PIT...X series)	6
a) zero setting plug-ins (TZZ and TZM types)	6
b) Other plug-ins (TCCO, TBT, TCS types)	7
- Types without plug-ins (PRT series)	7
- Pulse power potentiostats	
(PRT...X series)	8
- Modular potentiostatic units	
(PIT...MOD series)	9
ASA series	10
Accessories and ancillary equipment	12

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OPERATING PRINCIPLES

1. Use as a potentiostat

The basic function of a potentiostat is to maintain the potential of an electrode at a controlled (imposed) value, fixed or variable, in spite of variations of the other factors in the electrochemical system in which it is included (fig. 1).

The potential of that electrode, E_t , called the working electrode (or controlled electrode), is measured against another electrode, E_r , called the reference electrode, which is placed in the electrochemical cell, very close to the working electrode. A third electrode, E_a , called the auxiliary electrode (or counter-electrode) allows the electrical current necessary to produce the wanted potential difference to flow through the circuit. These three electrodes are schematically represented as shown on the adjacent diagram.

The working electrode and the auxiliary electrode are connected to the output terminals of the potentiostat. Current through this circuit is automatically controlled so that the potential difference between the working electrode and the reference electrode takes the desired value. This control process is carried out in the potentiostat by means of a differential amplifier, A_d , one input of which, e_1 , is connected to the reference electrode and the other input, e_2 , to a voltage source called the pilot voltage (or control voltage). This amplifier drives a power amplifier, A_p , which controls the output current of the potentiostat in such a manner that the potential difference between the working electrode and the reference electrode remains equal to the applied control voltage E_c . This voltage may vary according to any program within the rise time limits of the potentiostat employed.

2. Other uses




USE AS A CURRENT SOURCE (galvanostatic operation, fig. 2)

A resistor R of appropriate value is series connected with the load circuit Z ; the common point between elements Z and R is connected to the e_1 input of the amplifier. The potentiostat output current, i , which flows through the load circuit Z and the sampling resistor R , creates between the resistor terminals a voltage drop, iR , which is equal to the control voltage E_c , as a result of the control process. If the E_c voltage as well as the resistor R values are properly selected, the current - within voltage and currents limits stated for each model - is constant whatever the variations of the load impedance Z and of the mains supply voltage may be.

The current value, i , may be varied according to any program, by action on the value of the voltage E_c (or of the resistor R). Such a working mode is involved, in particular, in programmed chronoamperometry and galvanostatic techniques.

USE AS A VOLTAGE SOURCE (fig. 3)

This operating mode is used to control the voltage across the cell in a two-electrode system (both output terminals and e_1 input of the potentiostat being directly connected to the load circuit Z). Within voltage and current limits stated for each model, output voltage is then equal to control voltage, whatever the current drain and the mains supply voltage variations may be. As in the previous case, it is possible to vary the voltage value E_c according to any programme.

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PIT and PRT SERIES

Fast rise time potentiostats of the PRT and PIT series are fitted with a direct coupled dc amplifier with a transistor output stage.

Their electronic design is such that the polarity of control voltage, output voltage, and output current, can independently be reversed without switching or reversal of connections. In addition to the conventional use of any potentiostat (as controlled potential electrolysis current generator), these instruments are thus particularly useful for the manual or automatic plotting of current/voltage curves of electrochemical systems, for potentiokinetic studies as well as for three-electrode polarography and programmed voltammetry.

Their minimum rise time is very short : 0.1 to 1 μ s for PIT series, 1 to 10 μ s for PRT series, which is quite suitable for pulse techniques.

The basic electronic circuitry is identical for PRT, PRT...X and PIT...X potentiostat types. Instruments of the first two series are fitted with the same amplifier and differ only in their physical construction (the PRT...X being designed to accept plug-ins) ; instruments of the PIT series include particular electronic arrangements giving them a specially short response time. In each category, several models are available, covering a large range of output and current values, as shown in the following table :

The various potentiostats of the PRT and PIT series have some common characteristics, hereafter listed. They can be classified in four main groups : PRT ; PRT...X ; PIT...X ; PRT...MOD.

The first group contains instruments as simple as possible in their design and physical features, and moderately priced. Instruments of the PRT...X and PIT...X groups accept plug-ins which determine some characteristics and/or perform complementary functions. Potentiostats of the PRT...MOD (modular) group are representative of a new trend in electrochemical instrumentation ; they are mainly intended to build up complex and versatile systems.

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COMMON CHARACTERISTICS

● Pilot control voltage continuously adjustable from - 5 to + 5 V by means of a built-in 10-turn helical potentiometer. The potentiometer shaft, accessible through the knob, is machined so as to allow it to be driven at a constant speed by a SYNCHROVOLT plug-in auxiliary motor. The PIT...X models are to be externally controlled by a function generator, or by a potentiometer (which can be connected to the ± 5 V internal regulated voltage source intended for this purpose).

▲ Socket for connecting an external control generator (pilot scanner), for time-voltage programming, or, more generally, for controlling the output voltage by any variable phenomenon as translated into an electrical signal. Direct connection of the "cold" terminal of the external generator to the potentiostat ground (thus avoiding "ground loops" and parasitic signals at hf or mains frequency).

- Maximum amplitude of control voltage : ± 10 V ;

- Frequency response : dc to approximately 500 kHz (several megahertz for PIT series).

● Regulation accuracy typically better than $\pm 0,5$ mV (for 0 to 100 % of load and ± 15 % of mains voltage variations). For instruments fitted with a ZEROSTAT plug-in, typical accuracy is better than ± 50 μ V.

● Maximum noise and hum : less than 50 μ V rms referred to input.

● Amplifier input impedance : of the order of 10^6 M Ω .

● Built-in limiting circuits insuring full protection against overloads, short-circuits and misuse.

▲ Built-in center-zero meter for output voltage monitoring with ± 500 mV range push-button and overload protection.

■ Built-in center-zero ammeter for output current monitoring ; socket for external ammeter.

▲ Socket for connecting an external millivoltmeter, with high frequency decoupling for correct pulse operation (for the accurate measurement of the true working electrode potential).

◆ Built-in phase correction circuit allowing the potentiostat frequency response to be adapted to the characteristics of the electrochemical system under study in order to obtain optimum rise time and to avoid parasitic oscillations when the double layer capacity is high (1). On PRT and PRT...X potentiostats, this adjustment is made by means of a four-position bandwidth limiting switch : on the first time position, minimum response time is from 1 to 2 μ s ; on the three others, the response time is respectively about 20 μ s, 200 μ s and 2 ms. On PIT...X models, the correction is made by a continuously adjustable R-C circuit. External correction network can be used with all types of PRT and PIT potentiostats.

● Power requirements : 110 to 250 V ± 15 %, 40 to 60 Hz.

● Physical features : standard 19" rack mounting cabinet, two-tone grey enamelled.

▲ Except for PRT 10-05 S type.

■ Except for PRT 10-05 S and PIT types.

◆ For PRT 10-05 type, the time constant switch has only two positions : minimum (1 to 2 μ s) ; approximately 1 ms.

(1) The following papers deal with conditions to be fulfilled in order for a system containing a potentiostat and an electrochemical cell to be stable :

J. Schoen, K.E. Staubaek, Regelungstechnik, 1954, 2, p. 71 Sqg.
Some factors in potentiostat design, A. Bewick, M. Fleischmann, M. Liber, Electrochimica Acta, 1959, 1, p. 83 - 106.
The design and performance of potentiostats, A. Bewick et M. Fleischmann, Electrochimica Acta, 1963, 8, p. 89 - 106.
Facteurs déterminant le temps de réponse d'un potentiostat, J. Tacussel, Electrochimica Acta, 1966, 11, p. 437 - 42.

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PARTICULAR CHARACTERISTICS

In addition to individual characteristics (numerical values) and remarks relative to these various models of potentiostats, which are presented in the above table (page 3), the following particularities are to be noted :

Plug-in models (PRT...X and PIT...X series)

Potentiostats of PRT and PIT series including the letter X in their catalog number are to be provided with one or two interchangeable plug-ins. The type of these plug-ins depends upon zero adjusting mode and complementary characteristics desired.

Plug-in type potentiostats PRT 20-2 X, PRT 20-10 X, PRT 10-20 X, PIT 20-2 X, PIT 20-6 X respectively replace the discontinued models PRT 20-2 Z, PRT 20-10, PRT 7-20, PIT 20-2, PIT 20-6. The new amplifier with which they are fitted includes a symmetrical input stage with a pair of tubes working in electrometric conditions, followed by loading and coupling stages which drive a complementary transistor output stage. All semiconductors used are professional silicon types. Various advantages result from this new design : reduction of size and of power requirement ; improvement in regulation (approx. 10 times) and stability (approx. 3 times) ; reduced warm-up time ; pilot circuit input suitable for signal generators having one output terminal grounded.

a. Zero setting plug-ins

TZZ PLUG-IN :

It includes a carrier frequency amplifier (ZEROSTAT circuit) insuring automatic drift correction of the potentiostat direct coupled amplifier. The stability thus obtained is typically better than $\pm 50 \mu V$ per day (non cumulative drift). The time constant of the R - C smoothing circuit included in the output filter of the demodulator stage can be given two values in a 10/1 ratio, in order to obtain the optimum response for pulse operation.

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T ZM PLUG-IN :

This plug-in allows manual zero adjustment to be made. In practice, for electrochemical systems with high double layer capacitance, it can be advantageous to eliminate the supplementary RC time constant (of an order of magnitude of some seconds) introduced in the negative feedback loop by the ZERO-STAT circuit. A better transient response is thus obtained for relatively slow phenomena (i.e. having a duration from 100 μ s to 1 s), which is interesting, particularly in programmed voltammetry and in oscilloscopic polarography. The T ZM plug-in is fitted with two zero adjustment potentiometers (coarse and fine) ; the drift of the potentiostat amplifier, when this plug-in is inserted, is typically ± 1 mV per day.

When ordering, catalog number of the instrument will be completed by the suffix "M" or "Z" depending upon whether the desired plug-in is the T ZM or T ZZ type (for instance, a PRT 20-2 X potentiostat fitted with a T ZZ plug-in will be designated as PRT 20-2 X - Z).

b. Other plug-ins

As previously mentioned, potentiostats of the PRT...X series also include, on the front panel, a housing for another (optional) plug-in. This plug-in can be a TCCO type, which includes a simplified oscilloscope internally connected to the potentiostat output terminals in order to display interference signals (mains frequency residue or H.F. oscillations), and a continuously adjustable stabilization (phase correction) network. In the case when TCCO plug-in is not ordered, the instruments are normally supplied with a TBT blank plug-in. The housing has also been arranged to accept a function generator plug-in, presently in course of development.

In addition to the zero adjustment plug-in (TZZ or T ZM type), potentiostats of the PIT...X series can accept (and are normally delivered with) a special plug-in (type TCS) for stability correction. These potentiostats are generally used in connection with an external broad-band oscilloscope ; therefore no provision has been made for a CRT plug-in such as TCCO.

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PRT 10-05 L

Types without plug-ins (PRT series)

The following models have only slightly been modified in their external appearance : PRT 10-05 L (ex PRT 500 L), PRT 10-05 S (ex PRT 500 S), PRT 20-2 A (which replaces PRT 20-2 and PRT 2000 types). Internally, however, they have benefited from all electronic improvements developed for PRT...X and PIT...X types. Adjacent illustrations are those of PRT 10-05 models ; that of PRT 20-2 A is on the right top of the cover.

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Potentiostats of the PIT 20-2 X and PIT 20-6 X types, previously mentioned with regard to the plug-ins they can be fitted with, have specially been designed for ultrafast pulse techniques in electrochemical systems requiring currents of medium or high value (fuel cells, reactions in fused salts, etc.).

They include a high gain amplifier with only one voltage amplifying stage equipped with a video pentode, a buffer stage, and a complementary output stage fitted with high frequency power transistors. The amplifier transmitting the control voltage pulse from the pilot generator is also fitted with a single tube associated with transistors, according to the same basic electronic configuration as the main amplifier, which results in very low drift.

The pulse response of the PIT...X series instruments allows rise and fall times of 50 to 100 ns to be obtained ; they can also be controlled by sinewave signals, the frequency of which can be as high as 10 MHz, for small amplitudes.

Their main characteristics (some of which have already been mentioned above) are summarized hereafter :

- High impedance symmetrical input circuit ($\approx 10^{12} \Omega$).
- Open loop gain : with T2M plug-in : $3 \cdot 10^3$; with T2Z plug-in : greater than 10^6 (typical dc values). Gain bandwidth product : 20 MHz min.
- Adjustable plug-in phase correction circuit.
- Maximum amplitude of pilot signal : greater than $\pm 10 \text{ V}$.
- Response time (99 %) : less than 100 ns at full resistive load.
- Low impedance auxiliary output for the connection of an oscilloscope.

Several publications had been devoted to the electronic design of these instruments and to the experimental results obtained through their use (2) (3) (4) (5).

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3) A transistorized potentiostat with very fast rise time, J. Tacussel, Electrochimica Acta, 1966, 11, p. 449-55.

4) La méthode de la préimpulsion potentiostatique, J. Tacussel, Electrochimica Acta, 1966, 11, p. 443-7.

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We have recently undertaken design and production of modular units which can be associated together and with one or several basic units in order to constitute a complete electronic apparatus, particularly adapted to the solution of a specific problem, but easily modified for other applications by adding and/or rearranging modules (patents pending). Numerous uses are possible : Potentiostatic and potentiokinetic techniques. - Corrosion studies. - Controlled potential coulometry. - Automatic titrations. - Three-electrode polarography. - Oscilloscopic polarography. - Programmed voltammetry. - Chronoamperometry, etc.

For corrosion studies or differential polarography, the association of PRT 20-002 MOD units with a MAS 24-02 (mains supply) or a MAB 24 (storage battery supply) unit results in a "multiple" potentiostat, each channel having the following characteristics :

- Output voltage : 0 ± 20 V.
- Maximum output current : ± 20 mA.
- Maximum output power : 0,4 W.
- Input stage fitted with field effect transistors (the instrument being otherwise entirely fitted with high grade silicon transistors).
- Response time : approximately 10 μ s.
- Control voltage continuously adjustable from -5 to +5 V by means of a built-in 10-turn helical potentiometer.
- Socket for connection of an external pilot voltage source (permissible control voltage : ± 10 V).
- Regulation accuracy better than ± 1 to ± 2 mV (typical value : 0,5 mV) for 0 to 100 % load and ± 15 % of mains supply variation (or for a voltage variation from 22 to 27 V in case of battery supply).
- Noise and ac residual voltage : less than 50 μ V referred to input.
- Socket for connection of an external milliammeter.

Two more powerful models (output currents of, respectively, ± 100 and ± 300 mA) are in course of development.

OTHER PRT...X SERIES POTENTIOSTATS IN PREPARATION

Three high power models will soon be added to the PRT...X series : PRT 20-25 X, PRT 20-50 X and PRT 20-100 X (with an output voltage of ± 20 V and a maximum current of, respectively ± 25 , ± 50 and ± 100 A) and by a high output voltage model : PRT 200-1 (± 200 V, ± 1 A).

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These instruments are fitted with a carrier frequency amplifier, working at a frequency of 2 to 5 kHz ; as a consequence, their response time is of the order of a few milliseconds. The power output stage feeds the load circuit by means of a transformer and of silicon rectifiers. By proper choice of transformer ratio and of rectifier ratings, high output voltages or currents can easily be obtained and, if necessary, specific instruments can be built which, for a given output power, can afford combinations of voltage and current values different from those indicated in the following table.

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The various types of 500 W potentiostats differ only in their power amplifier, the adapting transformer and the output circuit rectifiers, which are located in a rack-mounting unit. The other elements (carrier frequency amplifier, control circuits, power supplies, overload protection devices, ventilating unit) also housed in rack mounting cases, are assembled in a standard 19" rack. It is thus possible for the user, with a single basic system and several power amplifier units, to dispose of different output voltage and current combinations (10 V and 50 A ; 20 V and 25 A ; 50 V and 10 A) each one being suitable to a given application. Other output units can be designed on special order.

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TYPE 100 W

COMMON TECHNICAL CHARACTERISTICS

• Control voltage continuously adjustable from 0 to ± 10 V by means of a built-in 10-turn helical potentiometer. The potentiometer shaft, accessible through the knob, is machined so as to allow it to be driven at a constant speed by a SYNCHROVOLT plug-in auxiliary motor.

• Socket for connecting an external control generator (pilot scanner) for time-voltage programming, or, more generally, controlling the output signal by any variable phenomenon as translated into an electrical signal.

- Maximum control voltage : ± 10 V.
- Frequency response : dc to approximately 100 Hz.

• Regulation accuracy better than ± 5 mV (typical value : 1 to 2 mV) for 0 to 100 % of load and $\pm 15\%$ of mains voltage variations.

• Maximum noise and hum : less than 2 mV.

• Amplifier input impedance : greater than $10^{11} \Omega$.

• Built-in offset voltage source : 0 ± 10 V.

• Full protection against short circuits, overloads and misuse.

• Built-in indicator for output voltage monitoring, protected against overloads.

• Supply : ac 110 to 250 V $\pm 15\%$; 40 to 60 Hz.

• Physical features : standard 5-unit cases for 19" rackmounting, two-tone grey enamelled.

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