# POTENTIOSTATS

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PRT 40-1X  
+ TZZ  
+ TCCO

PRT 20-2X  
+ TZZ  
+ TP - PRT

*Fast risetime*  
**POTENTIOSTATS**  
*with plug-ins (PRT... X series)*

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The fast response potentiostats of the PRT...X series described in this leaflet are entirely transistorized and equipped with a direct coupled DC amplifier. They can be fitted with one or two interchangeable plug-ins, which determine the zero adjustment mode and some complementary characteristics.

Their output power covers a large range of voltage and current values ; for certain particular applications, it is possible, upon special order, to design and build non-standard models.

In addition to their basic usage, which is to maintain at a determined value, fixed or adjustable, the potential of a working electrode, independently of the other factors likely to vary within the electrochemical system in which it is placed, these potentiostats can be used as sources of stabilized voltages (fixed, adjustable or programmable) and, with a sampling resistor inserted in the external circuit, as current generators. In the latter case, the current can be either fixed, or controlled, by variation of the pilot voltage, or by proper adjustment of the sampling resistor value.

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The electronic design of potentiostats of the PRT series is such that the polarity of control voltage, output voltage, and output current, can independently be reversed without switching or reversal of connections. In addition to the conventional use of any potentiostat (as controlled potential electrolysis current generator), these instruments are thus particularly useful for the manual or automatic plotting of current/voltage curves of electrochemical systems, for potentiokinetic studies as well as for three-electrode polarography and programmed voltammetry.

The amplifier they are fitted with includes a symmetrical input stage equipped with high transconductance field effect transistors, followed by loading and buffer stages which drive a complementary transistor output stage. All semiconductors used (field effect and bipolar transistors, diodes, integrated circuits) are silicon professional types. The electronic circuits are characterized by a very short warm-up time (2 to 3 seconds), a high zero stability (better than  $\pm 1$  mV), a low thermal dissipation and power consumption.

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\* This type is fitted with a digital display panel meter for measurement of output voltage and current, as well as of working electrode potential.

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- Pilot control voltage continuously adjustable from - 5 to + 5 V by means of a built-in 10-turn helical potentiometer. The potentiometer shaft, accessible through the knob, is machined so as to allow it to be driven at a constant speed by a SYNCHRO-VOLT plug-in auxiliary motor.

- Socket for connecting an external control generator (pilot scanner), for time-voltage programming, or, more generally, for controlling the output voltage by any variable phenomenon as translated into an electrical signal. Direct connection of the "cold" terminal of the external generator to the potentiostat ground (thus avoiding "ground loops" and parasitic signals at RF or mains frequency).

- Maximum amplitude of control voltage :  $\pm 20$  V ;
- Frequency response : dc to approximately 500 kHz.

- Regulation accuracy typically better than  $\pm 0,5$  mV (for 0 to 100 % of load and  $\pm 15$  % of mains voltage variations). For ins-

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truments fitted with a ZEROSTAT plug-in typical accuracy is better than  $\pm 50 \mu$ V.

- Maximum noise and hum : less than 50  $\mu$ V rms referred to input.
- Amplifier input impedance : of the order of  $10^{12} \Omega$ .
- Built-in limiting circuits insuring full protection against overloads, short-circuits and misuse.
- Built-in center-zero meter for output voltage monitoring, with  $\pm 500$  mV range push-button and overload protection.
- Built-in center-zero ammeter for output current monitoring ; socket for external ammeter.
- Socket for connecting an external millivoltmeter, with high frequency decoupling for correct pulse operation (for the accurate measurement of the true working electrode potential).
- Built-in phase correction circuit allowing the potentiostat frequency response to be adapted to the characteristics of the electrochemical system under study in order to obtain optimum rise time and to avoid parasitic oscillations when the double layer capacity is high. This adjustment is made by means of a four-position bandwidth limiting switch : on the first position, minimum response time has the value indicated in the above table ; on the three others, the response time is respectively about 20  $\mu$ s, 200  $\mu$ s and 2 ms. External correction network can be used.
- Power requirements : 110 to 250 V  $\pm 15\%$ , 40 to 60 Hz (200 to 250 V if power requirement is higher than 500 VA).
- Physical features : standard 19 in rack mounting cabinet, grey enamelled ; front panel of anodized aluminum.

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### a - Zero adjustment plug-ins

#### TZZ PLUG-IN

This plug-in includes a carrier frequency amplifier (ZEROSTAT circuit) insuring automatic drift correction of the potentiostat direct coupled amplifier. The stability thus obtained is typically better than  $\pm 50 \mu$ V per day (non cumulative drift). The time constant of the circuit included in the output filter of the demodulator stage can be given two values in a 10/1 ratio.

#### TZM PLUG-IN

This plug-in allows manual zero adjustment to be made. In practice, for electrochemical systems with high double layer capacitance, it can be advantageous to eliminate the additional RC time constant (which has an order of magnitude of some seconds) introduced in the negative feedback loop by the ZEROSTAT circuit. A better transient response is thus obtained for relatively slow phenomena (i.e. having a duration from 100  $\mu$ s to 1 s), which is interesting, particularly in programmed voltammetry and in oscilloscopic polarography. The TZM plug-in is fitted with two zero adjustment potentiometers (coarse and fine) ; the drift of the potentiostat amplifier, when this plug-in is inserted, is typically  $\pm 1$  mV per day.

■ When ordering, catalog number of the instrument will be completed by the suffix "M" or "Z" depending upon whether the desired plug-in is of the TZM or TZZ type (for instance, a PRT 20-2 X potentiostat fitted with a TZZ plug-in will be designated as PRT 20-2 X-Z).

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TCCO

## b - Other plug-ins


As previously mentioned, potentiostats of the PRT... X series also include, on the front panel, a recessed area which can accommodate another plug-in :

### TCCO PLUG-IN

This unit includes a simplified oscilloscope, internally connected to the potentiostat output terminals in order to display interference signals (mains frequency residue or RF oscillation), and a continuously adjustable stabilisation (phase correction) network.

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TP - PRT

### TP-PRT PLUG-IN

This plug-in is intended to generate a time dependant pilot signal. This signal is constituted by the sequence, single or recurrent, of three straightline segments, the slope of which can be made positive or negative. By proper adjustment, it is possible to generate a sequence of only two segments. In the non repetitive mode, the sequence can be triggered either by means of a push-button (one-shot manual control) or by an external electrical signal.

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The adjustment of each segment parameter (slope, limits) and of the offset voltage (pedestal, or starting potential) is made by rotary switches and ten-turn helical potentiometers with turns-counting dials. The main characteristics are as follows :

- Maximum output voltage :  $\pm 5$  V ;
- Four-position decade attenuator ;
- Offset voltage :  $\pm 5$  V, adjustable by a ten-turn potentiometer (pedestal).
- Slope of each segment : adjustable from  $\pm 10^2$  V/s to  $\pm 10^5$  V/s by a four-decade attenuator and a ten-turn potentiometer. Corresponding durations for the maximum voltage sweep :  $10^3$  to  $10^{-4}$  s ;
- Limit points adjustable, for each segment, from  $+5$  to  $-5$  V respective to the starting voltage (pedestal), by a ten-turn potentiometer ;
- Built-in power supply (automatically connected to mains by means of the plug-in connector) ; power requirements : approximately 12 VA.

### TBT PLUG-IN

This blank plug-in is normally supplied when the instrument is not fitted with either a TCCO or TP-PRT plug-in ; it can be used for the storage of accessories and cables.

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