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**Wulf Electroscope,
 Spare Sensitive System,
 Projection Attachment**

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The Wulf electrostatic projector is an apparatus which is used for many purposes, ranging from the simplest electrostatic experiments to measurements on radioactive substances. Each experiment can be rendered visible to a large circle of students by means of the projection attachment or a Leitz projector PRADO or „pradovit“, or optical accessories.

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1. Description

The case of the **Wulf electrostatic projector** (546 00) is sealed at the front and at the rear by glass plates, and has two openings, one at the top and the

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Notes

1. The five-digit numbers quoted in parentheses refer to the catalogue numbers of the respective apparatus.
2. The data DK... refer to the experiment descriptions published as a collection of leaflets in the „LEYBOLD PHYSICS LEAFLETS“. Published to date:

	Cat. No.
1st series (72 sheets)	599 21
2nd series (52 sheets)	599 22
3rd series (49 sheets)	599 23
4th series (50 sheets)	599 24

3. The specifications and illustrations are not binding in every detail for the design of the apparatus. It is our policy always to keep our manufacturing programme right up to date so that it makes full allowance for the latest knowledge acquired in all scientific and technical fields.

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other at one side. An adjustable carbon plate, the so-called counter-electrode (E) (Fig. 2), is fitted in the lateral opening and insulated against the case. The electrode is held in its respective position by the screw (A), and the adjustment is rendered reproducible by means of the setting ring (B). For cleaning, the insulator and counter-electrode can be removed after slackening 2 grub screws.

A large insulator which can be taken out is fitted in the top opening, and this is provided with a boring to carry the system (S). The bottom of the case is provided with a 4 mm boring for earthing (K_4) and two grooves for locating the electroscope in the rods of a Leitz PRADO or "pradovit" assembled as a demonstration projector.

The bottom of the case is also provided with a shaft, 10 mm ϕ and 10 cm long, which can be screwed off so that the apparatus can be fastened in an extension case supplied with it and secured with the latter in a stand base (300 02) or in a boss. The projection attachment can be fitted on this shaft, too. The shaft is provided with a 4 mm boring for earthing (K_3).

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The system is assembled as follows (Fig. 3):

One end of a thin aluminium leaf is secured directly to the system carrier, the other via a fine quartz loop. It can be seen extremely well when projected, since it is turned at an angle of 180° between its ends.

As a result of its special arrangement, which can be seen in Fig. 3, the loop allows the centre part of the leaf to

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deflect up to 20 mm. The deflection for a certain charge depends upon the position of the counter-electrode. The sensitivity of the Wulf electro-scope increases somewhat quadratically to the laid-on voltage until a labile position is finally reached. With each further increase in voltage, the leaf strikes the counter-electrode and discharges.

These characteristics of its system render the Wulf electro-scope suitable for special measuring methods (see Section 3). Other advantages of the system are its almost inertialess adjustment, its reproducibility, and its insensitiveness to vibrations.

A spare system in an aluminium box is supplied with each electro-scope. Additional **spare systems** can be obtained under Cat. No. 546 50. A capacitor plate with a 4 mm plug is also supplied with the Wulf electro-scope. When carrying out electrostatic experiments, this plug should be inserted in socket K₁.

The **projection attachment** (546 02) is used for projection if no Leitz projector is being used, and consists of two short stand rods located axially to each other and held by a ring-shaped boss with socket, and an adjusting screw. A projecting lens is slid over the one column, and an adjustable cross clamp over the other to secure the lamp housing.

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2. Assembly for projection

a) Projection by means of the projection attachment

The projection attachment (546 02) is slid on the shaft of the electro-scope, and the electro-scope is

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upon the position of the counter-electrode. With microscopic observation, and when the counter-electrode is approached to the system, exact measurements of up to 20 V can be made.

b) Circuit with auxiliary charge

The sensitivity of the leaf always reaches a maximum just prior to the leaf deflecting against the counter-electrode. Therefore, the loaded, insulated system is also suited for measuring small voltages of about 1 V. The system (K_1) should be loaded either by contact or by influence, and the voltage which is to be measured should be laid on either between (K_2) and (K_3) or (K_4), respectively, or with the counter-electrode charged between (K_1) and (K_3) or (K_4), respectively. With this circuit, it is possible, for example, to demonstrate an induced voltage pulse in a coil electrostatically.

c) Drop method

This method is not designed to measure a fixed deflection but to observe the leaf striking against the counter-electrode. Since the charge which is taken away is the same for each deflection, the number of deflections in a certain space of time gives the measurement of the amperage. The method is particularly suited for demonstrating the ionization of the air by radioactive rays, by the photo-effect, or by flame gas. The amperage is a measurement for the number of ions formed. Detailed instructions on this circuit are given, among others, in the Directions for Use of the ionization chamber (546 25), the Geiger counter (546 28), and the pair of spiral electrodes (546 29).

4. Precautionary measures

The system of the Wulf electroscopes is almost insensitive to vibrations, so that it is not necessary to take any special precautions when transporting the apparatus. On the other hand, we would draw your attention to several points which experience has shown to be the cause of damage to the system.

- When taking the system out of the Wulf electroscopes, the system may be easily damaged due to the counter-electrode being close to it. These manipulations should, therefore, only be carried out with special care and when the counter-electrode has been withdrawn.
- Also when mounting and removing the supplementary appa-

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secured in a stand base (300 02) either directly (Fig. 1) or via a stand rod, 25 cm long (300 41), by means of a boss or clamp (301 01 or 06) (Fig. 4). To illuminate the apparatus, place the lamp housing (450 60) with fitted single-lens condenser (460 17) in the cross-clamp and slide it right up to the electro-scope. The lens, electro-scope and lamp should all be located at the same level; the projection screen can be placed 2 to 8 metres away. The well-defined image of the system is produced by sliding the lens. Finally the lamp filament should be placed in the vertical, and the lamp should be moved in the housing until the lens is uniformly illuminated.

b) Projection by means of a Leitz projector (Figs 5 and 6)

One of the following appliances is required as a projector: Leitz PRADO 250 or 500 assembled as a demonstration projector (Cat. No. 442 50 or 51), or the Leitz "pradovit" n 12 (442 91) or n 24 (442 93) with horizontal attachment (442 58) and holder thereto (442 95). Besides, the condenser for 20 to 30 cm focal length is used.

c) Projection by means of optical accessories (Fig. 7)

The Wulf electro-scope is assembled with the aid of 1 stand base (300 02), 2 rods, 25 cm long (300 41), and 4 Leybold bosses (301 01), as shown in Fig. 7.

The lamp-house with lamp (450 60/51) and the single-lens condenser (460 17) are suited for illumination; for projection, a lens $f = 100$ mm or $f = 200$ mm (460 03 or 04) can be used, dependent upon the distance of the projec-

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tion wall. This assembly can also be made on the small optical bench (460 41/42).

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3. Operation

The system is in the correct position when the space between the carrier and the leaf can be clearly seen.

The Wulf electrostatic system is connected to the mains by means of three sockets (see Fig. 2): K_1 is connected to the system, K_2 to the counter-electrode, and K_3 or K_4 , respectively, to the case.

The Wulf electrostatic system is used for three measuring methods:

a) Circuit without auxiliary charge

The counter-electrode is earthed by connecting K_2 to K_3 or K_4 , respectively. The voltage which is to be observed (e. g. that of the sliding parallel-plate capacitor 544 22) is laid on between the system (K_1) and the earth. The system can also be charged either by contact or by influence.

When the charge reaches a certain value, dependent upon the position of the counter-electrode, the leaf of the system will deflect against the counter-electrode, and will be discharged again. The deflection of the leaf can also be utilised for charging by influence. A rubbed rod made of glass, hard rubber or plastic material should be placed close up until the leaf deflects. The counter-electrode should then be withdrawn a little, and the rod removed.

With the circuit without auxiliary charge the measuring range is between 200 and 1000 V, dependent

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ratus, the system may be drawn out of the case and be damaged. In order to avoid this, the system should always be held with one hand by its insulator.

- c) When bringing the counter-electrode closer, care should be taken that the system is in the correct position.
- d) When carrying out electrostatic experiments in which the leaf deflects frequently, the counter-electrode must be earthed, otherwise it will become charged with the result that the newly charged leaf will be drawn to the case and not to the counter-electrode, thus being exposed to damage.
- e) In all experiments where a productive source of current is used (power unit, capacitor), a safety resistor of at least 100 k Ω must at all costs be connected up in series. Well suited for this purpose is the safety resistor 100 M Ω (536 22).

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5. Installation of a spare system

The installation of a spare system is by no means a difficult process. Withdraw the counter-electrode, slightly turn the insulator with the system and remove it from the case, unscrew the top pin (Fig. 8), and remove the support of the damaged system. Only after having completed these preparations, take the spare system out of the box, and then screw it tightly into the boring of the insulator. When placing it into position, take care not to knock it against the casing because the protruding quartz loop is very sensitive. The system should be projected and, as described above, located in the correct position. The support of the damaged system should be screwed tightly in the case provided, and can be sent either to our works at Cologne or to one of our agents to have a new leaf fitted.

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6. Experiments

A number of supplementary apparatus are supplied with the Wulf electro-scope thus enabling a large number of experiments to be carried out, for example:

Zinc plate (546 31) and the pair of spiral electrodes (546 29) for demonstrating the photo-electric effect DC 535.215.1; a

and for control of a non-self-sustained discharge in air by means of a grid DC 537.545.2; b

Glass tube and plastic rod (541 01/04) for demonstrating the properties characteristic of electric charge DC 537.221; a

Capacitor plate (546 14) for determination of the saturation current and the dose of x-rays DC 537.531.9; a, DC 537.531; b and c

Test capacitor (546 23) for determination of the capacity of the Wulf electro-scope DC 537.723.5; c

The Wulf electro-scope as a current-measuring instrument DC 537.721; a

Ionization chamber (546 25) for measuring the half-life period of thorium emanation DC 539.163.1; a

Observation of ionization by collection in an ionization chamber DC 537.563; a, 539.16.08; k

for experiments on the extinction of α - and β -rays DC 539.163.2; b

for demonstration of ionization by means of radioactive substances DC 539.164.92; a

and for determination of the saturation current DC 539.164.92; c

Telescopic cylinder (546 27) and plate with radioactive preparation (546 44) for measuring the range of α -rays DC 539.163.2; a

Geiger counter (546 28) for counting α - and β -particles DC 539.16.08; a

Geiger-Müller counter for β - and γ -rays (559 10) with counter tube holder (559 11) for observation of β - and γ -rays DC 539.16.08; b

Determination of the activity of a radium source DC 539.16.08; k

Derivation of a frequency distribution from observations of the background counting rate DC 539.162.088; a

The relation between number of counts and frequency distribution DC 539.162.088; b and others.

Furthermore we would draw your attention to the book "Das Wulf-Elektroskop" by Th. Wulf, which specifies approximately 90 experiments and which is extremely suited for detailed tuition (Cat. No. 546 46).

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