

Auto-indexing program

- uses as starting input the previously stored list of 25 or less non-indexed reflections
- investigates the crystal without any further input about the crystal
- determines correct indices for these reflections
- calculates orientation matrix of the crystal
- does least-squares refinement of the matrix components
- calculates metrical tensor (Niggli matrix)
- determines normalized reduced unit cell
- prints pairs of perpendicular directions in direct lattice.

Matrix transformation program

- transforms unit cell and transformation matrix
- transforms orientation matrix
- transforms indices in reflection list on request
- saves transformation on specific file; old unit cell recalled by single command.

Reflection centering programs

- work with an input list of reflection positions, up to 25 reflections are possible
- use program controlled slits for fast running centering program
- produce an accuracy of centering which is better than 0.005° (this program is incorporated in the reflection search program and automatically centers the found reflections)
- prepare a new list of optimized reflection positions which are stored at same location as input list
- use the gravity center for optimization. This avoids $K\alpha_1/K\alpha_2$ split up problem
- contain special program provided to do additional, very accurate determination of 2θ value based on measurements at $+$ and -2θ
- contain special program provided to do scan around any axis or combination of scans with print out of reflection profile and information about reflection centering
- achieve very accurate determination of cell constants and check alignment of the goniometer.

Orientation Matrix programs

- provide different starting input possibilities
 - calculate orientation matrix using unit cell parameters and 2 reflection positions

```
"
LS
 1 -1.001  2.000 -2.997  0.0361
 2 -1.001  1.997 -3.998  0.0337
 3 -3.000  1.998 -2.000  0.0239
 4 -1.002  0.999 -4.002  0.0320
 5 -3.001  1.001 -1.999  0.0300
 6 -2.001  0.998 -2.999  0.0377
 7  5.998 -0.002 -0.998  0.0225
 8  0.002  2.006  4.001  0.0700
 9 -0.001  6.000 -0.002  0.0198
10 -0.002 -6.002  0.001  0.0242
11 -0.003  0.998  6.003  0.0328
12  0.003 -1.000 -5.997  0.0271
13  0.000  1.994  8.002  0.0414
```

RECIPROCAL AXIS MATRIX

```
 0.082611  0.023330 -0.007429
-0.052394 -0.035914 -0.091665
-0.000895  0.077627 -0.001481
```

NIGGLI-VALUES

```
148.1774  163.4779  157.8939
-33.7179 -55.9565 -49.1129
VOLUME= 1605.6960
```

CELL PARAMETERS

```
12.1728  12.7858  12.5656
102.1149 111.4587 108.3945
-0.209869 -0.365829 -0.315555
```

DIRECT AXIS MATRIX

```
11.510980  0.007187 -6.582241
-0.870474 -0.207156 -10.330611
-3.862276 12.784171 -2.801172
```

SIGMA DIRECT AXIS MATRIX

```
0.003272  0.001364  0.002561
0.005523  0.002303  0.004323
0.003876  0.001616  0.003034
```

SIGMA CELL PARAMETERS

```
0.0034  0.0035  0.0021
0.0199  0.0206  0.0318
0.000339 0.000335 0.000527
```