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Electron beam deflection tube

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The electron beam deflection tube allows investigations on the behavior of cathode rays in both electric and magnetic fields. Experiments with electrostatic or magnetic deflection of the electrons yield values corresponding to the order of magnitude of the specific electron charge e/m and the electron velocity v respectively.

The electrostatic deflection of the cathode rays is performed by means of two capacitor plates mounted in the tube. The (virtually) homogeneous magnetic field of the pair of Helmholtz coils (555 06) is used for the magnetic deflection.

Literature:

Physics Experiments, Volume 2, "Electricity-Electronics" (599 932).

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1 Safety Notes

- Plug pins (a_2) and (b_1) - cf. Fig. 1 - are to be covered safe to touch by means of two-way plug adapters (501 641) (high tension).
- Do not increase the filament voltage above 6.5 V (danger of damage to the heating coil!).
- Do not move the tube when the filament is incandescent.
- Strong mechanical stress exerted on the plastics caps glued over the tube wall, by pressure, pulling or shock should be avoided.
- The plug pins of the tube should only be connected using **one** cable as several superposed plugs might cause an inadmissible mechanical stress.
- Use a high-voltage proof heating voltage source in the experiment 3.2.1.

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2 Technical data; description

Directly heated tungsten cathode	
filament voltage U_f	6 V AC/DC
filament current I_f	approx. 1.35 A
anode voltage U_A	1 to 5 kV DC
distance between capacitor plates	approx. 5.4 cm
diameter of glass bulb	approx. 13 cm
total length of tube	approx. 30 cm


Fig. 1

Elektronenstrahl-Ablenkröhre (555 12) in Experimentierständer (555 05)
Electron beam deflection tube (555 12) mounted on the universal stand (555 05)

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Operational parts (Fig. 1)

- (a) Electron gun, consisting of directly heated tungsten incandescent cathode and cylinder-shaped anode
- (a₁) Pair of sockets connected to the incandescent cathode
- (a₂) Plug pin connected to the anode
- (b) Capacitor plates for electrostatic deflection
- (b₁) Plug pins connected to the capacitor plates
- (c) Fluorescent screen with graticule divided in centimeters

3 Operation

3.1 Operating apparatus:

Filament voltage: 6 V; 1,5 A, e. g.
tapped from transformer 6 V, 12 W 562 73
or
from high-voltage proof voltage source, e. g.
demountable transformer, consisting of
coil with 500 turns 562 21
coil with 50 turns 562 18
U-core with yoke and clamping device 562 11/12

Important:

A high-voltage proof heater voltage source is necessary for the experiments 3.2.1 (cathode potential of 1 to 5 kV to earth.)

Anode voltage:

1 to 5 kV DC; 2 mA; smoothed, continuously adjustable, e. g. tapped from the high-voltage power supply, 10 kV 522 37

Voltage for the capacitor plates:

1 to 5 kV DC; smoothed; adjustable,
e. g. tapped from high-voltage power supply 522 37

Voltage for pair of Helmholtz coils 555 06:

0 to 10 V DC; 1 A, smoothed, continuously adjustable,
e. g. tapped from
DC power supply, 0 to ± 15 V, 1.2 A 522 30

3.2 Examples of experiments

3.2.1 Electrostatic deflection

An electron passing through the (homogeneous) electric field of a plate capacitor (capacitor voltage U_p , distance between the plates d) with the velocity v moves on a parabolic trajectory (cf. Fig. 2).


Fig. 2

$$y = \frac{1}{2} \frac{e}{m} \frac{E}{v^2} x^2 \text{ mit } E = \frac{U_p}{d} \quad (\text{I})$$

$$y = \frac{1}{2} \frac{e}{m} \frac{E}{v^2} x^2 \text{ with } E = \frac{U_p}{d} \quad (\text{I})$$

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By means of an experiment assembly according to Fig 3. it is possible

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1. to confirm the relation (II) (approximately). For this purpose related pairs of values of x and y are plotted graphically with x^2 as abscissa and y as ordinate and the resulting curve confirms the proportionality $y \sim x^2$ from (II);
2. to determine v from (I) when x and y are found experimentally and $\frac{e}{m}$ is known.


Fig. 3

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3.2.2 Electromagnetic deflection; measuring example for determining e/m and v

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An electron of mass m and charge e moving perpendicularly to a magnetic field $\mu_0 H$ receives a centripetal force $\mu_0 H e v$ forcing it into a circular path. The centrifugal force $(mv^2)/r$ is equivalent to this centripetal force:

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For the velocity v of the electrons which is determined by the anode voltage U_A the following equation holds true:

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From (I) and (II) it follows for the specific charge e/m of the electron:

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
U_A can be measured directly; r and H are determined from the experimentally found data and the geometrical dimensions of the tube and the pair of Helmholtz coils in the following way:

Determination of r (see Fig. 4)

The radius of curvature r of the electron beam visible on the fluorescent screen follows from the relation that can be taken from Fig. 4.

$$r^2 = x^2 + (r - y)^2.$$

It follows therefrom

$$y = \frac{x^2 + y^2}{2r} \quad (\text{IIIa})$$

by way of calculation of r . Pairs of values of x and y are plotted graphically with $x^2 + y^2$ as abscissa and y (see Fig. 4) as ordinate. The proportionality constant $1/2r$ is determined from the straight line best fitting the points.

Fig. 4

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Determination of $B = \mu_0 \cdot H$

The magnetic flux density B is either determined experimentally (e.g. with the tangential B-probe, 516 60, and Teslameter, 562 62) or calculated according to Biot-Savart's law. According to this, the following applies for the field strength H of the (virtually) homogeneous magnetic field of the pair of Helmholtz coils:

$$H = \frac{n R^2 I}{(R^2 + a^2)^{3/2}} \quad (\text{IIIb})$$

n = number of turns per coil
 R = radius of the coil
 a = half mean distance of the coils
 I = current per coil

For the pair of Helmholtz coils used here, with

$n = 320$ (given),
 $R = 6.8$ cm (measured) and
 $a = 3.4$ cm (measured)

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Assemble arrangement according to Fig. 5;
connect the coils in series so that the direction of current flow is the same in both coils.


Fig. 5

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Graphical determination of r according to (IIIa) from several pairs of coordinates x, y (see Fig. 6 and table 1)

$$y = \frac{1}{2r} [x^2 + y^2]$$

x	cm	4	5	6	7	8
y	cm	0.4	0.7	1.0	1.4	1.8
$x^2 + y^2$	cm ²	16.2	25.5	37.0	51.0	67.2
$\tan \alpha = \frac{1}{2r}$	$\frac{1}{\text{cm}}$	0.028 (from diagram 6)				
r	cm	17,9				


Fig. 6

For e/m it follows with

$U_A = 4000 \text{ V}$,

$H = 10,5 \cdot 10^2 \text{ Am}^{-1}$

$r = 0.179 \text{ m}$

$\mu_0 = 1,25 \cdot 10^{-6} \text{ Vs A}^{-1} \text{ m}^{-1}$

$$\frac{e}{m} = \frac{2 \cdot 4000}{(1,25 \cdot 10^{-6} \cdot 10,5 \cdot 10^2 \cdot 0,179)^2} = 1,45 \cdot 10^{11} \text{ As kg}^{-1}$$

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$U_A = 4000 \text{ V}$,

$H = 10,5 \cdot 10^2 \text{ Am}^{-1}$,

$r = 0,179 \text{ m}$,

$\mu_0 = 1,25 \cdot 10^{-6} \text{ Vs A}^{-1} \text{ m}^{-1}$

$$\frac{e}{m} = \frac{2 \cdot 4000}{(1,25 \cdot 10^{-6} \cdot 10,5 \cdot 10^2 \cdot 0,179)^2} = 1,45 \cdot 10^{11} \text{ As kg}^{-1}$$