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High tension power supply

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1. Description

The high tension power supply meets the requirements for a reliable source of variable high-tension D.C. for carrying out electrostatic experiments, for operating ionization chambers, Geiger counters, Geiger-Müller counters and gas discharge tubes, for post-acceleration in the school cathode ray oscilloscope (575 15), and for other uses. If not stated otherwise, it is operated from the A. C. mains, 220 volts, 50 cycles *), and supplies two equal series-connected D. C. voltages which at the same time can be regulated continuously from 0 to approximately 3000 volts each with a rotary potentiometer. Built-in safety resistors of approx. 1 MΩ which are in series with the sockets + and — limit the D. C. so that the high tension side can be short-circuited without

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*) The instrument may be obtained also for 110, 130, 150 or 240 volts.

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Notes

- The five-figure numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.

- The data DC ... refer to the experiment descriptions published as a collection of leaflets in the "LEYBOLD PHYSICS LEAFLETS". Published to date:

	Cat. No.
1st series (72 sheets)	599 21
2nd series (52 sheets)	599 22
3rd series (49 sheets)	599 23
4th series (50 sheets)	599 24
(translation prepared)	599 24

- The specifications and illustrations are not binding in every detail for the design of the apparatus. It is our policy always to keep our manufacturing programme right up to date so that it makes full allowance for the latest knowledge acquired in all scientific and technical fields.

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causing any damage. A large band-pass filter renders the D.C. voltage practically free from humming.

The high tension power supply is located in a casing measuring (170 × 220 × 200) mm fitted at the top with a carrier handle. The front plate is fitted with the three high tension sockets, the rotary knob of the potentiometer, one earth socket, a control lamp and the mains switch. The symbols beside the mains switch are:

Dot inside circle:
The apparatus is switched on.

Dot outside circle:
The apparatus is switched off.

The rear of the high tension power supply is fitted with a DIN 49 490 two-pin apparatus plug for connecting a mains supply cable. The cable with "SCHUKO" safety plug is supplied with the apparatus. The high tension itself has no definite potential against earth. Anyone of the three high tension sockets can be earthed.

Undesirable mains voltage variations can be eliminated by using the magnetic voltage stabilizer (522 60) in front of the high tension power supply. The latter will then deliver an essentially constant high tension.

2. Technical data

The basic circuit diagram of the high tension power supply is illustrated in Fig. 2.

Outputs:

No-load voltage 0 ... 2 × approx. 3,000 volts.

Short-circuit current approx. 2.5 milliamps.

The table below shows the voltage delivered between sockets 0 and + (or —) when the potentiometer is turned up fully, as a function of the current drawn by the load. Besides, the respective ripple and noise is indicated as a proportion of the hum amplitude to the terminal voltage.

Current	Terminal voltage	Ripple + noise
	approx.	
0 mA	3,000 V	< 0.05 %
0.5 mA	2,000 V	approx. 0.1 %
1 mA	1,300 V	approx. 0.5 %
1.5 mA	600 V	approx. 1 %

Fuse:

Fuse, rated value of 500 milliamps., 20 mm long, in rear plate of casing.

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3. Note

The casing, which actually is free of voltage, must be earthed in a reliable manner to protect whoever is carrying out the experiment. This is either possible by means of the third strand of the mains cable (e. g. safety contact of the SCHUKO plug) or by means of the earth socket on the front plate of the casing.

Owing to the large band-pass filter of the unit, the high tension sockets still carry high voltage even after the unit has been switched off. This high voltage does not fade away until several minutes have elapsed. If necessary the discharge should be accelerated by connecting the three high tension sockets in series. Owing to the installation of the large band-pass filter, the high voltage only decreases gradually when the rotary potentiometer is turned down.

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u. v. m.

4. Applications

Descriptions of experiments where the high tension power supply is used:

Demonstration of the photoelectric effect by means of the Wulf electroscope (DC 535.215.1; a, 3rd series, sheet 2).

Saturation current in the case of ionization by X-rays (DC 537.531.9; a, 3rd series, sheet 28).

Control of a non self-sustained discharge in air by means of a grid (DC 537.545.2; b, 3rd series, sheet 35).

Demonstration of α - and β -particles with point counter and electroscope (DC 539.16.08; a, 3rd series, sheet 38)

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