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Geiger Counter

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
Fig. 1

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1. Description

The counter is used in conjunction with the Wulf electroscope (546 01 A) to detect individual α -rays emanated from weakly radioactive preparations.

The apparatus is designed as a metal cylinder, approximately 5 cm long, with one side open. The other side is sealed by means of an insulator. The insulator is fitted with a fine-pointed needle connected to a 4 mm socket. For assembling, the cylinder is equipped with a short rod and a plug. Other accessories to the apparatus are a needle with a very weak radium preparation, stored in a small glass tube, and a small book of aluminium foils used to seal the opening of the counter.

2. Assembly

The leaf in the system of the Wulf electroscope serves to indicate the α -particles. It should be adjusted in

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Notes

- The five-figure numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.
 - The data DC . . . refer to the experiment descriptions published as a collection of leaflets in the "LEYBOLD PHYSICS LEAFLETS".
Published to date:
- | | Cat. No. |
|----------------------------------|----------|
| 1st series (72 sheets) | 599 21 |
| 2nd series (52 sheets) | 599 22 |
| 3rd series (49 sheets) | 599 23 |
3. The data and illustrations are not binding in every detail for the design of the apparatus. It is our sole aim always to adapt our manufacturing programme to the most recent knowledge gained in all scientific and technical fields.

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such a manner that the glow discharge emitted by a single α -particle either brings about a perceptible enlargement of the deflection or causes the leaf to strike against the counter electrode. Among the many circuits possible, the one illustrated in Fig. 2 has proved to be particularly efficient.

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N is the high-tension power supply (522 37). The high-tension capacitor (538 82) charged to approximately 3000 volts by means of either the Van de Graaff generator (541 70) or the demonstration transformer (562 11/12/14/16), with the rectifier valve RFG 5 (586 40 and 586 41) can also be used instead. V is the high-tension voltmeter 4 kV (204 39), W the safety resistor, e. g. the resistance coil 100 Mohms (536 22), and Z the counter. The voltage must be adjusted to exactly the right value. It should be increased until the leaf of the electroscope, set at sensitivity, deflects constantly, even without radioactive preparation. Then it should be reduced with the voltage until the leaf does not deflect more than twice a minute. After the voltage has been regulated, the preparation for the counter should be placed in front of the counter's opening. It is so weak that the individual α -particles can be indicated separately, thus rendering it possible to observe the irregular temporal sequence of emission. When the electroscope is fully sensitive, each particle causes the leaf to strike, and when the influence plate is withdrawn a little the deflection is enlarged perceptibly.

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3. Special Hints

The tip of the counter must be cleaned very carefully to remove the finest of fibres and dust. If this is not done, the

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voltage set according to the indications will be too low so that the counter does not function. The best way of cleaning the tip is to place it in the mouth and to draw it between the lips. Furthermore, the insulator of the counter must be protected against dust, and should also be cleaned.

To prevent the prepared counter from becoming contaminated again, one of the supplied aluminium foils can be placed over the opening of the counter, this being effected by coating the periphery of the opening with a little grease and pressing it against the foil. When sealed in this manner, the well-cleaned counter remains ready for use for a few weeks. The foils, which are only a very few μ thick, do not absorb the α -rays to any great extent. In order to protect it, the preparation should only be stored in the small glass tube provided for this purpose.

4. Experiments

Demonstration of α - and β -particles with point counter and electroscope DC 539.16.08; a, 3rd series, sheet 38

Detailed descriptions of experiments can also be found in the book by Th. Wulf: „Das Wulf-Elektroskop“ (in German, 546 46).

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