




Surveyor UV/Vis Plus Detector

Hardware Manual

60053-97108 Revision D January 2009

© 2009 Thermo Fisher Scientific Inc. All rights reserved.

Surveyor and ChromQuest are registered trademarks and LightPipe is a trademark of Thermo Fisher Scientific Inc. in the United States.

Kel-F is a registered trademark of 3M Company.

The following are registered trademarks in the United States and other countries: Microsoft and Windows are registered trademarks of Microsoft Corporation. Delrin, Teflon, and Tefzel are registered trademarks of E. I. du Pont de Nemours and Company.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Thermo Fisher Scientific Inc. provides this document to its customers with a product purchase to use in the product operation. This document is copyright protected and any reproduction of the whole or any part of this document is strictly prohibited, except with the written authorization of Thermo Fisher Scientific Inc.

The contents of this document are subject to change without notice. All technical information in this document is for reference purposes only. System configurations and specifications in this document supersede all previous information received by the purchaser.

Thermo Fisher Scientific Inc. makes no representations that this document is complete, accurate or error-free and assumes no responsibility and will not be liable for any errors, omissions, damage or loss that might result from any use of this document, even if the information in the document is followed properly.

This document is not part of any sales contract between Thermo Fisher Scientific Inc. and a purchaser. This document shall in no way govern or modify any Terms and Conditions of Sale, which Terms and Conditions of Sale shall govern all conflicting information between the two documents.

Release history: Revision A October 2005, Revision B April 2006, Revision C March 2008, Revision D January 2009.

For Research Use Only. Not regulated for medical or veterinary diagnostic use by U.S. Federal Drug Administration or other competent authorities.



Regulatory Compliance

Thermo Fisher Scientific performs complete testing and evaluation of its products to ensure full compliance with applicable domestic and international regulations. When the system is delivered to you, it meets all pertinent electromagnetic compatibility (EMC) and safety standards as described in the next section by product name.

Changes that you make to your system might void compliance with one or more of these EMC and safety standards. Changes to your system include replacing a part or adding components, options, or peripherals not specifically authorized and qualified by Thermo Fisher Scientific. To ensure continued compliance with EMC and safety standards, replacement parts and additional components, options, and peripherals must be ordered from Thermo Fisher Scientific or one of its authorized representatives.

Surveyor UV/Vis Plus Detector

EMC Directive 89/336/EEC amended by 92/31/EEC and 93/68/EEC

EMC compliance has been evaluated by Underwriters Laboratories Inc.

EN 55011	1998	EN 61000-4-3	2002
EN 61000-3-2	1995, A1; 1998, A2; 1998, A14; 2000	EN 61000-4-4	1995, A1; 2001, A2; 2001
IEC 61000-3-2	2000	EN 61000-4-5	1995, A1; 2001
EN 61000-3-3	1995	EN 61000-4-6	1996, A1; 2001
IEC 61000-3-3	1994	EN 61000-4-11	1994, A1; 2001
EN 61326-1	1997	ENV 50204	1995
EN 61000-4-2	1995 A1; 1998 A2; 2001	CISPR 11	1999, A1; 1999, A2; 2002
FCC Class A, CFR 47 Part 15 Subpart B: 2003			

Low Voltage Safety Compliance

This device complies with Low Voltage Directive EN 61010-1:2001.

FCC Compliance Statement

THIS DEVICE COMPLIES WITH PART 15 OF THE FCC RULES. OPERATION IS SUBJECT TO THE FOLLOWING TWO CONDITIONS: (1) THIS DEVICE MAY NOT CAUSE HARMFUL INTERFERENCE, AND (2) THIS DEVICE MUST ACCEPT ANY INTERFERENCE RECEIVED, INCLUDING INTERFERENCE THAT MAY CAUSE UNDESIRE OPERATION.
--



CAUTION Read and understand the various precautionary notes, signs, and symbols contained inside this manual pertaining to the safe use and operation of this product before using the device.

Notice on Lifting and Handling of Thermo Scientific Instruments

For your safety, use a team effort to lift or move Thermo Fisher Scientific LC instruments that are too heavy or bulky for one person to lift or move. In compliance with international regulations, Thermo Fisher Scientific is notifying you that the Surveyor UV/Vis Detector weighs 19.5 kg (43 lbs).

Notice on the Proper Use of Thermo Scientific Instruments

In compliance with international regulations: Use of this instrument in a manner not specified by Thermo Fisher Scientific could impair any protection provided by the instrument.

Notice on the Susceptibility to Electromagnetic Transmissions

Your instrument is designed to work in a controlled electromagnetic environment. Do not use radio frequency transmitters, such as mobile phones, in close proximity to the instrument.

For manufacturing location, see the label on the instrument.

WEEE Compliance

This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2002/96/EC. It is marked with the following symbol:



Thermo Fisher Scientific has contracted with one or more recycling or disposal companies in each European Union (EU) Member State, and these companies should dispose of or recycle this product. See www.thermo.com/WEEERoHS for further information on Thermo Fisher Scientific's compliance with these Directives and the recyclers in your country.

WEEE Konformität

Dieses Produkt muss die EU Waste Electrical & Electronic Equipment (WEEE) Richtlinie 2002/96/EC erfüllen. Das Produkt ist durch folgendes Symbol gekennzeichnet:



Thermo Fisher Scientific hat Vereinbarungen mit Verwertungs-/Entsorgungsfirmen in allen EU-Mitgliedsstaaten getroffen, damit dieses Produkt durch diese Firmen wiederverwertet oder entsorgt werden kann. Mehr Information über die Einhaltung dieser Anweisungen durch Thermo Fisher Scientific, über die Verwerter, und weitere Hinweise, die nützlich sind, um die Produkte zu identifizieren, die unter diese RoHS Anweisung fallen, finden sie unter www.thermo.com/WEEERoHS.

Conformité DEEE

Ce produit doit être conforme à la directive européenne (2002/96/EC) des Déchets d'Équipements Électriques et Électroniques (DEEE). Il est marqué par le symbole suivant:



Thermo Fisher Scientific s'est associé avec une ou plusieurs compagnies de recyclage dans chaque état membre de l'union européenne et ce produit devrait être collecté ou recyclé par celles-ci. Davantage d'informations sur la conformité de Thermo Fisher Scientific à ces directives, les recycleurs dans votre pays et les informations sur les produits Thermo Fisher Scientific qui peuvent aider la détection des substances sujettes à la directive RoHS sont disponibles sur www.thermo.com/WEEERoHS.

CAUTION Symbol	CAUTION	VORSICHT	ATTENTION	PRECAUCION	AVVERTENZA
	Electric Shock: This instrument uses high voltages that can cause personal injury. Before servicing, shut down the instrument and disconnect the instrument from line power. Keep the top cover on while operating the instrument. Do not remove protective covers from PCBs.	Elektroschock: In diesem Gerät werden Hochspannungen verwendet, die Verletzungen verursachen können. Vor Wartungsarbeiten muß das Gerät abgeschaltet und vom Netz getrennt werden. Betreiben Sie Wartungsarbeiten nicht mit abgenommenem Deckel. Nehmen Sie die Schutzabdeckung von Leiterplatten nicht ab.	Choc électrique: L'instrument utilise des tensions capables d'infliger des blessures corporelles. L'instrument doit être arrêté et débranché de la source de courant avant tout intervention. Ne pas utiliser l'instrument sans son couvercle. Ne pas enlever les étuis protecteurs des cartes de circuits imprimés.	Descarga eléctrica: Este instrumento utiliza altas tensiones, capaces de producir lesiones personales. Antes de dar servicio de mantenimiento al instrumento, éste deberá apagarse y desconectarse de la línea de alimentación eléctrica. No opere el instrumento sin sus cubiertas exteriores quitadas. No remueva las cubiertas protectoras de las tarjetas de circuito impreso.	Shock da folgorazione. L'apparecchio è alimentato da corrente ad alta tensione che può provocare lesioni fisiche. Prima di effettuare qualsiasi intervento di manutenzione occorre spegnere ed isolare l'apparecchio dalla linea elettrica. Non attivare lo strumento senza lo schermo superiore. Non togliere i coperchi a protezione dalle schede di circuito stampato (PCB).
	Chemical: This instrument might contain hazardous chemicals. Wear gloves when handling toxic, carcinogenic, mutagenic, or corrosive or irritant chemicals. Use approved containers and proper procedures to dispose waste oil.	Chemikalien: Dieses Gerät kann gefährliche Chemikalien enthalten. Tragen Sie Schutzhandschuhe beim Umgang mit toxischen, karzinogenen, mutagenen oder ätzenden/reizenden Chemikalien. Entsorgen Sie verbrauchtes Öl entsprechend den Vorschriften in den vorgeschriebenen Behältern.	Chimique: Des produits chimiques dangereux peuvent se trouver dans l'instrument. Portez des gants pour manipuler tous produits chimiques toxiques, cancérogènes, mutagènes, ou corrosifs/irritants. Utiliser des récipients et des procédures homologuées pour se débarrasser des déchets d'huile.	Química: El instrumento puede contener productos químicos peligrosos. Utilice guantes al manejar productos químicos tóxicos, carcinógenos, mutágenos o corrosivos/irritantes. Utilice recipientes y procedimientos aprobados para deshacerse del aceite usado.	Prodotti chimici. Possibile presenza di sostanze chimiche pericolose nell'apparecchio. Indossare dei guanti per maneggiare prodotti chimici tossici, cancerogeni, mutageni, o corrosivi/irritanti. Utilizzare contenitori aprovo e seguire la procedura indicata per lo smaltimento dei residui di olio.
	Heat: Before servicing the instrument, allow any heated components to cool.	Hitze: Warten Sie erhitzte Komponenten erst nachdem diese sich abgekühlt haben.	Haute Temperature: Permettre aux composants chauffés de refroidir avant tout intervention.	Altas temperaturas: Permita que los componentes se enfríen, ante de efectuar servicio de mantenimiento.	Calore. Attendere che i componenti riscaldati si raffreddino prima di effettuare l'intervento di manutenzione.
	Fire: Use care when operating the system in the presence of flammable gases.	Feuer: Beachten Sie die einschlägigen Vorsichtsmaßnahmen, wenn Sie das System in Gegenwart von entzündbaren Gasen betreiben.	Incendie: Agir avec précaution lors de l'utilisation du système en présence de gaz inflammables.	Fuego: Tenga cuidado al operar el sistema en presencia de gases inflamables.	Incendio. Adottare le dovute precauzioni quando si usa il sistema in presenza di gas infiammabili.
	Eye Hazard: Eye damage could occur from splattered chemicals or flying particles. Wear safety glasses when handling chemicals or servicing the instrument.	Verletzungsgefahr der Augen: Verspritzte Chemikalien oder kleine Partikel können Augenverletzungen verursachen. Tragen Sie beim Umgang mit Chemikalien oder bei der Wartung des Gerätes eine Schutzbrille.	Danger pour les yeux: Des projections chimiques, liquides, ou solides peuvent être dangereuses pour les yeux. Porter des lunettes de protection lors de toute manipulation de produit chimique ou pour toute intervention sur l'instrument.	Peligro par los ojos: Las salicaduras de productos químicos o partículas que salten bruscamente pueden causar lesiones en los ojos. Utilice anteojos protectores al manipular productos químicos o al darle servicio de mantenimiento al instrumento.	Pericolo per la vista. Gli schizzi di prodotti chimici o delle particelle presenti nell'aria potrebbero causare danni alla vista. Indossare occhiali protettivi quando si maneggiano prodotti chimici o si effettuano interventi di manutenzione sull'apparecchio.
	General Hazard: A hazard is present that is not included in the above categories. Also, this symbol appears on the instrument to refer the user to instructions in this manual. When the safety of a procedure is questionable, contact your local Technical Support organization for Thermo Fisher Scientific San Jose Products.	Allgemeine Gefahr: Es besteht eine weitere Gefahr, die nicht in den vorstehenden Kategorien beschrieben ist. Dieses Symbol wird im Handbuch außerdem dazu verwendet, um den Benutzer auf Anweisungen hinzuweisen. Wenn Sie sich über die Sicherheit eines Verfahrens im unklaren sind, setzen Sie sich, bevor Sie fortfahren, mit Ihrer lokalen technischen Unterstützungsorganisation für Thermo Fisher Scientific San Jose Produkte in Verbindung.	Danger général: Indique la présence d'un risque n'appartenant pas aux catégories citées plus haut. Ce symbole figure également sur l'instrument pour renvoyer l'utilisateur aux instructions du présent manuel. Si la sûreté d'un procédé est incertaine, avant de continuer, contacter le plus proche Service Clientèle pour les produits de Thermo Fisher Scientific San Jose.	Peligro general: Significa que existe un peligro no incluido en las categorías anteriores. Este símbolo también se utiliza en el instrumento para referir al usuario a las instrucciones contenidas en este manual. Cuando la certidumbre acerca de un procedimiento sea dudosa, antes de proseguir, póngase en contacto con la Oficina de Asistencia Técnica local para los productos de Thermo Fisher Scientific San Jose.	Pericolo generico. Pericolo non compreso tra le precedenti categorie. Questo simbolo è utilizzato inoltre sull'apparecchio per segnalare all'utente di consultare le istruzioni descritte nel presente manuale. Quando è in dubbio la misura di sicurezza per una procedura, prima di continuare, si prega di mettersi in contatto con il Servizio di Assistenza Tecnica locale per i prodotti di Thermo Fisher Scientific San Jose.

CAUTION Symbol	CAUTION	危險警告	危險警告
	Electric Shock: This instrument uses high voltages that can cause personal injury. Before servicing, shut down the instrument and disconnect the instrument from line power. Keep the top cover on while operating the instrument. Do not remove protective covers from PCBs.	電撃: この計測器は高電圧を使用し、人体に危害を与える可能性があります。保守・修理は、必ず作業を停止し、電源を切ってから実施して下さい。上部カバーを外したままで計測器を使用しないで下さい。プリント配線板の保護カバーは外さないで下さい。	電撃: 儀器設備使用會造成人身傷害的高伏電壓。在維修之前，必須先關儀器設備並切除電源。務必要在頂蓋蓋上的情況下操作儀器。請勿拆除PCB保護蓋。
	Chemical: This instrument might contain hazardous chemicals. Wear gloves when handling toxic, carcinogenic, mutagenic, or corrosive or irritant chemicals. Use approved containers and proper procedures to dispose waste oil.	化学物質: 危険な化学物質が計測器中に存在している可能性があります。毒性、発がん性、突然変異性、腐食・刺激性などのある薬品を取り扱う際は、手袋を着用して下さい。廃油の処分には、規定の容器と手順を使用して下さい。	化學品: 儀器設備中可能存在有危險性的化學物品。接觸毒性致癌、誘變或腐蝕／刺激性化學品時，請配帶手套。處置廢油時，請使用經過許可的容器和程序。
	Heat: Before servicing the instrument, allow any heated components to cool.	熱: 熱くなった部品は冷えるのを待ってから保守・修理を行って下さい。	高溫: 請先等高溫零件冷卻之後再進行維修。
	Fire: Use care when operating the system in the presence of flammable gases.	火災: 可燃性のガスが存在する場所でシステムを操作する場合は、充分な注意を払って下さい。	火災: 在有易燃氣體的場地操作該系統時，請務必小心謹慎。
	Eye Hazard: Eye damage could occur from splattered chemicals or flying particles. Wear safety glasses when handling chemicals or servicing the instrument.	眼に対する危険: 化学物質や微粒子が飛散して眼を傷つける危険性があります。化学物質の取り扱い、あるいは計測器の保守・修理に際しては防護眼鏡を着用して下さい。	眼睛傷害危險: 飛濺的化學品或顆粒可能造成眼睛傷害。處理化學品或維儀器設備時請佩戴安全眼鏡。
	General Hazard: A hazard is present that is not included in the above categories. Also, this symbol appears on the instrument to refer the user to instructions in this manual.	一般的な危険: この標識は上記以外のタイプの危険が存在することを示します。また、計測器にこの標識がついている場合は、本マニュアル中の指示を参照して下さい。	一般性危險: 說明未包括在上述類別中的其他危險。此外，儀器設備上使用這個標誌，以指示用戶本使用手冊中的說明。
	When the safety of a procedure is questionable, contact your local Technical Support organization for Thermo Fisher Scientific San Jose Products.	安全を確保する手順がよくわからない時は、作業を一時中止し、お近くのサーモエレクトロンサンローゼプロダクトのテクニカルサポートセンターにご連絡ください。	如對安全程序有疑問，請在操作之前與當地的菲尼根技術服務中心聯繫。

Contents

	Preface	xi
	Related Documentation	xi
	Safety and Special Notices	xi
	Contacting Us	xii
Chapter 1	Introduction	1
	Functional Description	2
	LightPipe Flowcell	3
	Specifications	5
Chapter 2	Installation	7
	Installation Checklist	8
	Unpacking and Inspecting Instrument	8
	Making Initial Instrument Preparations	9
	Checking the Power Setting and Fuses	10
	Making the Initial Back Panel Connections	11
	Ethernet Connection	11
	Analog Output Connections	11
	Analog Output Settings	13
	Unit ID Settings	13
	Remote Communications Connections	13
	System Interconnect Cable Connections	15
	Connecting and Installing the Flowcell	19
	Powering On the Detector for the First Time	20
Chapter 3	Maintenance	21
	Recommendations for Routine Maintenance	22
	Tools Required	22
	Cleaning Detector External Surfaces	22
	LightPipe Flowcell	23
	Removing the LightPipe Flowcell	23
	Storing the LightPipe Flowcell	26
	Installing the LightPipe Mounting Assembly	26
	Installing the LightPipe Flowcell	27

	Standard Analytical Flowcell	29
	Removing the Standard Flowcell Assembly	29
	Installing the Standard Analytical Flowcell	31
	Cleaning Flowcells.	32
	Cleaning with Organic Solvents	32
	Cleaning with Nitric Acid	33
	Replacing the Lamps	34
	The Deuterium Lamp	35
	The Tungsten Lamp	37
Chapter 4	Troubleshooting.	41
	Detector-Related Problems	42
	Error Messages	47
Chapter 5	Accessories and Replaceable Parts	49
	Replaceable Parts	49
	Optional Flowcells for the UV/Vis Detector	50
Appendix A	Surveyor Firmware Upgrade Utility.	51
	Index	55

Preface

This *Surveyor UV/Vis Plus Detector Hardware Manual* provides you with information on how to set up and maintain your Surveyor UV/Vis Plus detector.

Related Documentation

In addition to this guide, Thermo Fisher Scientific provides the following documents for the Surveyor UV/Vis Plus Detector:

- *Surveyor Plus Preinstallation Requirements Guide*
- *Surveyor Plus Getting Connected Guide*
- *Surveyor Plus Getting Started with ChromQuest Guide*
- *Surveyor Plus Getting Started with Xcalibur Guide*

Safety and Special Notices

Make sure you follow the precautionary statements presented in this guide. The safety and other special notices appear in boxes.

Safety and special notices include the following:



CAUTION Highlights hazards to humans, property, or the environment. Each CAUTION notice is accompanied by an appropriate CAUTION symbol.

IMPORTANT Highlights information necessary to prevent damage to software, loss of data, or invalid test results; or might contain information that is critical for optimal performance of the system.

Note Highlights information of general interest.

Tip Highlights helpful information that can make a task easier.

Contacting Us

There are several ways to contact Thermo Fisher Scientific for the information you need.

❖ To contact Technical Support

Phone	800-532-4752
Fax	561-688-8736
E-mail	us.techsupport.analyze@thermofisher.com
Knowledge base	www.thermokb.com

Find software updates and utilities to download at mssupport.thermo.com.

❖ To contact Customer Service for ordering information

Phone	800-532-4752
Fax	561-688-8731
E-mail	us.customer-support.analyze@thermofisher.com
Web site	www.thermo.com/ms

❖ To copy manuals from the Internet

Go to mssupport.thermo.com and click **Customer Manuals** in the left margin of the window.

❖ To suggest changes to documentation or to Help

- Complete a brief survey about this document by clicking the link below. Thank you in advance for your help.



- Send an e-mail message to the Technical Publications Editor at techpubs-lcms@thermofisher.com.

Introduction

The Surveyor UV/Vis Plus Detector (see [Figure 1](#)) is a member of the Surveyor Plus™ family of LC instruments.

This chapter provides an introduction to the Surveyor UV/Vis Plus Detector.

Contents

- [Functional Description](#)
- [LightPipe Flowcell](#)
- [Specifications](#)

Figure 1. Surveyor UV/Vis Plus Detector




Functional Description

The Surveyor UV/Vis Plus Detector is a full-featured, time-programmable, variable-wavelength UV/Vis (ultraviolet / visible) absorbance detector. It operates in either the single or dual wavelength mode in the UV or visible ranges.

The UV/Vis detector is a bench top unit for inclusion into the Surveyor LC System. It is controlled remotely over an Ethernet communications link from a PC using the ChromQuest or Xcalibur data systems.

The UV/Vis detector consists of a dual light source, an optical bench (see [Figure 2](#)), dual photodiode detectors, and several printed circuit boards (PCBs).

Figure 2. The Surveyor UV/Vis Plus Detector optical bench



To provide a complete spectrum of ultraviolet and visible light, the detector uses a deuterium lamp for the UV range (190 to 365 nm) and a tungsten lamp for the visible range (366 to 800 nm). The lamps are protected by a cover with a special safety interlock to reduce the possibility of human exposure to harmful UV light. A computer-controlled shutter selects which lamp is used for sample analysis. The light from the selected lamp goes to a holographic grating that disperses the light into a spectrum. The grating can be turned to select which portion of the spectrum and which wavelength of light is used. The alignment of the grating is controlled by a stepper motor for precise wavelength control. Absorbance linearity is increased by a system of baffles that reduces stray light in the optical bench.

From the grating, the light travels to a fiber optic beamsplitter, where a portion of the light is redirected to a reference photodiode. The signal from the reference photodiode is used to compensate for changes in lamp intensity and helps to maintain consistent performance as the lamp performance varies. From the beamsplitter, most of the light passes through a lens and then on through the flowcell and onto the sample photodiode. The light falling on the sample photodiode is then processed by the analog PC board and converted into data that can be analyzed by the computer.

The available flowcell options include:

- Standard analytical flowcell with a 10 mm pathlength
- LightPipe flowcell with a 50 mm pathlength

Other flowcells can also be used in the detector. See [Table 3](#).

LEDs on the front of the module keep the user informed of power, communications, run, and lamp status.


LightPipe Flowcell

The Surveyor UV/Vis Plus Detector can be ordered with a 50 mm LightPipe flowcell.

The internal bore of the LightPipe flowcell is 50 mm long and has a volume of 10 μL . The 5 cm optical pathlength produces a signal five times that of a conventional 1 cm flowcell, while the low internal volume of 10 μL minimizes band broadening. A special low refractive index coating of the internal bore ensures a high optical throughput and minimizes short-term noise.

Mobile phase enters the LightPipe flowcell through a port in the bottom of the flowcell and exits through a port in the top of the flowcell. See [Figure 3](#). Directing the flow from the bottom to the top of the flowcell reduces the development of air bubbles.

Figure 3. Schematic of LightPipe flowcell, showing the internal flowpath



Sensitive optical fibers are exposed at both ends of the flowcell as shown in [Figure 4](#). Thermo Fisher Scientific ships the LightPipe flowcell with end caps to protect these optical fibers. See [Figure 5](#). Leave these protective end caps in place until you install the LightPipe flowcell. Replace the protective end caps for storage if you remove the LightPipe flowcell from the detector.



CAUTION Do *not* touch the ends of the LightPipe flowcell. Doing so could damage the sensitive optical fibers. If you must grasp the ends of the LightPipe flowcell, wear clean, talc-free gloves.

Figure 4. End of LightPipe flowcell, showing the sensitive optical fiber





Figure 5. LightPipe flowcell with protective end caps



Specifications

Light Source:	Deuterium and tungsten lamps, re-aligned
Wavelength Range:	190 nm to 800 nm • D₂ Lamp: 190 to 365 nm • W Lamp: 366 to 800 nm
Wavelength Accuracy:	±1 nm at 254 nm and 656 nm
Spectral Bandwidth:	6 nm
Absorbance Range:	0.0005 AUFS to +3.0 AUFS
Short Term Noise* :	• 6 µAU/cm at 254 nm, 1 s rise time, single wavelength mode (50 mm LightPipe flowcell) • 20 µAU/cm at 254 nm, 1 s rise time, single wavelength mode (10 mm standard analytical flowcell)
Drift*:	• ≤ 1 mAU/h at 254 nm and 546 nm (50 mm LightPipe flowcell) • ≤ 0.1 mAU/h at 254 nm and 546 nm (10 mm standard analytical flowcell)
Warm-up:	90 min to meet noise and drift specifications
Linearity:	deviation ≤ 5% up to 2.0 AU at 256 nm
Scan Rate:	4, 5, 6.67, 10, or 20 Hz (user selectable)
Rise Time:	0, 0.1, 0.2, 0.5, 1, 2, 3, 4, or 5 s (user selectable)
Cell Dimensions:	• 10 mm, 15 µL standard analytical • 50 mm, 10 µL LightPipe
Cell Pressure Rating:	1000 psi
Analog Outputs (2):	Unattenuated at 10 mV/AU, 100 mV/AU, or 1.0 V/AU
Remote Controls:	Start, Zero, RJ45 Ethernet interface for ChromQuest or Xcalibur
Dimensions:	18 cm (7.1 in.) × 36 cm (14.2 in.) × 47 cm (18.5 in.) (<i>h</i> × <i>w</i> × <i>d</i>)
Weight:	19.5 kg (43 lbs)
Power Requirements:	100/115 or 230 V ac; 50/60 Hz, 225 VA max.
Fuses:	• 2.5 A at 100/115 V ac (slo-blo) • 1.25 A at 230 V ac (slo-blo)
Operating Temperature:	+10 to +30 °C
Storage Temperature:	- 40 to +70 °C
Operating Humidity:	5% to 95% non-condensing relative humidity
Product Certification:	CE, TUV-C/US, FCC (EMI)

* According to ASTM E1657-94 "Standard Practice for Testing Variable-Wavelength Photometric Detectors Used in Liquid Chromatography" (for 50 mm flowcell, 5 nm bandwidth, 2 sec rise time, and MeOH at 1 mL/min.).

Installation

This chapter describes the initial installation of your Surveyor UV/Vis Plus Detector, including connection to other chromatographic instrumentation. The Installation Checklist, on the back of this page, is an abbreviated version of this chapter and can be used as a quick reference of how to conduct a successful installation. Make a copy of the checklist and fill it out when the installation is complete. Include the completed checklist in your maintenance records.

Note Perform the installation in the sequence presented on the Installation Checklist and detailed in this chapter.

Contents

- [Installation Checklist](#)
- [Unpacking and Inspecting Instrument](#)
- [Making Initial Instrument Preparations](#)
- [Checking the Power Setting and Fuses](#)
- [Making the Initial Back Panel Connections](#)
- [Connecting and Installing the Flowcell](#)
- [Powering On the Detector for the First Time](#)

Installation Checklist

The following installation checklist is a brief summary of the steps that you need to complete for the proper installation of your Surveyor UV/Vis Plus Detector:

- ☐ Unpack and inspect your instrument ([page 8](#)).
- ☐ Read the safety notices in the preface of this manual.
- ☐ Make the initial instrument preparations ([page 8](#)).
- ☐ Check the power settings and fuses ([page 10](#)).
- ☐ Make the initial back panel connections ([page 11](#)).
- ☐ Connect and install the flowcell ([page 19](#)).
- ☐ Power on the detector for the first time ([page 20](#)).
- ☐ Download the Surveyor UV/Vis Plus Detector firmware, if required ([Appendix A, “Surveyor Firmware Upgrade Utility.”](#))

This Surveyor UV/Vis Plus Detector was installed by:

(Name)

(Date)

Unpacking and Inspecting Instrument

Carefully remove the detector from the shipping container and inspect both the detector and the packaging for any signs of damage. If you find any damage, save the shipping materials and immediately contact the shipping company.

The shipping container should contain the following items:

- Surveyor UV/Vis Plus Detector
- Accessory kit

Detectors ordered with a standard 10 mm flowcell arrive with the flowcell installed in the detector. Detectors ordered with the 50 mm LightPipe are shipped with the LightPipe packaged in the accessory kit.

The accessory kit (P/N 60053-62004) contains the items listed in [Table 1](#). Carefully check to make sure you have received all the items listed on the packing list. If any items are missing, contact your Thermo Fisher Scientific representative immediately.

Table 1. Surveyor UV/Vis accessory kit

Item	Quantity	Part Number
Cable, Ext. Run/Auto Zero	1	6040-0102
8-pin cable connector, 3.81 mm pitch, minicombico	2	00004-02511
Cable, RJ45 Ethernet communication, CAT-5, shielded	1	70111-63302
Cable, analog SS DET	2	6040-0103
Nut, 1/16 in., 1/4-28	2	3256-0024
Ferrule, 7/4-28 ETFE	2	3256-0025
Nut and Ferrule set tube	1	3256-0022
LC Test Mix, 1 mL	1	A3469-020
Funnel, Surveyor	1	F3012-010
Tubing, convoluted, detector drain	1	F5034-010
Screwdriver, Hex, ball end, 3/32-in, black	1	00725-00040
Screwdriver, Hex, ball end, 1/8-in, black	1	00725-00042
Screwdriver, Hex, ball end, 9/64-in, black	1	00725-00044
Surveyor UV/Vis Plus Detector Hardware manual	1	60053-97108

Making Initial Instrument Preparations

Place the detector on a bench top as close as possible to the chromatographic column outlet (minimizing the length of tubing necessary for connection to the flowcell inlet). If you are using other Surveyor modules, place the UV/Vis detector on top of the Surveyor Autosampler Plus. Be sure to place the HPLC system in a draft-free location away from an open window, air conditioner vents, or other circulating air source. A stable room temperature (± 0.5 °C) is necessary for applications requiring maximum detection sensitivity. Allow at least 15 cm (6 in.) of clear space between the back panel of the detector and any wall or obstruction. This provides access to the back-panel connectors and a free flow of cooling air.


You will need the following tools for installation:

- narrow-tip screwdriver (2 mm wide)
- #2 Phillips screwdriver
- 1/4 in. open end wrench
- 1/8 in. hex wrench (supplied)
- 3/32 in. hex wrench (supplied)
- 9/64 in. hex wrench (supplied)

Checking the Power Setting and Fuses

The detector is shipped with the voltage and fuses preset for your location. To verify the power setting, look through the cut-out window on the power entry module located at the lower right on the detector back panel. See [Figure 6](#).

Figure 6. Power entry module showing selection of the voltage



CAUTION Do not connect the power cord to your instrument without first verifying that the voltage is properly set for your location! Never connect the detector to, or operate the detector with, an electrical line source with power fluctuations greater than 10% above or below the nominal rated line voltage!

❖ To check the fuses and change the voltage setting, if necessary

1. Place the tip of the narrow-blade screwdriver in the small slot to the left of the fuse holder in the power entry module and push to the right. The fuse holder should pop out.
2. Pull the fuse holder out of the power entry module.
3. Verify that the fuse holder contains the appropriate size fuses: T2.5 A for 100/115 V, or T1.25 A for 230 V.
4. If necessary, slide out the voltage selection PC board and re-install it so that the desired voltage label is upright and readable. This voltage must agree with the incoming line voltage.
5. Slide the fuse holder back into the power entry module until it snaps in.

IMPORTANT Do not connect the power cord or turn on the detector yet. Wait until you perform the procedure in [‘Powering On the Detector for the First Time’](#) on [page 20](#) to connect the power cord.

Making the Initial Back Panel Connections

Use the cables provided in the installation kit to make the connections to the back panel of the detector, as described in this section. [Figure 7](#) shows the back panel of the detector.

This section contains the following topics:

- [Ethernet Connection](#)
- [Analog Output Connections](#)
- [Analog Output Settings](#)
- [Unit ID Settings](#)
- [Remote Communications Connections](#)
- [System Interconnect Cable Connections](#)

Ethernet Connection


Connect your Ethernet switch to the Ethernet connector outlet of the detector using the supplied shielded Ethernet cable (P/N 70111-63302).

Analog Output Connections

Two analog signal cables (twin axial computer cables) are provided in the installation kit to use for connecting the analog outputs from the UV/Vis detector to other data collection devices.

The analog signal cables have three wires protruding from the ends of the shielded cable. Two of these wires are electrically insulated. They are used for carrying an analog signal to data collection devices. Typically, the wire with the clear insulation is connected to the positive analog output, and the wire with the black insulation is connected to the signal ground (sometimes referred to as the negative signal).

Figure 7. UV/Vis detector back panel, showing the location of the connectors and power entry module




The third wire is not insulated. This bare wire is used to ground the cable shielding. The cable shielding reduces signal noise caused by radio frequency interference. The cable shielding is most effective if the bare wire is grounded at just one end.

At the detector back panel there are three different analog voltage outputs per channel: CHA – 10 mV, CHA – 100 mV, CHA – 1 V, CHB – 10 mV, CHB – 100 mV, CHB – 1 V; and a single ground per channel: GND A and GND B.

Each of the analog signal wire ends in the signal cables have a 1/4-in. of bare, tinned wire. The ends are stripped and soldered to allow electrical contact and to prevent fraying. To make each electrical connection, do the following: Insert the end of the wire into the appropriate terminal in the 8-pin terminal connector. Hold the wire in place while you tighten the small terminal set screw firmly onto the wire. Insert the terminal connector into the contact pins numbered 9 through 16 at the back of the instrument. See [Figure 8](#).

IMPORTANT Do not connect the detector ground terminals to any earth ground on your data system computer. This incorrect connection will lead to an increased noise level and a subsequent decrease in sensitivity.

Figure 8. UV/Vis detector outputs full-scale voltage connections



Analog Output Settings

The analog outputs are controlled by the wavelength, bandwidth, rise time, and zero functions of the detector. These outputs are compatible with data collection systems using any of the three different voltages (10 mV, 100 mV, or 1 V) by selecting the appropriate terminal of the analog output terminal connector (Figure 8).

Unit ID Settings

The Surveyor UV/Vis Plus Detector is shipped with the unit ID preset to 01 using the two rotary switches located on the back panel. The unit ID must match the Stack ID setting specified for the detector within the chromatography data system used to control it. You specify the Stack ID for the detector when you configure it for communication with your chromatography data system. The range of values for the unit ID is 01 to 99. The value of 00 is reserved for special service functions.

Remote Communications Connections

The Surveyor UV/Vis Plus Detector has the remote communications inputs RUN and ZERO, and the remote communications output EVENT.

RUN

The RUN input receives an inject signal from the autosampler, or from a manual injector, and sends a run-start signal to the PC-based software.

During installation, you must connect this UV/Vis detector RUN input to your autosampler or injector output. If your LC system contains a Surveyor Autosampler, you install the system synchronization wiring harness (system interconnect cable) to accomplish this task.

For Thermo Scientific autosamplers (other than the Surveyor Autosampler), connect the UV/Vis detector RUN input to the autosampler Inject Out output, and connect the UV/Vis detector GND terminal to the autosampler ground terminal.

For other autosamplers or injectors, first determine whether the device uses a TTL signal or a relay (a contact-closure), to remotely trigger detectors. Refer to the autosampler or injector documentation to determine whether a TTL signal or relay is used. Then wire the device to the UV/VIS detector, by doing one of the following:

- If the autosampler or injector uses a TTL signal to trigger the run, connect the TTL trigger to the detector RUN terminal, and connect the signal ground terminal (or TTL return terminal) to the detector GND terminal.
- If the autosampler or injector device output relay consists of a single, normally-open terminal and a common ground terminal, connect the normally-open terminal to the detector RUN terminal, and connect the common ground terminal to the detector GND terminal. If both of the relay contacts are electrically isolated from each other and from ground (the device simply closes a contact between the two relay terminals), then connect either one of the relay terminals to the detector RUN terminal, and connect the other relay terminal to the detector GND terminal.

Note The detector ground (GND) terminals (pins 1 and 7) on the back of the detector are tied to a single digital ground. You can use either ground terminal for digital input or output return connections.

ZERO

The UV/Vis detector ZERO connector on the back panel of the detector is used to zero the detector signal output from a remote device (generally at the start or end of each sample run).

The UV/VIS detector can be zeroed remotely with either a TTL low signal or with a contact closure.

During installation, you must connect this UV/Vis detector ZERO input to your remote device output. When the Surveyor UV/Vis Plus Detector is used with the Surveyor Autosampler, use the system interconnect cable to make this connection.

First determine whether the external device uses a TTL signal or a relay to remotely trigger detectors. Refer to the injector documentation of the external device to determine whether a TTL signal or relay is used. Then wire the device to the UV/VIS detector, using the appropriate procedure:

- If the external device uses a TTL signal to zero the detector, connect the TTL trigger to the detector ZERO terminal, and connect the signal ground terminal (or TTL return terminal) to the detector GND terminal.
- If the external device uses a relay to zero the detector, connect the normally-open terminal to the detector ZERO terminal, and connect the common ground terminal to the detector GND terminal. If the relay contacts are electrically isolated from each other (and from ground), then connect either one of the relay terminals to the detector ZERO terminal, and connect the other relay terminal to the detector GND terminal.

EVENT

The UV/Vis detector EVENT output connector is used to trigger an external device such as a fraction collector.

If the external device is triggered by a contact closure, connect the UV/Vis detector EVENT terminal (pin 8) to the positive pin on the external device input, and connect one of the UV/Vis detector GND terminals (pin 1 or 7) to the external device negative pin.

Note The external device Input terminal might not have markings indicating positive and negative polarity. In this case, connect the UV/VIS EVENT terminal to one of the pins, and connect a UV/VIS GND terminal to the other pin.

If the external device is triggered using a TTL signal, connect the +5 V Output (pin 2) to the positive Input terminal of the external device, and connect the UV/Vis detector EVENT Output (pin 8) to the negative Input terminal of the external device.

System Interconnect Cable Connections

The system interconnect cable coordinates the timing of the Surveyor plus modules during an injection sequence.

There are two versions of the system interconnect cable:

- The older version has five combicon connectors. The detector connector is identified by the DET tag on its adjacent cable as shown in [Figure 9](#).
- The newer version has seven combicon connectors, with three connectors labeled DETECTOR as shown in [Figure 10](#).

Figure 9. 5-connector version of the system interconnect cable





Figure 10. 7-connector version of the system interconnect cable




Connecting the 5-Connector System Interconnect Cable

❖ To connect the 5-connector interconnect cable

1. To connect the autosampler, plug the A/S connector into the left, 8-pin socket on the back panel of the autosampler. See [Figure 11](#).

Figure 11. Back panels of the Surveyor stack, showing the 5-connector interconnect cable




2. To connect the pump, do one of the following:
 - To connect the Surveyor LC Pump Plus, plug the connector identified by the blue sticker with the text LC PUMP to the 8-pin receptacle on the back panel of the pump. There is a blue label with the text LC PUMP above the receptacle.
 - To connect the Surveyor MS Pump Plus, plug the connector identified by the yellow sticker with the text MS PUMP to the 8-pin receptacle on the back panel of the pump. There is a yellow label with the text MS PUMP above the receptacle.
3. To connect a Surveyor UV/Vis Plus Detector, plug the DET connector into the left, 8-pin socket on the back panel of the detector.

Connecting the 7-Connector System Interconnect Cable

❖ To connect the 7-connector system interconnect cable to the LC modules

1. To connect the Surveyor Autosampler Plus plug the A/S connector into the left, 8-pin socket on the back panel of the autosampler as shown in [Figure 13](#).
2. To connect the pump, do the following:
 - To connect the Surveyor LC Pump Plus, plug a connector labeled PUMP into the 8-pin socket on the back panel of the pump.
 - To connect the Surveyor MS Pump Plus, do the following:
 - Connect the end of the adapter cable labeled PUMP to one of the PUMP connectors of the 7-connector interconnect cable.
 - Plug the end of the adapter cable labeled MS PUMP into the 8-pin socket on the back panel of the MS pump. The adapter cable for the MS pump is shown in [Figure 12](#).

Figure 12. Adapter cable for the Surveyor MS Pump Plus




3. To connect the Surveyor UV/Vis Plus Detector, plug a detector connector into the left, 8-pin socket on the back panel of the detector.

2 Installation

Making the Initial Back Panel Connections

Figure 13. Back panels of the modules of a Surveyor Plus LC system, showing the 7-connector system interconnect cable



Connecting and Installing the Flowcell

The Surveyor UV/Vis ships in any of three configurations, as listed in [Table 2](#).

Table 2. Shippable configurations of the Surveyor UV/Vis Plus Detector

Part Number	Configuration
Surveyor-UV0	with no flowcell
Surveyor-UV1	with 10 mm Standard Analytical flowcell
Surveyor-UV5	with 50 mm LightPipe flowcell

The standard analytical flowcell comes already installed in your detector. For more information on the standard analytical flowcell, see ‘[Standard Analytical Flowcell](#)’.

The LightPipe flowcell comes packed in a small, separate box within the shipping carton for the detector. This small box contains the LightPipe flowcell (with protective caps on each end) and a plastic bag containing the inlet and outlet tubing and FingerTight fittings. Unpack the flowcell and install it as described in ‘[LightPipe Flowcell](#)’.

Optional flowcells can be installed on the Surveyor UV/Vis Plus Detector configurations, according to the directions included in their packaging. See [Table 3](#) for a list of available optional flowcells.

Table 3. Optional flowcells for the UV/Vis detector

Part Number	Flow cell Type	Pathlength (mm)	Volume (μL)
9550-0103S	Kel-F Bio Cell	6	9.0
9550-0101S	3 mm Semi-Prep	3	4.5
9550-0100S	6 mm Analytical	6	9.0
9550-0150-01S	High Pressure Microbore	2	.25
9550-0197S	Titanium Bio Cell	6	9.0
9550-053S	3 mm Microbore	3	1.2
9550-0234S	10 mm Analytical (Standard on Surveyor-UV1)	10	15
9550-0265S	6 mm Microbore	6	5.0
9550-0263S	Cuvette Holder		
9550-0267S	High Pressure Standard	6	9.0
803237	50 mm LightPipe (Requires F3005-010 LightPipe mounting assembly)	50	10

Most of the flowcells listed in [Table 3](#) require flowcell cover P/N 60053-40005. The 50 mm LightPipe (P/N 803237) requires flowcell cover P/N 60053-40004. The High Pressure Microbore flowcell requires flowcell cover P/N 60053-40006.

Powering On the Detector for the First Time

After you have installed the flowcell, the detector is ready to be powered on.

❖ To power on the detector for the first time

1. Ensure that the power switch at the front of the unit is in the Off position (released or out position).

IMPORTANT If you are using a previous version of the LightPipe flowcell that has natural PEEK end caps, never turn the power on without the flowcell cover in place. When the instrument is powered on, it performs its Start-up diagnostics. Ambient light can affect these diagnostics.

IMPORTANT Do not operate the instrument in the Run or Preview Run mode before you install the flowcell cover, which minimizes thermal fluctuations.

2. Attach the power cord to the power entry module on the back panel of the detector and connect it to the power source.
3. Turn the power on by pushing the power button in to engage it.

When the instrument is powered on, it performs its Start-up diagnostics. After a few seconds, the LEDs on the front of the detector turn solid green. If one or more LEDs flash green or amber, or turn solid amber, there might be a problem. The COMM LED remains amber until the software is started and initiates communication with the detector.

See [Chapter 4, “Troubleshooting,”](#) for a complete description of LED status indications. Call your Thermo Fisher Scientific Service Representative if you require assistance.

Maintenance

This chapter describes the maintenance procedures that you must perform to ensure optimum performance of the UV/Vis detector.

The performance of your Surveyor UV/Vis Plus Detector depends on the maintenance of all parts of the detector. It is your responsibility to maintain your UV/Vis detector by properly performing the maintenance procedures on a regular basis.

Procedures described in this chapter are to be performed only by personnel trained in the service and maintenance of this instrument. If you have any questions on proper maintenance, or would like to arrange for a preventive maintenance program, please contact your Thermo Fisher Scientific representative.

Contents

- [Recommendations for Routine Maintenance](#)
- [Tools Required](#)
- [Cleaning Detector External Surfaces](#)
- [LightPipe Flowcell](#)
- [Standard Analytical Flowcell](#)
- [Cleaning Flowcells](#)
- [Replacing the Lamps](#)

Recommendations for Routine Maintenance

Table 4 lists recommendations for routine maintenance of the Surveyor UV/Vis Plus Detector. Use this table as a basis for developing your maintenance program in accordance with the practices at your company.

Table 4. Recommended Routine Maintenance

Procedure	Interval
Exterior Detector Housing Cleaning	As needed
Flowcell Cleaning	As needed*
Lamp Replacement	
a. (D ₂) Deuterium	a. Every 1000 h or as required
b. (W) Tungsten	b. Every 2500 h or as required

* Good laboratory practices dictate that the flowcell should be flushed with clean solvent after every use. This will greatly reduce the frequency with which you need to clean it.

Tools Required

The following tools are required for the maintenance procedures in this chapter:

- narrow-tip, flat-bladed screwdriver (2 mm wide)
- #2 Phillips screwdriver
- 1/4 in. open end wrench
- 1/8 in. hex wrench (supplied)
- 3/32 in. hex wrench (supplied)
- 9/64 in. hex wrench (supplied)

Cleaning Detector External Surfaces

The external surfaces of the detector should be kept clean and dry. To clean the outside of the detector, wipe it with a dust-free cloth or a damp cloth (moistened with water only) to remove dirt or stains.

LightPipe Flowcell

This section describes how to remove the LightPipe flowcell from, and install the flowcell in, the UV/Vis detector. Contact your Thermo Fisher Scientific Service Representative if you have any questions regarding LightPipe maintenance or replacement.

This section contains the following topics:

- [Removing the LightPipe Flowcell](#)
- [Storing the LightPipe Flowcell](#)
- [Installing the LightPipe Mounting Assembly](#)
- [Installing the LightPipe Flowcell](#)




CAUTION Do not disassemble the LightPipe housing or tighten the screws on the housing. Doing so will damage the LightPipe flowcell. Thermo Fisher Scientific cannot be held responsible for any damage done to the LightPipe by attempts to disassemble the housing or to tighten the screws. Please contact your Thermo Fisher Scientific Service Representative with any questions regarding LightPipe maintenance or service.

Removing the LightPipe Flowcell

❖ To remove the LightPipe from the mounting assembly

1. Open the front doors of the detector.
2. Using a 1/8 in. hex wrench (supplied), unscrew the LightPipe cover retention screw. Then remove the LightPipe cover. See [Figure 14](#).

Figure 14. Surveyor UV/Vis Plus Detector with the LightPipe and lamp covers on






CAUTION Do not touch the surface of the photodiode. When removing the photodiode mount, be sure to place the mount where the photodiode will not get dirty or scratched.

3. Loosen and remove the two photodiode mount retention screws using the 3/32 in. hex wrench (supplied). See [Figure 15](#).
4. Carefully pull the photodiode mount straight back and off the LightPipe mounting assembly. Rest the photodiode mount out of the way inside the detector housing or on the benchtop.

Figure 15. Photodiode assembly, showing two mounting screws



Photodiode assembly
mounting screws

5. Loosen the thumbscrew holding the retention block to the LightPipe retention bar. See [Figure 16](#).
6. Lift the retention bar shown in [Figure 17](#) until the bar is clear of the alignment notch in the LightPipe.
7. Carefully remove the LightPipe from the LightPipe mounting assembly without touching the ends of the LightPipe.



CAUTION Do not touch the ends of the LightPipe flowcell. This could damage the sensitive fiber optics. See [Figure 17](#). If you must grasp the ends of the LightPipe flowcell to remove it from the assembly, use clean, talc-free gloves.

Figure 16. LightPipe mounting assembly, showing retaining block





Figure 17. LightPipe mounting assembly, showing retention bar



Storing the LightPipe Flowcell

❖ To store the LightPipe flowcell

1. Remove the FingerTight fittings and any attached tubing.
2. Immediately replace the protective end caps on the ends of the LightPipe to protect the sensitive fiber optics. See [Figure 18](#).

Figure 18. LightPipe flowcell with protective end caps, ready for storage




Installing the LightPipe Mounting Assembly

If you are replacing a standard analytical flowcell with a LightPipe flowcell, you will need to install the LightPipe mounting assembly in the detector.

❖ To install the LightPipe mounting assembly

1. Remove the current flowcell. See [“Removing the Standard Flowcell Assembly”](#) on [page 29](#) for instructions on removing the standard flowcell.
2. Carefully align the two holes in the mounting assembly with the mounting screw holes in the detector housing.
3. Ensure that you can swing the retention bar forward into the body of the mounting assembly. See [Figure 19](#).
4. Using the 9/64-in. hex driver (supplied), screw the two LightPipe mounting assembly screws into the detector. Do not overtighten the screws.

Figure 19. LightPipe mounting assembly shown attached to detector



Installing the LightPipe Flowcell

The LightPipe flowcell assembly contains the parts listed in [Table 5](#).

Table 5. Parts for LightPipe flowcell assembly

Part Number	Description
803237	Flow cell Assembly, with inlet/outlet tubing and fittings (5 cm LightPipe)
2522-0285	FingerTight PEEK Ferrule Nuts
803260	Inlet tubing, with insulation, PEEK 1/16 x 0.005 in. ID (RED)
703950	Outlet tubing, PEEK 1/16 x 0.01 in. ID (Blue)

❖ To install the LightPipe in the mounting assembly

1. If the photodiode mount is attached to the end of the LightPipe mounting assembly, detach it as follows.
 - a. Loosen and remove the two photodiode mount retention screws ([Figure 15](#)) by using the 3/32 in. hex wrench (supplied).
 - b. Carefully pull the photodiode mount straight back and off the LightPipe mounting assembly. Place the photodiode mount out of the way inside the detector housing or on the benchtop.

2. Remove the protective end caps from the ends of the LightPipe.



CAUTION Do not touch the ends of the LightPipe. This could damage the sensitive optics. If you must grasp the ends of the LightPipe to remove it from the assembly, use clean, talc-free gloves.

3. Using the supplied fittings, attach the inlet tubing and the outlet tubing to the flowcell as follows:
 - a. Connect the insulated red PEEK™ inlet tubing to the port in the bottom of the flow cell.
 - b. Connect the blue PEEK outlet tubing to the port in the top of the flowcell.
 4. Being careful not to touch the ends of the flowcell, hold the flowcell with the alignment slot facing upward as you insert it into the mounting assembly.
 5. Secure the flowcell in the mounting assembly as follows:
 - a. Pull the retention bar down over the alignment slot in the flowcell (Figure 17).
 - b. Align the retaining block with the notch in the base of the assembly.
 - c. Slide the thumbscrew onto the retention bar. Then, push the retaining block firmly against the flowcell and tighten the thumbscrew finger tight (Figure 16).
 6. Reattach the photodiode mount as follows:
 - a. Carefully align the holes in the photodiode mount with the holes in the LightPipe mounting assembly. Be careful not to touch the photodiode while manipulating the photodiode mount.
 - b. Screw in and tighten the two photodiode mount retention screws (Figure 15) using the 3/32 in. hex wrench (supplied).
- IMPORTANT** To ensure optimum performance, do not run the Surveyor UV/Vis Plus Detector without replacing the LightPipe cover. The LightPipe cover minimizes thermal fluctuations. If you are using a previous generation LightPipe flowcell that has natural PEEK end caps, install the flowcell cover before you power on the detector.
7. Guide the tubing through the slots in the top and bottom of the flowcell compartment cover. Then secure the cover by screwing in the LightPipe cover retention screw (see Figure 14) with the 1/8 in. hex wrench (supplied).
 8. Connect the free end of the inlet tubing to the outlet of the LC column.
 9. Close the front doors of the detector.
 10. Before you start the pump flow, ensure that the outlet tubing drains to the waste reservoir.

Standard Analytical Flowcell

This topic describes the removal and installation of the standard analytical flowcell. For other flowcell problems, such as a cracked window or leaks that occur in locations other than at the inlet/outlet fittings, contact your Thermo Fisher Scientific service representative.

This section contains the following topics:

- [Removing the Standard Flowcell Assembly](#)
- [Installing the Standard Analytical Flowcell](#)



CAUTION Flow cells are factory-assembled units that should not be disassembled for any reason.

Removing the Standard Flowcell Assembly

To clean or replace the flowcell, you need to remove it from the UV/Vis detector.


❖ To remove the flowcell

1. Turn off the power to the detector power and disconnect the power cord from the back panel of the detector.
2. Open the front doors of the detector.

The flowcell assembly is located behind the flowcell cover.

3. Remove the flowcell cover by unscrewing the flowcell cover retention screw and pulling it straight out and away from the flowcell assembly. See [Figure 20](#).

Figure 20. Surveyor UV/Vis Plus Detector with the flowcell cover installed



4. Loosen and remove the two photodiode mount retention screws. See [Figure 21](#).
5. Carefully pull the photodiode mount straight back and off the flowcell assembly. Rest the photodiode mount out of the way inside the detector housing or on the benchtop. See [Figure 22](#).
6. To free the tubing from the UV/Vis detector, remove the thumbscrew that holds the tubing retainer bracket to the UV/Vis detector.
7. Using the 9/64 in. hex driver (supplied), unscrew the two flowcell assembly mounting screws from the detector.
8. Carefully pull the assembly out of the UV/Vis detector and toward you to remove it from the detector.
9. If the flowcell needs to be completely removed from the system, disconnect the flowcell inlet tube from the column and disconnect the flowcell outlet tubing from the waste reservoir.

[Figure 22](#) shows the detector with the flowcell assembly removed.

Figure 21. Photodiode mount and tubing retainer bracket





Figure 22. Surveyor UV/Vis Plus Detector with flowcell assembly removed



Installing the Standard Analytical Flowcell

❖ To install the standard analytical flowcell

1. With the inlet tube on the bottom left and the flowcell slot on top, slide the flowcell assembly into the detector and under the alignment pin.
2. Replace the two flowcell assembly mounting screws.
3. Align the photodiode mount with the flowcell assembly, so that the two holes with nylon washers align with the photodiode mounting posts on the flowcell assembly.
4. Install the two photodiode mount retention screws using the 9/64 in. hex driver.
5. Connect the inlet tubing to the chromatographic column and the outlet tubing to the waste reservoir.
6. Replace the tubing retainer bracket and tighten the thumbscrew just enough to hold the tubing in place.
7. Replace the standard analytical flowcell cover by guiding the tubing through the slots in the flowcell compartment cover and then securing the cover with the flowcell cover mounting screw.
8. Connect the power cord to the power entry module on the detector back panel.

Cleaning Flowcells

The LightPipe and the standard analytical flowcell need to be cleaned to maintain optimum performance. Good laboratory practice dictates that you flush flowcells with clean solvent after every use. This practice keeps cell contamination at a minimum, but occasionally the exterior and/or interior surfaces of the flowcell can become contaminated. Flow cell contamination is usually caused by precipitation of salts or by solubility problems. Signs of a contaminated flowcell are increased baseline noise, signal spiking, erratic or drifting baselines, low light intensity and/or increased backpressure.

This section contains the following topics:

- [Cleaning with Organic Solvents](#)
- [Cleaning with Nitric Acid](#)

Cleaning with Organic Solvents

If you suspect that your flowcell needs to be cleaned, start with the following procedure using organic solvents.

❖ To clean the flowcell with organic solvents

1. Disconnect the column from the chromatographic system to avoid column degradation. Connect the flowcell inlet directly to the chromatographic pump.
2. Flush the cell with water to prevent a reaction between the last solvent used in the chromatographic system and the cleaning solvent that is going to be used.

Note Ensure that the cleaning solvents you plan to use is miscible with the solvent already present in the flowcell and pump. Isopropanol is a good choice for most applications. If the last solvent in the pump was an aqueous buffer solution, be sure to pump 25 to 40 mL of HPLC-grade water (or equivalent) through the system to remove any salts *before* you flush the system with the cleaning solvents. This wash will help avoid salt precipitation problems.



CAUTION Thermo Fisher Scientific does not recommend using a syringe to force solvent through the flowcell. Pressurizing the syringe could cause a leak or rupture, resulting in a dangerous and uncontrolled spraying of solvent.

3. Flush the flowcell with 40 to 50 mL of cleaning solvent (HPLC-grade water, methanol, or isopropanol).
4. Flush the cell with water to prevent a reaction between the cleaning solvent and the mobile phase that is used in your application.

Cleaning with Nitric Acid

Methanol or isopropanol is generally sufficient for cleaning a flowcell. However, if the flowcell remains contaminated after flushing it with organic solvents, perform the following procedure using nitric acid.



CAUTION Nitric acid is a strong oxidizing acid, and it can react vigorously with alcohols (especially methanol). Be sure to wear protective clothing and eye protection and adhere to your company's safety procedures for the proper handling and disposal of corrosive acids. Flush the flowcell with water to remove all traces of alcohol before flushing it with nitric acid.

❖ To clean the flowcell with nitric acid

1. Completely remove the flowcell assembly from the detector housing by following the procedure in “[Removing the LightPipe Flowcell](#)” on [page 23](#) or “[Removing the Standard Flowcell Assembly](#)” on [page 29](#). (This will prevent possible leaks from harming the mechanical and electronic components of the detector.)
2. Ensure that the column is removed from the chromatographic system to avoid column degradation. The flowcell inlet should be connected directly to the chromatographic pump.
3. This step is very important! Flush the cell with water to prevent a reaction between the last solvent used in the chromatographic system and the cleaning solution that is going to be used.
4. Prepare a 20% (vol/vol) solution of nitric acid in HPLC-grade water.



CAUTION Thermo Fisher Scientific does not recommend using a syringe to force acid solutions through the flowcell. Pressurizing the syringe could cause a leak or rupture, resulting in a dangerous and uncontrolled spraying of acid.



CAUTION Before you pump cleaning solutions through the flowcell:

- Ensure that the LC column has been removed from the chromatographic system.
- Ensure that water was the last solvent in the LC pump and solvent reservoir.

5. Pump the nitric acid solution through the flowcell using the chromatographic pump.
6. After you have finished the cleaning procedure, pump another 25 to 40 mL of water through the flowcell to remove all traces of nitric acid. Monitor the pH of the outlet stream of the flowcell to ensure that the acid has been completely flushed out.
7. See “[Installing the LightPipe Flowcell](#)” on [page 27](#) or “[Installing the Standard Analytical Flowcell](#)” on [page 31](#) to reinstall your flow cell.



CAUTION Flow cells are factory-assembled units that must not be disassembled under any circumstances. If you encounter contamination problems that cannot be remedied by using this cleaning procedure, contact your Thermo Fisher Scientific representative to arrange for repair or replacement of the flowcell.

Replacing the Lamps

The Surveyor UV/Vis Plus Detector has two lamps, one tungsten (W) and one deuterium (D₂). These lamps need to be properly cared for to maintain the optimum performance of your detector. Keep them clean and free from grease or fingerprints. Oils from your fingers that are left on the glass in the form of fingerprints can cause the lamps to fail. The lamps need to be warmed up before use. Typically, the warm up time needs to be at least 1.5 h to get the best performance from your detector.

Note When the instrument is turned on, the lamps are automatically ignited. If the lamps have been turned off, you need to load a method and wait for the lamps to warm up before you proceed with a run.

As the lamps age, there is a reduction in light output that results in increased baseline noise. If the noise level on the output signal of your detector is increasing and cleaning the flowcell does not help, you should change the appropriate lamp, using the procedures in this section. Be sure to replace only one lamp at a time to make operation verification easier.

Note The deuterium lamp has a lifetime of approximately 1000 h and the tungsten lamp has a lifetime of approximately 2500 h. Lamp lifetime varies depending upon the application.

To access the lamps, open the front doors of the detector. The deuterium and tungsten lamps are located in the housing on the left. The lamps are supplied pre-aligned in their individual assemblies to make them easy to install and align.

IMPORTANT Never loosen the screws that hold the lamp to its assembly or attempt to rotate or move the lamp up or down in the assembly. Either of these actions can cause a loss of alignment and degrade the system performance.

This section contains the following topics:

- [The Deuterium Lamp](#)
- [The Tungsten Lamp](#)

The Deuterium Lamp

The deuterium (D₂) lamp typically requires a warm-up time of twenty to thirty minutes. However, for applications that demand great sensitivity, you need to allow a warm-up period of 1.5 h.

IMPORTANT The intensity of the deuterium lamp falls off very slightly over a period of time after it is turned on. You should wait at least 1.5 h for the lamp to stabilize after a cold start before collecting data in the spectral range of the deuterium lamp.

The deuterium lamp lifetime is usually at least 1000 h. The deuterium lamp assembly is equipped with a chronometer that tracks the total hours of lamp operation. To read the chronometer, note the position of the “gap” in the mercury tube against the graduated background.

IMPORTANT The surfaces of both lamps must be kept free of fingerprints and smudges. For this reason, the lamps should only be handled with clean, talc-free gloves. If either lamp requires cleaning, clean with a lint-free lens paper moistened with methanol or isopropanol before replacing the lamp cover.

This topic contains the following subtopics:

- [Removing the Deuterium Lamp](#)
- [Installing the Deuterium Lamp](#)

Removing the Deuterium Lamp

❖ To remove the deuterium lamp


1. Disconnect the power cord from the detector back panel and ensure that the instrument is turned off.
2. Open the front doors of the detector.



CAUTION Because intense UV light can damage your eyes, the detector automatically turns off the lamps when the lamp cover is removed. Because the lamps run hot, turn off the power and allow sufficient time (approximately 30 minutes) for the lamps to cool before you remove the lamp cover.


3. To remove the lamp cover, use an Allen wrench to unscrew the retention screw. See [Figure 23](#). Then pull the cover straight back to expose the lamp assemblies.

Figure 23. Surveyor UV/Vis Plus Detector showing the lamp cover retention screw



4. Unplug the deuterium lamp lead from the detector. See [Figure 24](#). Take care not to twist the connector as you gently pull it out.
5. Loosen the two thumbscrews that hold the lamp assembly in place and pull the assembly straight out.

Figure 24. Lamp housing removed




Installing the Deuterium Lamp

❖ To install a new deuterium lamp

1. Hold the deuterium lamp assembly so that the leads are at the top. Slide the assembly onto the alignment pin shown in [Figure 25](#).
2. Securely fasten the assembly in place with the two thumbscrews and aluminum standoffs.
3. Connect the lamp lead to the upper terminal in the lamp compartment.
4. Replace the lamp housing and secure it with the screw. Close the front doors of the detector.
5. Connect the power cord to the back panel of the detector.

Figure 25. Deuterium lamp assembly alignment pin



The Tungsten Lamp

The tungsten (W) lamp typically requires only 15 min of warm-up time. Its lifetime is approximately 2500 h. An increase in noise can indicate that the lamp needs to be replaced.

This topic contains the following subtopics:

- [Removing the Tungsten Lamp](#)
- [Installing the Tungsten Lamp](#)

IMPORTANT The surfaces of both lamps must be kept free of fingerprints and smudges. For this reason, the lamps should only be handled with clean, talc-free gloves. If either lamp requires cleaning, clean with a lint-free lens paper moistened with methanol or isopropanol before replacing the lamp cover.

Removing the Tungsten Lamp

❖ To remove the tungsten lamp


1. Turn the detector off and disconnect the power cord from the back panel of the detector.



CAUTION To avoid burning yourself, always allow sufficient time (approximately 30 min) for the lamp to cool before removing it.

2. Remove the lamp housing by removing the lamp cover retention screw (see [Figure 23](#)) and pulling the cover straight back to expose the lamp assembly.
3. Unplug the tungsten lamp lead from the detector, taking care not to twist the connector as you gently pull it out. See [Figure 26](#).
4. Loosen the thumbscrew and the aluminum standoff that hold the lamp assembly in place.
5. Pull the assembly straight out to remove it from the detector.

Figure 26. Detector with lamp cover removed, showing the tungsten lamp assembly



Tungsten lamp lead


Tungsten lamp assembly
thumbscrew

Installing the Tungsten Lamp

❖ To install the tungsten lamp

1. Hold the lamp assembly so that the leads are at the top. Slide the assembly onto the two alignment pins shown in [Figure 27](#). The alignment pins are located on either side of the detector monochromator aperture.
2. Securely fasten the assembly in place with the thumbscrew and aluminum standoff.
3. Connect the lamp lead to the lower terminal in the lamp compartment.
4. Replace the lamp housing and fasten it securely with the screw.
5. Connect the power cord to the back panel of the detector.

Figure 27. Tungsten lamp removed



Alignment pins

Troubleshooting

This chapter provides information for troubleshooting possible detector and chromatographic system problems. The information is organized in a table of symptoms, possible causes, and remedies. Because many of the problems attributed to the detector might actually be due to other components in the chromatographic system, references and potential solutions to these types of problems have also been included. Some common error messages and numeric codes are included to aid in troubleshooting errors.

Contents

- [Detector-Related Problems](#)
- [Error Messages](#)

Detector-Related Problems

Table 6 lists detector-related problems along with suggestions for corrective action.

Table 6. Troubleshooting detector-related problems (Sheet 1 of 5)

Symptom	Cause	Remedy
1. Spikes on baseline.	a. Continuous gas bubbles in the flowcell.	a. Degas the mobile phase. Connect a backpressure device to the flowcell (check backpressure rating).
	b. Immiscible solvent bubbles following mobile phase changeover.	b. Flush the flowcell with isopropanol, and then with mobile phase.
	c. Electrical interference.	c. Check the electrical lines for good connections and/or interference from broadcast radiation. Check for ground loops.
	d. Extremely large fluctuations in voltage on POWER line.	d. Remove systems (for example, ovens) that cause voltage fluctuations, isolate the detector to a “quiet” circuit, or use an uninterruptible power supply (UPS).
	e. Ambient light.	e. Install a flowcell cover. The flowcell cover must be installed before power On and must not be removed during Run and Preview Run operation.

Table 6. Troubleshooting detector-related problems (Sheet 2 of 5)

Symptom	Cause	Remedy
2. Random noisy baseline.	a. Contaminated flowcell.	a. Flush the flowcell with cleaning solvents as described in “Cleaning Flowcells” on page 32 . Check for leaks.
	b. Leak in sample inlet line.	b. Check all of the fittings from the column outlet to the flowcell inlet for leaks.
	c. Bubble trapped in flowcell.	c. Increase the flow rate until the bubble is removed. Supply a backpressure device to the flowcell (check pressure rating to avoid rupturing flowcell).
	d. Leaking flowcell.	d. Replace the flowcell.
	e. Insufficient lamp warm-up.	e. Allow a 30 minute warm-up for normal operation and 1.5 hours for maximum sensitivity.
	f. Lamp aging or defective.	f. Replace the lamp.
	g. Ground loop problem between integrator and detector.	g. Check for proper cable connections to the detector output; do not ground at both ends of the cable.
	h. Flow cell or lamps dirty.	h. Clean the dirty component as described in “Cleaning Flowcells” on page 32 .
	i. Integrator input voltage does not match detector output voltage.	i. Verify that the integrator is connected to the appropriate Analog Output connections on the detector. See “Analog Output Connections” on page 11 . Check the attenuation setting on the integrator.

Table 6. Troubleshooting detector-related problems (Sheet 3 of 5)

Symptom	Cause	Remedy
3. Excessive baseline drift Drift.	a. Flow cell contaminated.	a. Flush the flowcell with cleaning solvents as described in “ Cleaning Flowcells ” on page 32 . Check for leaks.
	b. Mobile phase contamination.	b. Replace with fresh mobile phase made with high-purity solvents.
	c. Material bleeding from column.	c. Clean or replace the column.
	d. Leaks in system or flow cell.	d. Check all of the fittings for leaks. Replace the flowcell.
	e. Tiny bubble trapped in flowcell.	e. Increase the flow rate until the bubble is removed. Connect a backpressure device to the flowcell outlet (check backpressure rating to avoid rupturing flow cell).
	f. Large temperature fluctuations.	f. Remove the system from drafts. Thermostatically control the column temperature.

Table 6. Troubleshooting detector-related problems (Sheet 4 of 5)

Symptom	Cause	Remedy
4. No peaks, or peaks much smaller than expected.	a. Incorrect wavelength setting.	a. Check the wavelength setting. Make sure that the correct file is selected.
	b. Lamp not On, or defective.	b. Make sure the lamp is On. Run the detector diagnostic tests to check the lamp. Replace the lamp if necessary.
	c. Integrator input voltage does not match detector output voltage.	c. Verify that the integrator is connected to the appropriate Analog Output connections on detector. (See “ Analog Output Connections ” on page 11.) Check the attenuation setting on the integrator.
	d. Insufficient sample reaching the detector.	d. Check the entire system for leaks. Verify the sample injection volume.
5. Broad, tailing peaks.	a. Rise time is too large (too slow).	a. Lower the rise time selection.
	b. Poor connection at flowcell inlet.	b. Check the end of the inlet tubing for a clean, flat surface free of obstructions.
6. Detector will not power up.	a. Tripped circuit breaker at power outlet.	a. Resolve the problem, and then reset the circuit breaker.
	b. Blown detector fuse.	b. Resolve the problem, and then replace the fuse.
	c. Incorrect voltage selected.	c. Verify that the Voltage Selector is set for the correct incoming line-voltage. (See “ Checking the Power Setting and Fuses ” on page 10.)
	d. Power cord not connected.	d. Connect the power cord.

Table 6. Troubleshooting detector-related problems (Sheet 5 of 5)

Symptom	Cause	Remedy
7. Detector does not go into run upon injection.	a. Detector not receiving Start signal.	a. Check the connection to the Run contacts on back panel.
8. LED1 - Power	a. Green - power On.	a. No action is required.
	b. Off - no power.	b. Check the fuse, line voltage, and so on.
9. LED2 - Comm	a. Green - communications normal.	a. No action is required.
	b. Amber - no communications with PC system.	b. Check the Ethernet connections and start the software. Check the module configuration in the software. Reload the firmware.
	c. Blinking Amber - firmware download in process or under serial control (stack ID set to 00).	c. Complete the download and re-establish communications with the PC. Reset the stack address to value other than 00.
10. LED3 - Run	a. Flashing Green - run in progress.	a. No action is required.
	b. Solid Green - ready.	b. No action is required.
	c. Blinking Amber - an error has occurred while performing a run.	c. Refer to the software to determine the nature of the error, or begin the run again.
11. LED4 - Lamps	a. Green - one or both lamps are On.	a. No action is required.
	b. Amber - both lamps are Off.	b. Turn On lamps from the software and allow 1.5 hours for warm-up.

Error Messages

Error messages are posted to the Status page of the data system that you are using to control your detector.

Note For further information, document the log entry and contact Thermo Fisher Scientific technical support.

Possible error messages include, but are not limited to, those listed in [Table 7](#).

Table 7. Possible error messages (Sheet 1 of 2)

Error Number	Meaning
2001	System reset
2002	System error (bus/address)
2004	Lamp case open
2011	Sample light is too high
2012	Low light on D2 lamp
2013	Low light on W lamp
2014	Volatile RAM test failed
2015	Non-Volatile RAM test failed
2016	12 volts failed
2017	12 volt unregulated failed
2018	-12 volts failed
2019	Motor voltage failed
2020	W lamp voltage or driver failed
2021	D2 lamp or driver failed
2022	5 volts failed
2023	Reference diode failed
2024	Sample diode failed (Could be exposed to light.)
2025	Wavelength home test failed
2026	Second wavelength home failed
2027	Analog linearity test failed
2028	D2 lamp failed / not present
2029	W lamp failed / not present
2030	Shutter failed
2031	Unknown shutter / lamp failure

Table 7. Possible error messages (Sheet 2 of 2)

Error Number	Meaning
2201	Invalid file number
2202	Invalid lamp state
2203	No replay while in a run
2205	Invalid current time
2206	Invalid lamp turn on time
2207	Invalid lamp turn off time
2208	Invalid time from ready time
2209	Invalid lamp on option
2210	Invalid lamp off option
2212	Invalid program time
2213	Invalid wavelength for program
2215	Invalid program type
2216	Invalid output range
2217	Invalid rise time filter
2218	Invalid auto zero time
2219	Invalid output range
2220	Invalid rise time filter
2221	Invalid run wavelength
2226	Invalid auto zero time
2227	Invalid rise time filter
2228	Invalid run time
2232	Invalid runs per wavelength
2234	Invalid offset
2235	Invalid zero on lambda change value
2236	Invalid channel 2 selection

Accessories and Replaceable Parts

This chapter contains the lists of accessories and replaceable parts that you can order from Thermo Fisher Scientific.

Contents

- [Replaceable Parts](#)
- [Optional Flowcells for the UV/Vis Detector](#)

Replaceable Parts

The following items are replaceable parts for the Surveyor UV/Vis Plus Detector:

Assembly, power entry w/switch	F3030-010
Connector, external	2110-0354
Cable, converter power	60053-63002
Cable, Main - Analog PCB	60053-63010
Cable, Main - Converter PCB	60053-63012
Cable, deuterium lamp	6040-0167
Cable, tungsten lamp	6040-0168
Cable, display LEDs	F5025-010
Cuvette holder assembly	9550-0263S
Cover, lamps	60053-40042
Cover, LightPipe	60053-40041
Door, right, assembly	60053-60068
Door, left assembly	60053-60069
Fuse, slo-blo, for 230 V applications (T1.25 A)	5101-1857
Fuse, slo-blo, for 100/115 V applications (T2.5 A)	5101-1858
Housing assembly	F3050-040S
Kit, accessory	60053-62004
Lamp, tungsten	9551-0022
Lamp, D2	9551-0023
Transformer, main	9551-0187-04
Patch cable, RJ-45, Cat 5, shielded	70111-63302
PCB, analog	F9080-010
PCB, Converter	F9060-060
PCB, Motherboard	F9055-010
Photodiode, pre-amp assembly	60053-63031
Screwdriver, hex, ball-end, 1/8 in., black	00725-00042

5 Accessories and Replaceable Parts

Optional Flowcells for the UV/Vis Detector

Screwdriver, hex, ball-end, 9/64 in., black.	00725-00044
Screwdriver, hex, ball-end, 3/32 in., black.	00725-00040
Solenoid, lamp shutter	9551-0221
Switch, magnetic proximity	5110-0058
Transformer, hi-voltage	9551-0200

Optional Flowcells for the UV/Vis Detector

The following flowcells and flowcell covers can be ordered from Thermo Fisher Scientific:

Kel-F bio cell	9550-0103S
3 mm semi-prep	9550-0101S
6 mm analytical	9550-0100S
High pressure Microbore	9550-0150-01S
Titanium bio cell	9550-0197S
3 mm microbore.	9550-053S
10 mm analytical (standard on Surveyor-UV1).	9550-0234S
6 mm microbore	9550-0265S
Cuvette holder	9550-0263S
High pressure standard	9550-0267S
50 mm LightPipe (requires F3005-010 LightPipe mount)	803237

Most of the above cells require flowcell cover P/N 60053-40005. The 50 mm LightPipe (P/N 803237) requires flowcell cover P/N 60053-40041. The high pressure microbore flowcell requires flowcell cover P/N 60053-40006.

Surveyor Firmware Upgrade Utility

From time to time, upgrades to Surveyor UV/Vis Plus firmware might become available. If you need to upgrade the firmware, or to replace firmware that has become corrupted, you can do so by using the Surveyor Firmware Upgrade Utility. This appendix provides instructions on how to use the Surveyor Firmware Upgrade Utility.

❖ To download firmware using the Surveyor Firmware Upgrade Utility

1. From Windows Explorer locate the firmware files you intend to download (ROM*.bin and UV*.bin):
 - For ChromQuest, these files are located in the Firmware folder on the ChromQuest CD.
 - For Xcalibur, these files are located in the following folder:
[drive]:\Xcalibur\system\Surveyor Firmware
2. Set the Unit ID on the back panel of the detector to the maintenance setting:
 - a. Ensure that the UV/Vis power is turned off.
 - b. Make a note of the current settings of the rotary switches on the back UV/Vis detector. You will return them to these positions at the conclusion of the procedure.
 - c. Set the switches of the UV/Vis to 00.
3. Ensure that your Surveyor UV/Vis Plus Detector is connected to your data system computer with the standard Ethernet connection.
4. Turn the power to the detector on.

The COMM LED blinks amber to confirm that the Surveyor UV/Vis Plus Detector is configured for downloading the firmware.


5. Start the Surveyor Firmware Upgrade Utility by double-clicking the Surveyor Firmware Upgrade Utility application. The Surveyor Firmware Upgrade Utility window opens. See [Figure 28](#).

6. In the Surveyor Firmware Upgrade Utility window, initiate communication between the utility and the Surveyor UV/Vis Plus Detector:
 - a. In the Device Type area, select the **Surveyor UV/Vis** option. See [Figure 28](#).



CAUTION Do **not** interrupt the firmware download process. Do **not** turn off the power to the Surveyor UV/Vis Plus Detector or close the Surveyor Firmware Upgrade Utility while the detector is connected to the Surveyor Firmware Upgrade Utility. If you lose power to the detector during the firmware download process, contact your Thermo Fisher Scientific representative. The firmware board will have to be reprogrammed.


Figure 28. Surveyor Firmware Upgrade Utility, showing the Surveyor UV/Vis option selected



- b. Click **Connect**.

The status display at the bottom of the Surveyor Firmware Upgrade Utility window indicates that a connection has been made. See [Figure 29](#).

Figure 29. Surveyor Firmware Upgrade Utility showing a successful connection to the Surveyor UV/Vis Plus Detector




7. In the File Names area, select the appropriate firmware files for the Surveyor UV/Vis Detector:
 - a. Click **Browse** to the right of the ROM box, and then select the ROM file ("ROM*.bin"). The filename appears in the ROM box. See [Figure 29](#).
 - b. Click **Browse** to the right of the APP1 box, and then select the APP1 file ("UV*.bin"). The filename appears in the App1 box. See [Figure 29](#).

IMPORTANT You must select both a ROM file and an APP1 file. Your Surveyor UV/Vis Plus Detector will not work properly if the files are not downloaded together.

8. Click **Download**, and wait for the download to finish.

The utility notifies you when it completes the download (see [Figure 30](#)). Once the download is complete, you can safely close the Surveyor Firmware Upgrade Utility and turn off the power to the detector.

Figure 30. Surveyor Firmware Upgrade Utility, showing a successful download



9. Close the Surveyor Firmware Upgrade Utility. The connection with the Surveyor UV/Vis Plus Detector automatically terminates.
10. Reset the Unit ID on the back panel of your detector for normal operation:
 - a. Turn off the power to the UV/Vis detector.

IMPORTANT Turn the UV/Vis power off before adjusting the rotary switches.

- b. Set the rotary switches back to their original positions.
- c. Turn on the UV/Vis power to resume normal operation.

Index

Numerics

- 5-connector interconnect cable
 - connecting 16
 - figure 15
- 7-connector interconnect cable
 - connecting 17
 - figure 16

A

- absorbance range 5
- ac line problems 42
- accessory kit 8
- adapter cable for the Surveyor MS Pump Plus (figure) 17
- alignment pins for deuterium lamp assembly (figure) 37
- analog outputs
 - full-scale settings 13
 - specifications 5

B

- baseline problems
 - noise, causes of 43
 - spikes, causes of 42

C

- Cautions
 - Do not touch the ends of the LightPipe flowcell lamp assembly 34
- cell dimensions 5
- cell pressure rating 5
- chromatogram, output settings 13
- column bleed 44
- connections
 - analog output 11
 - Ethernet 11
 - remote communications inputs 13

D

- deuterium lamp
 - installing 37
 - lifetime 34
 - removing 35
- dimensions 5
- Document Survey button xii
- drift 5, 44

E

- electrical interference, causes of 42
- electromagnet compatibility iv
- error messages, reference 47
- errors in the manual, reporting xii

F

- FCC compliance iv
- Firmware downloads
 - cable connections 51
 - location of firmware files 51
 - Surveyor Firmware Upgrade Utility 51
- flow cell
 - bubbles in the 42, 43
 - cover (figure) 29
 - covers 50
 - installing the 31
 - leaks 43
 - location of 29
 - options 50
 - removing the 29
 - safety precautions when cleaning 33
 - schematic 3
- fuse holder removal 10
- fuse sizes 10

G

- ground loop problems, causes of 43

I

initial power on 20
 installation
 accessory kit 8
 checklist 8
 flow cell 31
 positioning the detector 9
 powering on 20
 tools for 9
 unpacking 8

L

lamp cover (figure) 36
 lamp housing (figure) 36
 lamps
 cleaning the 35, 37
 installing the deuterium lamp 37
 installing the tungsten lamp 39
 removing the deuterium lamp 35
 removing the tungsten lamp 38

LEDs

Comm 46
 errors 46
 Lamps 46
 Power 46
 Run 46

LightPipe

installing 27
 ready for storage 26
 removing 23
 removing from assembly 23

linearity 5

M

maintenance
 cleaning external surfaces 22
 cleaning flow cell 32
 installing flow cell 31
 overview 21
 removing flow cell 29
 routine maintenance recommendations 22
 maintenance tools 22

N

nitric acid, safety precautions 33
 Notes
 lamp lifetimes 34
 lamp warm up 34
 lamps, cleaning the 35, 37

O

operating humidity 5
 operating temperature 5
 optical bench (figure) 2
 optical fiber (figure) 4

P

parts, replaceable 49
 peak problems
 broadening 45
 size 45
 tailing 45
 photodiode mount (figure) 30
 positioning the detector 9
 power requirements 5
 product certification 5

R

regulatory compliance iii
 remote controls 5
 rise time
 incorrect setting, results of 45
 settings 5
 rotary switches, unit ID settings 13

S

safety standards iii
 scan rate 5
 short term noise 5
 specifications 5
 spectral bandwidth 5
 stack address (note) 13
 storage temperature 5
 Surveyor UV/Vis
 installation 8
 optical bench 2

T

temperature fluctuations, results of 44
 tools
 needed for installation 9
 needed for maintenance 22
 troubleshooting
 error messages 47
 table of common problems 42

tungsten lamp [34](#)
 installing [39](#)
 lifetime [34](#)
 removing [38](#)

U

unit ID (note) [13](#)
unpacking your detector [8](#)
UPS, uninterruptible power supply [42](#)
UV/Vis, troubleshooting [42](#)

V

voltage
 ac line problems [42](#)
 changing the setting [10](#)
 voltage selector [10](#)

W

warm-up [5](#)
wavelength accuracy [5](#)
wavelength range [5](#)
weight [5](#)

