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When the shutter is closed the generator will switch off, unless another shutter is still open in one of the other tube shield positions. Alternatively, the generator can be switched off manually by operating the GENERATOR 0 push-button. All lamps will extinguish. Before the generator is switched off the power to the X-ray tube should be reduced, by first turning the mA control to minimum, and then the kV control. Any generation of HT will be inhibited if the shutter does not fully close after a close shutter command.


Fig. 3.1. Opening the window shutter.

3.3. SHUTTER CONTROL

For camera work it is usual to control the exposure time of the film to X-rays by using a PW 1722 timer. Two timers can be fitted: the left-hand timer will control the shutter of either window 1 or window 2, and the other timer the shutter of window 3 or 4. This paragraph describes the automatic and manual control of the shutter operation.

When the camera is installed on its bracket in front of the required exit window of the tube shield, open the shutter disc and push the camera forward until the collimator enters the hole in the disc. Ensure that the microswitch button on the tube shield is actuated by the pin or arm of the camera.

NOTE: The space between the front of the collimator and the exit opening of the tube shield should be kept as small as possible to ensure that radiation does not escape. Take care that the filter disc is free to move, and is not damaged.

The *filter disc* is rotated to one of six marked positions, and for five positions a filter can be placed in the X-ray beam to be used as a β -filter for the various X-ray tube targets, as follows:

Filter - Zr, used for Mo radiation

Filter - Ni, used for Cu radiation

Filter - Fe, used for Co radiation

Filter - Mn, used for Fe radiation

Filter - V, used for Cr radiation

A filter is not used in the sixth position, and the window is open.

The *shutter disc* is pivoted on its centre and is spring loaded. When a camera or goniometer is fitted, the disc is rotated against the spring to allow the X-ray beam to pass.


- A white button, which is depressed by fitting a camera or a goniometer. The button actuates a microswitch to complete a safety circuit, and this enables the associated window shutter to be opened.
- A window shutter, operated by an electro-magnet, to allow the X-ray beam to pass into the camera or goniometer. The shutter operation is selected by the front panel switch and push-buttons, and the shutter will only open when the safety circuit is complete. If the shutter does not fully close after a close shutter command, then generation of HT to the X-ray tube is inhibited.
- A red lamp cover over two fail-safe lamps, which are lit to warn that the associated shutter circuit is activated (shutter is open). The shutter closes if a lamp fails.
- Two coloured indicator flags, one on each side of the shutter mechanism housing. The flags are mechanically coupled to the shutter, and a green flag appears when the shutter is closed, and a red flag when it is open.
- Four screw holes for mounting the camera bracket slide.

The point focus beam appears at windows 1 and 3.

The line focus beam appears at windows 2 and 4.


2.3. Rear panel

The layout of the rear panel is shown in figure 2.4.

RD 2541

Fig. 2.4. Layout of the rear panel.

RD 2939

Fig. 2.2. Front panel controls.


Fig. 2.1. Front view of PW 1720/00 X-ray generator.

PHILIPS

Table-top X-ray Generator PW 1720

Instruction Manual

9499 300 12011

771111

PHILIPS


Table-top X-ray Generator PW 1720

Instruction manual

9499 300 12011

771111

CONTENTS

	<i>Page</i>
1. INTRODUCTION	1-3
2. DESCRIPTION OF THE UNIT	2-1
2.1. Front panel controls	2-2
2.2. Tube shield	2-3
2.3. Rear panel	2-4
2.4. H.T. generator	2-5
2.5. X-ray tube	2-5
2.6. Accessories	2-5
3. OPERATION	3-1
3.1. Initial checks	3-1
3.2. Sample analysis	3-1
3.3. Shutter control	3-2
3.3.1 With timer control	3-3
3.3.2 Without timers	3-3
3.4. Switching off a motor when the shutter closes	3-3
3.5. Use of the window	3-4
3.6. Changing the X-ray tube	3-5
3.7. Positioning the tube shield	3-6
3.8. Diffractometry operation	3-6
4. SERVICE INFORMATION	4-1
4.1. General	4-1
4.2. Technical specification	4-1
4.3. Safety precautions	4-2
4.4. Installation	4-3
4.5. Installation of optional accessories	4-4
4.6. Water cooling system	4-7
4.7. Maintenance and adjustments	4-8
5. CIRCUITS DESCRIPTION	5-1
5.1. General	5-1
5.2. Switch on/off	5-1
5.3. H.T. generator	5-2
5.4. Safety circuit	5-2
5.5. Tube current control	5-4
5.6. Shutter control	5-8
6. PARTS LIST	6-1
6.1. Spare parts	6-1
6.2. PW 1728/20 Spare parts kit	6-6
APPENDIX 1 X-RAY DIFFRACTION TUBE POWER RATINGS	A-1
APPENDIX 2 PW 1720/25 H.T. GENERATOR WITHOUT TUBE SHIELD*	A-3

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