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
PW 1156 Metallurgical Kit with PW 1162 Stress Attachment

- ☐ For use with horizontal or vertical goniometer
- ☐ Increases usefulness of standard diffractometer to the metallurgist
- ☐ Accepts specimens as large as 10 cm (4 inches) in any dimension
- ☐ Allows precise measurement of residual stress
- ☐ Includes calibrated bending device for measuring additional stress and absolute stress
- ☐ Makes possible measurement of local stress through use of a pin hole collimator

These two attachments, in combination, provide the metallurgist with a simple, precise means of making X-ray diffraction studies on bulky specimens.


This kit includes a base plate to allow the goniometer to be moved back from the center line of the X-ray beam about 4.75 cm beyond its normal position, to accommodate bulky specimens. A fully adjustable extension fixture allows the slit assembly to be moved forward. A standard powder specimen holder with extra-long shaft is provided along with a large specimen holder. A special fixture is used to bring the specimen surface exactly into the goniometer axis. It consists of a centering pin that fits into the specimen holder shaft.

The stress attachment consists of a track and slide which connects the receiving slit-detector assembly to the goniometer arm. This replaces the block used for this purpose in the standard Metallurgical Kit. It is adjustable so that the receiving slit moves exactly along a radius of the goniometer. This unit allows the precise determination of the position of the diffraction line at two psi angles, and thus may be used to provide precise determinations of residual stress.

Also included is an attachment which allows the specimen holder shaft to be rotated rapidly to the desired psi angle, and another attachment which provides for calibration of additional stresses applied to the specimen while on the goniometer. Mounted on the specimen table, it is used to apply known stresses to a flat specimen, and thus determine the stress constant of the material.

For transmission diffractometry, the accessory offers a number of geometrical advantages. Intensity loss caused by the monochromator is partially compensated by larger apertures. Much smaller volumes are required for transmission specimens compared with the standard diffractometer set-up. For reflection patterns it is not necessary to remove the specimen holder. The specimen can be tilted over 90° , the crystal removed, and a normal arrangement for reflection work obtained, without extra specimen preparation and mounting.

OPERATING CONDITIONS

	Transmission	Reflection
Tube	Cu 50 kV — 40 mA	Cu 50 kV — 40 mA
Filter	None	Ni
Detector	Proportional + Discriminator	Proportional + Discriminator
Scale	4×10^2 cps	4×10^2 cps
Divergence slit	2°	1°
Scatter slit	None	1°
Receiving slit	0.3 mm	0.2 mm
Scanning speed	$\frac{1}{2}^\circ 2\theta/\text{min}$	$\frac{1}{2}^\circ 2\theta/\text{min}$
Paper speed	400 mm/h	400 mm/h
Time constant	4 s	4 s

PW 1075 Curved Crystal Focalizer

Typical applications of this accessory are studies of texture, measurement of large c-spacings, measurement of angular errors due to deviations of specimen surface, line profile studies, and investigations of limited quantities of material.

The accessory permits comparison of a transmission pattern with a reflect on pattern from a standard diffractometer. Effects of fluorescent radiation are eliminated.

The curved crystal focalizer consists of a flat plate of quartz (1011) elastically bent by two adjustment knobs on the crystal holder. The surface approximates a section of a logarithmic spiral.

The holder is mounted on the detector arm of the standard diffractometer, with the counter tube moved further out. Counter tube and specimen rotate in a 2 : 1 ratio as in the standard instrument. Focusing is adjusted by changing the curvature of the crystal. (Geometry is illustrated in figure 4).

☐ For use with vertical or horizontal goniometer

☐ For transmission diffractometry in a number of applications


Figure 4.


Figure 4

Sp — transmission specimen
M — curved crystal focalizer
F — virtual image of the source

α — angular aperture of the primary beam
CT — counter-tube
MFC — monochromator focussing circle
SS1 and SS2 — soler slits
DS — divergence slits
RS — receiving slit
FC — focussing circle
GC — goniometer circle

X-RAY OPTICAL ARRANGEMENT

The X-ray optical set-up has been designed for use with one of the line foci of a modern X-ray diffraction tube.

The angular aperture of the beam in the plane of Bragg focussing (sometimes called Brentano para-focussing) is indicated by the divergent lines. It is defined by a single divergence slit (D) which also limits the primary beam to the specimen area.

Apertures of various sizes are available from 2 minutes, for large d-spacing measurement, to 4 degrees for the back reflection region where the area of sample irradiated by the incident beam is small. A 1° aperture is normally used for general work in the front reflection region. The intensity increases with increasing angular aperture and the resolution is practically unaffected. The instrument is designed for flat specimens (S), 15 x 20 mm. Smaller specimens may be used by reducing the size of the divergence slit.

The receiving slit (R), in front of the Receiving Parallel Slit Assembly (RP), defines the width of the converging reflected beam from the sample to the detector. The resolution decreases slightly and peak intensity increases markedly when the slit width is increased. Various sizes are available.

Equally spaced thin metal foils (Parallel Slit Assembly) limit the divergence of the beam in the plane perpendicular to Bragg focussing to an angular aperture of 2 degrees $\pm$ 2 minutes.

Two sets are used; one between the X-ray source and specimen (P) and the other (PP) between the specimen and detector (SS). These slits make it possible to obtain high resolution and excellent line shape for the extended line source.

A scatter shield cylinder surrounds the specimen and a scatterer slit (H) is placed between the receiving parallel slit assembly and detector, to ensure the detector receiving only those rays diffracted from the specimen surface.



Step scanning

Step scanning with both goniometers becomes possible with:

- Step scanning device PW 1053/00, consisting of a mechanical unit with gear wheels for stepwise angular advancement of the goniometer with 0.01° , 0.02° and 0.05° (2°) per step
- Step scan control unit PW 1364/00, which provides the necessary start/stop signals for the stepwise advancement of the goniometers PW 1050 and PW 1380.

When the above two units are extended with the numerator PW 4253, printer control PW 4255 and the high speed digital printer PW 4202, the numerator produces a two digit sequence number by which analogue or digital measuring results or the individual steps can be identified.

Fig. 1 X-ray optical arrangement of both horizontal and vertical goniometers

TECHNICAL FEATURES

Precision of angular setting (2° Θ and ϕ)

$\pm 0.0025^\circ$

Reading accuracy of angular setting

$\pm 0.002^\circ$

Scanning

Unidirectional or oscillatory between

-38° to $+175^\circ$ (2°)

Scanning speeds

$1/8^\circ$, $1/4^\circ$, $1/2^\circ$, 1° , 2° (2°) per minute by using change gears PW 1059/00, $1/01$, $1/02$ which are included in the PW 1050/25

Scanning stops

Quickly adjustable limit switches for limiting scanning range to required value; also prevent overtravel at extreme of scale

Step scanning

Mechanical step scanning device

PW 1053/00 can be mounted so that goniometer advances at 0.01° , 0.02° or 0.05° per step, enables the measurement of fixed count line profile

Goniometer radius, i.e. also distance between focal spot tube and specimen

170 mm

Motor

Special synchronous fractional horse-power type

Alignment

Easily performed with the aid of fixtures, gauges and a standard silicon sample which are part of the equipment; with new set of alignment tools PW 1077 the alignment procedure is more accurate, simpler and shortened considerably.

PW 1050/25 Vertical Goniometer

- ☐ A wide-range high- and low-angle instrument
- ☐ High precision achieved by special worm drive system
- ☐ Unique, high-resolution optical arrangement provides extremely accurate measurement of diffraction angles and line intensities

The vertical PW 1050/25 wide-range goniometer is especially suited to standard powder diffractometry, in which the specimen and radiation detector are rotated and positioned in a fixed 1 : 2 ratio.

DESCRIPTION

Two concentrically-mounted shafts are driven in a 1 : 2 ratio by a precision worm wheel and worm drive. The radiation detector probe is attached radially to the worm wheel.

With each revolution of the worm drive, the detector moves 1° . A dial is fixed to the worm drive for direct readings of the angle two-theta. The reading accuracy of angular setting is $\pm 0.002^\circ$.

The specimen is placed in a holder carried by a hollow shaft positioned within the hollow bearing of the worm wheel and coupled to the wheel by a set of accurate herringbone gears. These gears rotate the specimen or crystal at exactly one-half the angular speed of the scanning arm carrying the detector probe. A set of interchangeable spur gears supplied with the goniometer allows a variety of scanning speeds.

Alignment time, using the Phillips alignment tool set, is less than 30 minutes. Full X-ray protection is provided during alignment.



Diffractometer set

When used for X-ray powder diffractometry, the goniometer PW 1050 must be completed with the PW 1049/10 diffractometer set. This set includes all accessories which are required for doing powder diffractometry with the PW 1050.

For the PW 1380 the diffractometer set is included in the typenumber PW 1380.

Radiation protection








The radiation protection shield is an important part of the diffractometer set. It completely shields off the X-ray path during operation of the diffractometer; shuts off X-ray beam when cover shield is removed.

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