CHAÎNE HPLC

FICHE N° 4041

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

Part No. 0993-8912

Notice

The information contained in this document is subject to change without notice.

PERKIN-ELMER MAKES NO WARRANTY OF ANY KIND WITH REGARD TO THIS MATERIAL, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. Perkin-Elmer shall not be liable for errors contained herein for incidental consequential damages in connection with furnishing, performance, or use of this material.

Copyright Information

Reproduction or publication in any form or format is prohibited without written permission of The Perkin-Elmer Corporation or any of its subsidiaries.

Copyright © 1995 The Perkin-Elmer Corporation.
All rights reserved.
Printed in the United States of America.

Release Information

Manual Part No.	Release	Release Date
0993-8912	C	September, 1995

Trademarks

Perkin-Elmer is a registered trademark of The Perkin-Elmer Corporation.
PE XPRESS is a trademark of the Perkin-Elmer Corporation.
Teflon, Tefzel, and Vespel are Trademarks of E.I. duPont de Nemours.

Registered names, trademarks, etc., used in this document, even when not specifically marked as such, are protected by law.

Company Name and Address

The Perkin-Elmer Corporation
Instrument Division
761 Main Avenue
Norwalk, Connecticut 06859-0010


Figure 1-2. Front Right Side View of Autosampler.



Figure 7-9. Sample Transfer Tubing Path

PERKIN ELMER

The Perkin-Elmer Corporation, 761 Main Ave., Norwalk, CT 06859-0012 U.S.A.
Tel: (203) 762-1000. Fax: (203) 762-6000
Bodenwerk Perkin-Elmer GmbH, Postfach 10 17 61, D-88647 Überlingen,
Germany. Tel: (07551) 810. Fax: (07551) 1612
Perkin-Elmer Ltd., Post Office Lane, Beaconsfield, Bucks HP9 1QA, England.
Tel: (0494) 676161. Fax: (0494) 679331/3

Part No. 0993-8910

Notice

The information contained in this document is subject to change without notice.

PERKIN-ELMER MAKES NO WARRANTY OF ANY KIND WITH REGARD TO THIS MATERIAL, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. Perkin-Elmer shall not be liable for errors contained herein for incidental consequential damages in connection with furnishing, performance, or use of this material.

Copyright Information

Reproduction or publication in any form or format is prohibited without written permission of The Perkin-Elmer Corporation or any of its subsidiaries.

Copyright © 1997 The Perkin-Elmer Corporation.
All rights reserved.
Printed in the United States of America.

Release Information

Manual Part No.	Release	Release Date
0993-8910	E	April 1997

Trademarks

Perkin-Elmer is a registered trademark of The Perkin-Elmer Corporation.
PE XPRESS is a trademark of The Perkin-Elmer Corporation.
Luer-Lok is a trademark of Becton Dickenson.
Teflon and Tefzel are trademarks of E.I duPont de Nemours.

Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are protected by law.



Figure 7.22. Exploded View of the Station Head Pump.

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