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1. Description

The apparatus is used for demonstrating gravitation and for determining the gravitation constant. It is a highly sensitive torsion balance and readings are taken by means of spot of light.

The rectangular metal casing is covered by glass plates so that the movable system can be seen and is fitted underneath with a stand column 12 mm. in diameter and on top with the carrier tube for the torsion wire. The carrier tube supports the adjustable and fixable torsion head to which the 25 cm. long torsion band made of bronze with a 0.01 x 0.15 mm. cross-section is fastened. This carries a L-shaped suspension with a concave reflector of approximately 30 cm. focal length fastened to it and two 15 g. leaden balls each 5 cm. from the axis of rotation. The space for this suspension in the casing has been kept as small as possible to avoid air currents, and is subdivided by plastic panes.

The suspension is locked by the two knurled-head screws underneath the casing. This is indispensable when transporting the apparatus. In the case of a permanently installed gravitation torsion balance a locking of the suspension is only then necessary when interference with the apparatus by unauthorized persons cannot be avoided.

Underneath the rectangular casing there is a rotary holder on which two 1.5 kg. leaden balls are placed. These balls can be approached to each of

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Notes

- The five-figure numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.
- The data DK ... refer to the experiment descriptions published as a collection of leaflets in the "LEYBOLD PHYSIKALISCHE HANDBÜCHER".
- The data and illustrations are not binding in every detail for the design of the apparatuses. It is our sole aim always to adapt our manufacturing programme to the most recent knowledge gained in all scientific and technical fields.

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