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1. Description

The counter is used in conjunction with the Wulf electroscope (546 01 A) to detect individual α -rays emanated from weakly radioactive preparations.

The apparatus is designed as a metal cylinder, approximately 5 cm long, with one side open. The other side is sealed by means of an insulator. The insulator is fitted with a fine-pointed needle connected to a 4 mm socket. For assembling, the cylinder is equipped with a short rod and a plug. Other accessories to the apparatus are a needle with a very weak radium preparation, stored in a small glass tube, and a small book of aluminium foils used to seal the opening of the counter.

2. Assembly

The leaf in the system of the Wulf electroscope serves to indicate the α -particles. It should be adjusted in

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Notes

- The five-figure numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.
- The data DC . . . refer to the experiment descriptions published as a collection of leaflets in the "LEYBOLD PHYSICS LEAFLETS".

Published to date:

	Cat. No.
1st series (72 sheets)	599 21
2nd series (52 sheets)	599 22
3rd series (49 sheets)	599 23

- The data and illustrations are not binding in every detail for the design of the apparatus. It is our sole aim always to adapt our manufacturing programme to the most recent knowledge gained in all scientific and technical fields.

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