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

THE DUNNING COLORIMETER

Colorimeter Box

The colored solutions in the Dunning Colorimeter have proved to be stable. If, however, at any time they appear to have faded, the ampules should be returned for examination and correction. A slight precipitate may sometimes be noticed, due to the action of the alkaline solution on the glass; this does not affect the color.

Delivered for \$5.00

HYNSON, WESTCOTT & DUNNING
BALTIMORE

Size
 $2\frac{1}{2} \times 3\frac{1}{2} \times 5$
inches

Complete Apparatus in Polished Wooden Box

The Dunning Colorimeter

FOR ESTIMATING THE QUANTITY OF PHENOLSULPHONE-
PHTHALEIN EXCRETED WHEN APPLYING THE
Rowntree and Geraghty Renal Functional Test

HYNSON, WESTCOTT & DUNNING

Diagnostic Agents and Appliances

BALTIMORE,

MARYLAND

THE DUNNING COLORIMETER

For estimating the quantity of Phenolsulphonephthalein excreted when applying the

ROWNTREE AND GERAGHTY RENAL FUNCTIONAL TEST

DIRECTIONS: The patient's bladder should be completely emptied and one exactly measured mil (1 c. c.) of sterile physiologically standardized phenolsulphonephthalein solution from an ampule injected, either intravenously or intramuscularly, as preferred; the former method requiring much less time for completing the test than the latter. *The ampules contain more than one mil*, to permit the withdrawal of the exact quantity required. To open ampules, scratch the neck just between the small bulb and the body of the ampule with the accompanying file, and break off the top with the forefinger and thumb.

INTRAVENOUS INJECTION: Collections of urine and estimations may be made in fifteen minutes, thirty minutes or one hour. The average normal eliminations are: 35 % to 45 % in fifteen minutes; 50 to 65 % in thirty minutes and 65 % to 80 % in one hour.

INTRAMUSCULAR INJECTIONS: The solution should be injected into the lumbar muscles and collections made of the urine excreted during the first hour and ten minutes, the dye appearing in the urine in about ten minutes after intramuscular injections. Collections of the urine excreted during the second hour should also be made. The Phenolsulphonephthalein excreted during each period should be separately estimated colorimetrically. The average normal eliminations of the intramuscular injections are 50 % the first hour and a total of 80 % in two hours.

ESTIMATIONS: To estimate the quantity of phenolsulphonephthalein excreted, using the Dunning Colorimeter: Dilute the specimen to about 200 mls (c. c.) with water and render it alkaline by the addition of 10 mls (c. c.) or more, if necessary, of a 5 % solution of sodium hydroxide, then farther dilute the alkaline urine with sufficient water to make it measure one litre. Enough of this dilution should be perfectly clarified, by filtration, to fill the open ampule. This ampule, containing the specimen, should then be placed in the center hole of the comparison box and compared with the sealed and marked test ampule until the one that most nearly matches the specimen in color is found; the number on the ampule indicating the percentage of phenolsulphonephthalein excreted. From this comparison, the phenolsulphonephthalein content of the specimen may be approximately estimated. More accurate comparisons may be made by using two test ampules at a time, one on either side of the specimen.

NOTE: *If, after adding the alkali to the specimen, the coloration is slight, showing small excretion of the dye, then the dilution should be carried only to 250 or 500 mls (c. c.) and the reading divided by 4 or 2, as the extent of dilution may demand.*

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