

PRÉSERVER
SAUVEGARDER
VALORISER

PULVERISATEUR - DEVILBISS

FICHE N° 12869

Période de fabrication : -

Fabricant : fabricant non renseigné

Domaines : Santé

Sous-domaines : Pharmacologie

Organisme : Université d'Auvergne, Faculté de médecine et de pharmacie

Ville : Clermont-Ferrand

Modèle :


Matériaux :


Description

C'est une boîte rectangulaire en carton qui s'ouvre par le dessus. Elle contient un flacon cylindrique en verre ainsi qu'une fine tige en inox. Il manque la poire et l'embout de la tige en or.

Utilisation

On devait fixer la poire sur la tige, puis appuyer dessus pour vaporiser la solution médicamenteuse contenue dans le flacon. La grandeur de la tige permettait d'atteindre le fond de la gorge ou du nez. Ce vaporisateur est utilisé pour le traitement des infections nasales et buccales mais aussi pour le traitement des maladies pour animaux.







To Purchasers of DeVilbiss Atomizers

DIRECTIONS AND GUARANTY

TO CLEAN—Hold finger over hole in spray point at same time compressing bulb vigorously, thus blowing back into the bottle any sediment that may have lodged in either the spray point or outlet tube.

If this does not suffice, spread spray tube apart at end—remove the tip—and clean out the small hole on each side. In replacing the tip, be sure the word “up” is on top.

TO STERILIZE—Pass spraying tip through a flame; or unscrew glass outlet tube and boil metal parts and bottle.

When spraying the nose, use nasal guard over spray tip for best results. When adjusting nasal guard, slip guard over spray tip until end of spray tip protrudes through opening.

When spraying throat the nasal guard should be removed in order that the spray tip may be correctly adjusted, to permit proper direction of the spray.

For Additional Directions See Reverse Side

Pour nous citer :

Base de la Mission nationale de sauvegarde et de valorisation du patrimoine scientifique et technique contemporain, PATSTEC, Pulverisateur - DeVilbiss (fabricant non renseigné), <https://www.patstec.fr/ressources/objets/detail?id=14976>, consulté le 2026-07-27